

 	<b>F722</b> <b>TECHNICAL SHEET</b>	<b>Edition</b>	<b>Rev</b>	<b>Editor</b>
		1/2014	01	RAQ-MKT

Ref. code: IMM120071

# INFANT TRANSPORT UNIT PLUS



DESCRIPTION
<ul style="list-style-type: none"> <li>• Newborn emergency transportation system</li> <li>• Complete with inner_xing belts</li> <li>• Grants a warm, protected and safe environment</li> <li>• Innovative foldable-bag system, PVC material allowing, when not used, to be taken everywhere, always ready for use</li> <li>• 3 transparent sides allowing complete and continuous monitoring of newborn by medical personnel</li> <li>• Internal mattress ensuring comfortable transport</li> <li>• 2 external handles for transportation</li> <li>• Arranged for the anchorage tot the ambulance stretcher (EN 1789:1999)</li> <li>• Dimension open: cm 90x38x34,5(h)</li> <li>• Dimensions closed: cm 90x38x15(h)</li> <li>• Supplied with transport bag</li> </ul>



<b>PRODUCT CLASSIFICATION</b>	Class I, accordingly to MDD 93/42/CEE and further updates
<b>PACKAGING</b>	Singulary packaged

PRODUCT CONFORMITY/REFERENCE NORMS	
MDD93/42/CEE	Directive concerning medical devices
ISO14971	Medical devices risk management application
ISO 15233-1	Symbols to be used with medical devices labels

Rev. 11/2017

Oscar Boscarol srl, via Enzo Ferrari 29, Bolzano – IT  
 Tel. +39 0471932893 – Fax +39 0257760140  
[www.boscarol.it](http://www.boscarol.it) - [info@boscarol.it](mailto:info@boscarol.it)