



Emergency Medical Systems







## INDEX

	PAGE
MEDICAL SUCTION UNIT	6
EMERGENCY BAGS & RUCKSACKS	24
IMMOBILIZATION	42
RESUSCITATOR BAGS	62
DIAGNOSTICS	72
REANIMATION	86
MEDICATION	96
OTHERS, ACCESSORIES, DIDACTICS	142
CONTACTS AND SELLING CONDITIONS	148

Dear Customer,  
herewith we inform you that due to our Ministry of Health instructions related to the advertisement of medical devices, diagnostics in vitro and surgical medical devices, the featuring of the pages hereafter are appointed to professionals. We also remind you that our products description only aims to inform and are not intended for promotional purpose.





## INTRODUCTION



Oscar Boscarol srl has always encouraged the research for vanguard solutions and technologies with compelling quality that might match the operators needs during emergency activities. This is the force that led our products to become a helping hand, an indispensable support and constant reference for all insiders. Today the Oscar Boscarol srl is a worldwide leader exporting half of its products to foreign markets.





## MADE IN ITALY

---

Our OB medical suction units are entirely manufactured in Italy at our headquarters in Bolzano in order to guarantee care and quality control, starting from the materials used up to the assembly of the device and the control of its correct functioning before customer delivery.

## HISTORY-RELIABILITY

---

We have been manufacturing medical suction units since the 90s and we are at the top of the international competition with over 40,000 units sold worldwide (Europe, Asia, India, South America and Africa).

## COMPLIANCE WITH REGULATIONS

---

All our devices comply with current regulations: our quality system was developed according to the ISO 9001 and ISO 13485 standards, by always placing first the need to produce qualitatively impeccable devices that comply with the mandatory requirements.

## DIRECTIVES, RULES AND REGULATIONS OF REFERENCE

---

Current regulatory references our suction units comply with.

**MDD 93/42/CEE** Council Directive concerning medical devices

---

**ISO 10079-1** Medical suction equipment - Part 1: electrically powered suction equipment

---

**IEC 60601-1** Medical electrical equipment. General requirements for basic safety and essential performance

---

**IEC 60601-1-11** Medical electrical equipment general requirements for basic safety and essential performance - Collateral standard: requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment

---

**IEC 60601-1-12** Medical electrical equipment general requirements for basic safety and essential performance - Collateral standard: requirements for medical electrical equipment and medical electrical systems used in the emergency medical services environment

---

**EN 1789** Medical vehicles and their equipment - Road ambulances

---

**ECE R10 REV. 5** Uniform provisions concerning the approval of vehicles with regard to electromagnetic compatibility

---





OB 3000

## OB 3000

FLOW RATE-MAX VACUUM	33 LPM - 800 mbar (80 kPa)
RATED VOLTAGE	11÷25 Vdc
AUTONOMY	> 70 minutes
IP RATING	IP 44
BATTERY TYPE	removable LIPO 5 Ah 11,1 V
WARRANTY	5 years
PERIODIC MAINTENANCE	2 years (automatic reminder)
RTCA DO160G	YES - AVIO version
DIMENSIONS	360×244(h)×110 mm
WEIGHT	2,6 kg





## OB 3000

### POWERFUL AND RELIABLE

Max. flow rate 33 liter/minute with double piston suction pump and max. vacuum rate 800 mbar. Inner microprocessor records production data and functioning tracing and a safety inspection reminder indicates that the suction unit has to be sent for the scheduled maintenance. Analogic fluorescent vacuum meter, with mmHg scale upon request.



### BATTERY AND RECHARGING SYSTEM

Battery can be easily replaced by user and battery status indicator with no need to turn the device on.

Upon request the device can be supplied with a wall bracket (conformed to EN 1789) and with a 2-poles main power supply male connector.

### WEIGHT AND DIMENSIONS

FA version: 1000 ml autoclavable canister with protection overflow valve and protection filter directly on the cover.

FM version: 1000 ml OB-J canister suitable for the use of Serres liners. Extreme lightwightness and reduced dimensions make of this unit an easily portable device.



- 
- |          |   |                                      |
|----------|---|--------------------------------------|
| BSU3000  | - | OB 3000 FA MEDICAL SUCTION UNIT      |
| BSU3020  | - | OB 3000 FM MEDICAL SUCTION UNIT      |
| BSU3000A | - | OB 3000 AVIO FA MEDICAL SUCTION UNIT |
| BSU3020A | - | OB 3000 AVIO FM MEDICAL SUCTION UNIT |
-



## OB MINI

FLOW RATE-MAX VACUUM	25 LPM - 850 mbar (85 kPa)
RATED VOLTAGE	11÷25 Vdc
AUTONOMY	> 100 minutes
IP RATING	IP 43
BATTERY TYPE	LIPO 5 Ah 11,1 V
WARRANTY	3 years
PERIODIC MAINTENANCE	2 years (automatic reminder)
RTCA DO160G	YES - AVIO version
DIMENSIONS	290×240(h)×95 mm
WEIGHT	1,8 kg



## OB MINI

### RELIABLE AND ULTRA LIGHTWEIGHT

Max. flow rate 25 liter/minute with double piston suction pump and max. vacuum rate 850 mbar. Inner microprocessor records production data and functioning tracing and a safety inspection reminder indicates that the suction unit has to be sent for the scheduled maintenance. Only 1,8 kg for a device complete with 500cc canister: maximum lightweight and compactness.



### JAR

500 cc autoclavable canister or upon request 1 l FA and 1 l FM version.  
FA version: 1000 ml autoclavable canister with protection overflow valve and protection filter directly on the cover.  
FM version: 1000 ml OB-J canister suitable for the use of Serres liners.



### OPTIONAL BAG SIDE

2-poles main power supply male connector.

Blue transport bag.

Lateral poche with velcro available upon request.

BSU360EU - OB MINI 500ML MEDICAL SUCTION UNIT BLUE EU  
BSU360JP - OB MINI 500ML MEDICAL SUCTION UNIT BLUE JP/USA  
BSU360UK - OB MINI 500ML MEDICAL SUCTION UNIT BLUE UK  
BSU370 - OB MINI AVIO 500ML MEDICAL SUCTION UNIT BLUE EU  
BSU370JP - OB MINI AVIO 500ML MEDICAL SUCTION UNIT BLUE JP  
BSU370UK - OB MINI AVIO 500ML MEDICAL SUCTION UNIT BLUE UK  
BOR89471 - MINI TYRPOL BLUE BAG  
BSU364EU - OB MINI FA MEDICAL SUCTION UNIT BLUE EU  
BSU364JP - OB MINI FA MEDICAL SUCTION UNIT BLUE JP  
BSU364UK - OB MINI FA MEDICAL SUCTION UNIT BLUE UK  
BSU368EU - OB MINI FM MEDICAL SUCTION UNIT BLUE EU  
BSU368JP - OB MINI FM MEDICAL SUCTION UNIT BLUE JP  
BSU368UK - OB MINI FM MEDICAL SUCTION UNIT BLUE UK





OB 2012

## OB2012

FLOW RATE-MAX VACUUM	30 LPM - 800 mbar (80 kPa)
RATED VOLTAGE	12 Vdc
AUTONOMY	> 60 minutes
IP RATING	IP 34d
BATTERY TYPE	12 V 4 Ah Pb/Ca
WARRANTY	3 years
PERIODIC MAINTENANCE	2 years
RTCA D0160G	NO
DIMENSIONS	370×290(h)×130 mm
WEIGHT	4,2 kg





## HARD AND RESISTANT

Max. flow rate 30 liter/minute with double piston suction pump and max. vacuum rate 800 mbar. Autonomy >60 minutes.



## JAR

FA version: 1000 ml autoclavable canister with protection overflow valve and protection filter directly on the cover.  
FM version: 1000 ml OB-J canister suitable for the use of Serres liners.



## RECHARGING SYSTEM AND BAG SIDE

Upon request the device can be supplied with a wall bracket (conformed to EN 1789) and with a 2-poles main power supply male connector. Red transport bag. Lateral pouche with velcro available upon request.

---

BSU150	- OB2012 FA MEDICAL SUCTION UNIT
BSU154	- OB2012 FA MEDICAL SUCTION UNIT + WALL BRACKET
BSU158	- OB2012 FA MEDICAL SUCTION UNIT COMPLETE
BSU100	- OB2012 LINER MEDICAL SUCTION UNIT
BSU104	- OB2012 LINER MEDICAL SUCTION UNIT + WALL BRACKET
BSU108	- OB2012 LINER MEDICAL SUCTION UNIT COMPLETE
BOR89469	- MINI TYRPOL RED BAG

---

Codes content summary table on p.13.



OB 1000

## OB1000

FLOW RATE-MAX VACUUM	25 LPM - 850 mbar (85 kPa)
RATED VOLTAGE	12 Vdc
AUTONOMY	> 45 minutes
IP RATING	IP 34d
BATTERY TYPE	12 V 3.4 Ah Pb/Ca
WARRANTY	2 years
PERIODIC MAINTENANCE	2 years
RTCA DO160G	YES - AVIO version
DIMENSIONS	335×290(h)×105 mm
WEIGHT	3,3 kg



## RELIABLE AND AVIO VERSION

Max. flow rate 25 liter/minute with double piston suction pump and max. vacuum rate 850 mbar. Autonomy >45 minutes. Available in AVIO version: conformed to RTCA DO 160G.



## JAR

FA version: 1000 ml autoclavable canister with protection overflow valve and protection filter directly on the cover.  
FM version: 1000 ml OB-J canister suitable for the use of Serres liners.



## RECHARGING SYSTEM AND BAG SIDE

Upon request the device can be supplied with a wall bracket (conformed to EN 1789) and with a 2-poles main power supply male connector. Red transport bag.

Lateral pouche with velcro available upon request.

---

BSU210	- OB1000 FA MEDICAL SUCTION UNIT
BSU216	- OB1000 FA MEDICAL SUCTION UNIT EU
BSU220	- OB1000 FA MEDICAL SUCTION UNIT COMPLETE
BSU232	- OB1000 AVIO FA MEDICAL SUCTION UNIT
BSU212	- OB1000 FM MEDICAL SUCTION UNIT
BSU226	- OB1000 FM MEDICAL SUCTION UNIT EU
BSU228	- OB1000 FM MEDICAL SUCTION UNIT COMPLETE
BSU230	- OB1000 AVIO FM MEDICAL SUCTION UNIT
BOR89469	- MINI TYRPOL RED BAG

---

Codes content summary table on p.12.

## OB1000 CODES INDEX



**OB1000**



**12 VCC CABLE**

BSU855: 2 poles  
recharging supplying  
cable



**WALL BRACKET**

BSU806: wall bracket



**POWER SUPPLY**

BSU897EU: 2 poles  
power supply 13,7V

### BSU210

OB1000 FA MEDICAL  
SUCTION UNIT



-

-

### BSU216

OB1000 FA MEDICAL  
SUCTION UNIT EU



-



UK and JP/USA connections  
available upon request

### BSU220

OB1000 FA MEDICAL  
SUCTION UNIT COMPLETE



-

### BSU232

OB1000 AVIO FA MEDICAL  
SUCTION UNIT



-

-

### BSU212

OB1000 FM MEDICAL  
SUCTION UNIT



-

### BSU226

OB1000 FM MEDICAL  
SUCTION UNIT EU



-



UK and JP/USA connections  
available upon request

### BSU228

OB1000 FM MEDICAL  
SUCTION UNIT COMPLETE



-

### BSU230

OB1000 AVIO FM MEDICAL  
SUCTION UNIT



-

-

**IT IS POSSIBLE TO REQUEST MAIN POWER SUPPLY ADAPTERS WITH PLUGS FOR UK AND JP/USA SOCKETS**



## OB2012 CODES INDEX



**OB2012**



**12 VCC CABLE**

BSU855: 2 poles  
recharging supplying  
cable



**WALL BRACKET**

BSU800: wall bracket



**POWER SUPPLY**

BSU897EU: 2 poles  
power supply 13,7V

### BSU150

OB2012 FA MEDICAL  
SUCTION UNIT



-

-

### BSU154

OB2012 FA MEDICAL  
SUCTION UNIT + WALL  
BRACKET



-

### BSU158

OB2012 FA MEDICAL  
SUCTION UNIT COMPLETE



UK and JP/USA connections  
available upon request

### BSU100

OB2012 LINER MEDICAL  
SUCTION UNIT



-

-

### BSU104

OB2012 LINER MEDICAL  
SUCTION UNIT + WALL  
BRACKET



-

### BSU108

OB2012 LINER MEDICAL  
SUCTION UNIT COMPLETE



UK and JP/USA connections  
available upon request

**IT IS POSSIBLE TO REQUEST MAIN POWER SUPPLY ADAPTERS WITH PLUGS FOR UK AND JP/USA SOCKETS**

## OB500 - STATIONARY MEDICAL SUCTION UNIT

### JAR

### STATIONARY SUCTION UNIT FOR AMBULANCE

Max. flow rate 30 liter/minute and max. vacuum rate 800 mbar. Cable for the connection to the vehicle lighter (12Vcc).  
Dimensions: 9x8,5x6 cm  
Weight: 4,6 kg

FA version: 1000 ml autoclavable canister with protection overflow valve and protection filter directly on the cover.  
FM version: 1000 ml OB-J canister suitable for the use of Serres liners.



### MODULE OR EMBEDDING

Available in MODULE version with external regulation device and in EMBEDDING underwall version.

- 
- BSU464 - OB500 FM MEDICAL SUCTION UNIT UNDERWALL EMBEDDING
  - BSU442 - OB500 FA MEDICAL SUCTION UNIT UNDERWALL EMBEDDING
  - BSU462 - OB500 LINER MEDICAL SUCTION UNIT UNDERWALL EMBEDDING
  - BSU414 - OB500 FM MEDICAL SUCTION UNIT MODULE
  - BSU402 - OB500 FA MEDICAL SUCTION UNIT MODULE
  - BSU412 - OB500 LINER MEDICAL SUCTION UNIT MODULE
-

## MANUAL MEDICAL SUCTION UNIT

---



### E-VAC MANUAL SUCTION UNIT

---

Functional, safe, portable, powerful and economic suction unit, it solely works by hand pressure, allowing to achieve a depression up to 450 mbar.

Its "GUN" shape permits to use it with just one hand during emergency activities. Disinfectable and washable 250 ml jar for secretion collection.

Dimensions: about 20x6x19 cm

Weight: about 310 g

CE medical device.

---

BSU660 - E-VAC MANUAL SUCTION UNIT

BSU668 - RIGID PROBE FOR EVAC

BSU664 - JAR FOR E-VAC

---



### STERILE NEWBORN SUCTION UNIT

---

Sterile newborn suction unit.

Lightweight, economic and efficient.

Suitable inside newborn reanimation cases and sets. Disposable.

---

400100 - STERILE NEWBORN SUCTION UNIT

---

## FM AND LINER - JARS WITH DISPOSABLE SERRES LINERS



### OB-J CANISTER

New canister OB-J for secretion collection made of polypropylene for medical use. The jar features one canister and one adapter that allows the use of the jar with the disposable liner produced by the company Serres (non included). Complete with "L" joint. The canister is vapour autoclavable at a max temperature of 121°C and at a pressure of 2bar (200kPa). The jar can also be disinfected accordingly to local hospital procedures.  
Capacity: 1000 ml.  
CE medical device.

BSU506 - OB-J CANISTER



### SERRES CANISTER

Serres secretion jar for Serres disposable bags, made of polycarbonate and complete with safety overflow valve.  
Capacity: 1000 ml.  
CE medical device.

57308 - SERRES CANISTER



### SERRES DISPOSABLE BAGS

1000 ml disposable bag with lid and integrated antibacterial filter. Available in set of 6-12-36 disposable bags to be used with OB-J or Serres jar.  
CE medical device.

BSU705- 6 SERRES DISP. BAGS  
BSU706- 12 SERRES DISP. BAGS  
BSU707- 36 SERRES DISP. BAGS



## FA - AUTOCLAVABLE JARS



### 500 ML CANISTER

Suction jar for medical suction units, made of autoclavable and shockproof Macrolon with pressure closing lid made of silicone-rubber. Complete with overflow valve. Capacity: 500 ml.

BSU504 - 500ML CANISTER



### OB-J FA CANISTER

Autoclavable suction jar (121° C - 2 bar max.) made of polypropylene for medical use. Complete of overflow security valve. The protection filter and the angular joint connector for the aspiration tube patient side are directly connected to the cover. Capacity: 1000 ml. CE medical device.

BSU500 - OB-J FA CANISTER



### PROTECTION FILTER

Protection filter with filtrating support made of PTFE material and polypropylene hermetic container. Maximum suitable pressure up to 1 bar (100kPa). Available in sets of 5-15-40 filters. CE medical device.

BSU730 - 5 FILTER SETS  
BSU732 - 15 FILTER SETS  
BSU734 - 40 FILTER SETS

## RECHARGING AND SUPPLYING SYSTEM FOR SUCTION UNIT



### WALL BRACKET

Wall bracket anchorage for Boscarol suction units with integrated recharging system. Anchorage system tested up to 10g capacity. Anchorage system tested up to 10g capacity, obligatory requirement according to norm EN 1789. Extremely compact and functional wall bracket, not protruding out of rescue vehicle wall, with fast release button. Cable for the connection to the vehicle electric system (12V power supply).

BSU806 - WALL BRACKET FOR OB1000  
BSU800 - WALL BRACKET FOR OB2012 AND OB 3000



### MAIN POWER SUPPLY ADAPTER 2 POLES 13,7 V

Main power supply adapter with 2 poles male connector. Suitable for all BSU suction units. Wide range of input voltage: 100-240 Vac 50/60 Hz. Output voltage 13,7 Vdc. Max. nominal current output 5 A. The suction unit can be turned-on and fully operational during the charging activity. It is an accessory of medical device. Available in EU, UK and JP version.

BSU897EU - 13,7 V 2P POWER SUPPLY EU  
BSU897UK - 13,7 V 2P POWER SUPPLY UK  
BSU897JP - 13,7 V 2P POWER SUPPLY JP



### RECHARGING-SUPPLYING CABLE 2P 12VDC

2 poles recharging and power supply cable for Medical Suction Units models. Supplied with universal car lighter plug.

BSU855 - SUPPLYING CABLE 2P 12VCC

## ACCESSORIES



### ASPIRATION PROBE

Jankauer rigid aspiration probe, complete of 180 cm long connecting tube with female connection. Patient aspiration using only one hand. Rigid canula instantly ready and easy to use. Vacuum regulation by the thumb of the hand. Suitable for every kind of medical secretion suction unit. Singularly sterile packaged. CE medical device.


126140107191 - RIGID ASPIRATION PROBE COMPLETE  
126140105021 - SET AIRWAY STERILE BULB  
126140104021 - CONNECTING TUBE MT 1,8 FOR JANKAUER



### SUCTION UNIT CATHETER

Transparent suction catheter, with silicone. 2 lateral holes, disposable in sterile package with colored connector, available in different sizes.

11214101003 - CATHETER BLACK CH10  
11214101004 - CATHETER WHITE CH12  
11214101005 - CATHETER GREEN CH14  
11214101006 - CATHETER ORANGE CH16  
11214101007 - CATHETER RED CH18  
11214101008 - CATHETER YELLOW CH20



### CATHETER JOINT FINGER TYP

Suction-catheter joint "FINGER TYP" supplied in sterile package. Hole for vacuum regulation. Available in set. CE medical device.

BSU750 - FINGER-TYP 5 PCS SET  
BSU752 - FINGER-TYP 15 PCS SET  
BSU754 - FINGER-TYP 50 PCS SET



### SILICONE PATIENT TUBE

Autoclavable medical silicone patient tube  
Length: about 130 cm

BSU902 - SILICONE PATIENT TUBE 130CM



## PRACTICAL AND RESISTANT

---

Our rescue rucksacks are designed to guarantee practicality and resistance. We personally check the quality of our productions to provide even better products.



## EMERGENCY RUCKSACKS



### EMERGENCY RUCKSACK VIPER 2.0 - BOR87000N

Made of red nylon in tear-proof and waterproof fabric. Equipped with double handle (upper or lateral). Large front opening pocket with 2 internal elastic bands. 2 small side pockets and a large one on the outside. Equipped with front, side and rear reflective strips. The part in contact with the back of the user is equipped with cushions. Large front opening pocket with 2 internal elastic bands. 3 external perimeter pockets. 8 removable coloured pouches with transparent front and Velcro strips. Internal central partition with zipper. Reinforced bottom and rubber feet. Supplied empty. Dimensions: about 46x26x50 cm. Weight: about 3,8 kg.



### EMERGENCY RUCKSACK VIPER TYRPOL - BOR87012

Made of red tearproof fabric Tyrrpol material. Equipped with upper handle and shoulder straps for transport. Wide frontal total-opening pocket with 2 internal elastic strips. 2 perimeter external pockets. Lateral and frontal yellow refractive strips. Supplied with 6 coloured removable pouches with transparent frontal window and Velcro-strips. Central bulkhead with zip, complete with elastics. Reinforced bottom basin and rubber feet. Dimensions: 34,5x19,5x46,5 cm. Weight: about 4,2 kg. Supplied empty.



### EMERGENCY RUCKSACK MINI VIPER 2.0 - BOR87002N

Red tearproof and water-proof material. Equipped with double upper and side handle for easier transport. Large front opening pocket with 2 internal elastic bands for the fixing of material. 2 small side pockets and a large one on the outside. Equipped with front, side and rear reflective strips. The part in contact with the back of the user is equipped with cushions. 4 removable coloured pouches with transparent front and Velcro strips. Reinforced bottom and rubber feet. Supplied empty. Dimensions: about 44x26x44cm. Weight: about 2,3 kg. Maximum weight. capacity: 30kg

## EMERGENCY RUCKSACKS AND TROLLEY

---

### EMERGENCY RUCKSACK MASTERMED 2.0 - BOR87020N

---



Made of red nylon in tear-proof and waterproof fabric. The two front pockets are made with transparent internal compartments. The front pocket is equipped with elastic bands. The shoulder straps are padded and equipped with reflective strips along their entire length.

8 coloured inner bags with transparent front, 4 of which with handle. The part in contact with the back of the user is equipped with cushions.

Dimensions: about 35x30x53 cm.

Weight: about 4,3 kg.

Maximum weight capacity: 40 kg.

### EMERGENCY RUCKSACK MAXIMED - BOR87010

---



Made of red tearproof fabric material. Padded straps. 2 frontal equipped pockets. 1 perimeter external pocket. 1 perimeter external pocket with 4 coloured removable pouches. 4 inner removable coloured pouches with transparent frontal. Elastics alongside the internal perimeter of the rucksack. Oxygen cylinder holder straps placed centrally.

Dimensions: about 48x30x52(h) cm – 54 lts.

Weight: about 4,5 kg. Supplied empty.

### ZAINO MASTER TROLLEY TYRPOL - BOR87220

---



Trolley Rucksack Master Module empty, suitable for carrying materials in between hospital departments. Made of red tearproof impermeable Tyropol material. 2 refractive strips. Pullout trolley handle.

Pullout padded shoulder-straps. 2 frontal equipped pockets. Inner red, rigid, shockproof basin with central oxygen cylinder support. 4 removable coloured pouches with transparent frontal.

Dimensions: about 33x35x59(h) cm – 43 lts.

Weight: about 5 kg. Supplied empty.

## EMERGENCY BAGS

---



### BAG 33 - BOR87300

---

Made of red tearproof and impermeable fabric material. 1 capacious outer pocket  
5 inner adjustable spaces. Inside provided with elastics and a documents pocket.  
1 transparent frontal pocket for association logo.

Dimensions: about 36x18x35(h) cm.

Weight: about 1,6 kg. Supplied empty.



### MULTICASE - BOR80390

---

Impermeable tearproof material. To be carried by handle and shoulder strap. 3 transparent pockets inside inner top. Inner rigid basin.

Dimensions: about 48x26x16h cm.

Weight: about 1,2 kg. Supplied empty.




### PROFI BAG MINI - BOR89000

---

Emergency empty Profi Bag. Made of impermeable and tearproof red fabric material. With handle and shoulder straps. Modular interior. 2 wide lateral external pockets. Internal pocket with transparent.

Dimensions: about 36x26x20(h) cm.

Weight: 0,7 kg.



### PROFI BAG - BOR89010

---

Emergency empty Profi Bag. Made of impermeable and tearproof red fabric material. With handle and shoulder straps. Modular interior. 1 frontal pocket. Arranged to hold 1-2 lts oxygen cylinder. 2 wide lateral external pockets. Internal pocket with transparent.

Dimensions: about 60x25x27(h) cm.


Weight: 1,4 kg.



## TRAUMA BAG

### TRAUMA BAG 10 - BOR89040

Easy to carry either by hand, on shoulder or as a backpack.  
4 external lateral pockets, 4 inner modular spaces and 1 inner removable pouch with transparent frontal.  
Dimensions: about 45x33x30(h) cm.  
Weight: about 2,2 kg. Supplied empty.



### TRAUMA BAG 20 - BOR89050

Easy to carry either by hand, on shoulder or as a backpack.  
4 external lateral pockets, 5 inner modular spaces and 1 inner removable pouch with transparent frontal.  
Dimensions: about 53x33x34 cm.  
Weight: about 2,9 kg. Supplied empty.



### TRAUMA BAG 30 - BOR89060

Easy to carry either by hand, on shoulder or as a backpack.  
5 external lateral pockets, 5 inner modular spaces and 2 inner removable pouches with transparent frontal.  
Dimensions: about 63x37x31 cm.  
Weight: about 2,3 kg. Supplied empty.



BOR89040 - TRAUMA BAG 10 RED  
BOR89050 - TRAUMA BAG 20 RED  
BOR89060 - TRAUMA BAG 30 RED



## EMERGENCY RUCKSACKS ATLS


### ATLS O2 - BOR81145

Inner, rigid basin with central oxygen bottle holder. 4 coloured removable pouches with transparent frontal.  
Dimensions of all 3 models:  
about 32x25x53(h) cm.  
Volume: 35 lts. Weight: about 3,3 kg.



### ATLS MODULE - BOR81146

Inner, rigid basin with 4 compartments and 4 pouches.  
4 coloured removable pouches with transparent frontal.



### ATLS WK - BOR81147

Inner, rigid basin with 4 compartments.  
1 lateral oxygen bottle holder.  
Equipped with removable partitioning modules, Velcro type.  
4 coloured removable pouches with transparent frontal.



BOR81145 - EMERGENCY RUCKSACK ATLS O2  
BOR81146 - EMERGENCY RUCKSACK ATLS MODULE  
BOR81147 - EMERGENCY RUCKSACK ATLS WK

## RESCUE BAGS

---



### MINI RESCUE BAG - BOR80290

---

Compact and practical emergency bag  
Hand and shoulder transportable.  
Made in tearproof and waterproof material. 2 windows for the insertion of identity cards/logos. Lower part with inner shockproof basin. Upper part with elastics.

Size: about 39 x 18 x 25h cm.

Weight: about 1 kg .




### RESCUE BAG MODULE - BOR80340

---

Compact and practical empty rescue bag  
Impermeable tearproof material. Easy to carry either by hand, on shoulder or as a backpack. 2 frontal transparent pockets for association logo. Bottom side with shockproof inner basin. 8 removable pouches with transparent. Upper part with inner elastics suitable for equipment fixing.

Dimensions: about 47x25x27(h) cm

Weight: about 2,4 kg.




### BABY CHILDREN REANIBAG - BOR89090

---

Compact but extremely large rescue bag.  
Made of tear-proof and waterproof nylon.  
The fabric has bright colours and textures suitable for use in the paediatric sector.  
Equipped with 4 reflective strips placed on the sides of the bag. Bag complete with handles, belt and shoulder straps for comfortable and safe transport. The 2 external side pockets have internal spaces and elastic bands. The inside of the bag is divided into 4 separate spaces that can be removed. In the upper part there is a removable bag with a transparent front.  
Dimensions: about 44x30x29 cm.  
Weight: about 2,1 kg. Supplied empty.

## DRUGS BAGS



### DRUGS BAG 14 PHIALS - BOR89336

Small red drugs bag suitable to be fixed on a belt. Made of tear-resistant and water-repellent material. Comfortably worn on a belt with a heavy-duty belt loop. 2 parts with adjustable elastics to hold 14 different sized phials securely in place. Central padding adds extra protection for fragile vials. Dimensions: about 13x12x6 cm. CE medical device.




### DRUGS BAG 28 PHIALS - BOR89330

Book-opening drugs case empty. Made in tearproof impermeable red fabric material. 2 rear Velcro-strips for the fixing in bags and emergency rucksacks. 2 removable plates with elastics suitable for 28 different sized phials. Central padding for ampoules record. Dimensions: about 25x12x7(h) cm. CE medical device.



### DRUGS BAG 60 PHIALS - BOR89332

Big book-opening drugs bag. 2 parts with adjustable elastics to feature 60 different sized phials. Central padding with transparent pocket for ampoules list and expiry date. Made in red, tear-resistant, water-resistant material. 2 external Velcro-strips for the fixing in bags or rucksacks. Dimensions: about 25x25x7h cm. CE medical device.



### DRUGS BAG EASY RAIL + ELASTICS - 60 PHIALS BOR89340

Book-opening drugs case with EASY RAIL system. Made of tearproof impermeable Tyrol material. EASY RAIL system allows a quick insertion and removal of phials. Suitable for 60 different sized phials in addition to little vials (48 easy rail rigid spaces and 12 elastics). Internal padded divider. Transparent frontal window. Dimensions: about 29x25x8 cm. CE medical device.

## OXYGEN BOTTLE HOLDERS

---



### **PADDED OXYGEN CYLINDER HOLDER 2 L - BOR80613**

---

Padded oxygen cylinder bag (up to 2 lts).  
Tearproof impermeable material.  
Equipped with handle and shoulder strap.  
Equipped with compartments, pockets and elastics for items fixing.  
Dimensions: about 51x22 cm.  
Weight: about 0,7 kg. Supplied empty.



### **PADDED OXYGEN CYLINDER HOLDER 5 L - BOR80614**

---

Padded oxygen cylinder bag (up to 5 litres).  
Tearproof impermeable material . External pockets for accessories. Equipped with handle and shoulder strap.  
Dimensions: about 75x35x17 cm.  
Weight: about 0,5 kg.



### **OXYGEN CYLINDER BAG 2 L**

---

Oxygen cylinder bag (up to 2 litres). Made in nylon impermeable tearproof material.  
Dimensions red cylinder bag: Ø 13x44 cm.  
Dimensions blue cylinder bag: Ø 10x36 cm.  
Weight: about 200 g. Supplied empty.

---

BOR80615 - BLUE OXYGEN CYLINDER BAG  
BOR80617 - RED OXYGEN CYLINDER BAG

---



## POUCHES

---



### BLUE DOUBLE POUCHES - BOR81042

---

Dimensions: about 25x25x8 cm.



### STANDARD BIG POUCHES

---

Dimensions: about 31x15x8 cm

---

BOR89452 - STD BIG POCHE RED  
BOR89454 - STD BIG POCHE GREEN  
BOR89456 - STD BIG POCHE YELLOW

---



### STANDARD POUCHES

---

Dimensions: about 27x12x9 cm

---

BOR81024 - STD POUCH PINK  
BOR81028 - STD POUCH RED  
BOR81030 - STD POUCH GREEN  
BOR81032 - STD POUCH YELLOW  
BOR81034 - STD POUCH BLUE  
BOR81036 - STD POUCH VIOLET

---

## POUCHES

---



### BIG POUCHES WITH HANDLE

---

Dimensions: about 30X12X13 cm circa

BOR89440 - GREEN BIG POUCH WITH HANDLE  
BOR89442 - RED BIG POUCH WITH HANDLE  
BOR89444 - BLUE BIG POUCH WITH HANDLE  
BOR89446 - YELLOW BIG POUCH WITH HANDLE

---



### POUCHES WITH HANDLE

---

Dimensions: about 23x14x11 cm

BOR89420 - YELLOW POUCH WITH HANDLE  
BOR89422 - BLUE POUCH WITH HANDLE  
BOR89424 - GREEN POUCH WITH HANDLE  
BOR89426 - RED POUCH WITH HANDLE

---



### VIPER POUCHES

---

Dimensions: about 19x12x6 cm  
Velcro-strips: female

BOR81058- YELLOW VIPER POUCH  
BOR81060- GREEN VIPER POUCH  
BOR81062- RED VIPER POUCH  
BOR81064- BLUE VIPER POUCH  
BOR81066- VIOLET VIPER POUCH  
BOR81068- BLACK VIPER POUCH  
BOR81070- ORANGE VIPER POUCH

---



### MINI VIPER POUCHES

---

Dimensions: about 19x12x6 cm  
Velcro-strips: male

BOR82000 - VIOLET MINI VIPER POUCH  
BOR82001 - ORANGE MINI VIPER POUCH  
BOR82002 - BLACK MINI VIPER POUCH  
BOR82006 - DARK BLUE MINI VIPER POUCH

---

## POUCHES

---



### MINI POUCHES

---

Dimensions: about 15x12x8,5 cm

---

BOR89458 - BLACK MINI POUCH  
BOR89460 - RED MINI POUCH  
BOR89462 - BLUE MINI POUCH  
BOR89464 - GREEN MINI POUCH

---



### RED PVC MINI POCHE - BOR89469

---

Dimensions: about 10x7x13 cm

Designed to be attached to the OB2012, OB1000 AND OB MINI Boscrol Suction Unit transport bags.



### FANTASY POUCH

---

Dimensions small pouch: about 28x12x9 cm

Dimensions double pouch: about 25x25x8 cm

Dimensions medium: about 31x15x7 cm

---

BOR81050 - SMALL FANTASY POUCH  
BOR81052 - MEDIUM FANTASY POUCH  
BOR81054 - DOUBLE FANTASY POUCH WITH ELASTICS

---

## POUCHES



### OLIVE GREEN POUCHES

Dimensions small pouch:

about 24x15x11,5 cm

Dimensions big pouch:

about 34,5x12x14 cm

---

MIL1370 - OLIVE GREEN BIG POUCH WITH  
HANDLE

MIL1372 - OLIVE GREEN SMALL POUCH WITH  
HANDLE

---



### HOLSTER PROFESSIONAL

Particularly functional belt case. Made in tearproof impermeable fabric material. Suitable to contain one sphygmomanometer. 3 inner spaces with closing.

Dimensions: about 34x11x8(h) cm.

Weight: about 100 g. Supplied empty.

---

BOR206 - BLUE HOLSTER PROFESSIONAL

---



### FIRST AID - BOR220

Made in tearproof impermeable red material. Book opening. Complete with 3 net pockets for emergency equipment. Waist wearing mode.

Dimensions: about 13x6x18(h) cm.

Weight: about 300 g. Supplied empty.



### FABRIC BAG - BOR84072

Fabric empty bag. Tearproof impermeable blue material. With shoulder strap.

Dimensions: about 29x16x16 cm.

Weight: about 210 g. Supplied empty.



## E-HOLSTER

---



### EMPTY E-HOLSTER

---

E-Holster belt made in tearproof and water-resistant red material. 5 compartments with pockets for gloves and for Pen Light. Belt loop of about 5,5 cm. Dimensions: about 14x3x19(h) cm. Weight: about 100 g. Supplied empty.

---

BOR84062 - BLUE EMPTY E-HOLSTER

BOR84063 - RED EMPTY E-HOLSTER

---



### COMPLETE E-HOLSTER

---

E-Holster belt made in tearproof and water-resistant red material. 5 compartments with pockets for gloves and for Pen Light. Belt loop of about 5,5 cm. Complete with: 1 first aid scissors, 1 anatomic forceps, Cocher forceps, 1 Pen Light, 1 Tourniquet. Dimensions: about 14x3x19(h) cm. Weight: about 200 g.

---

BOR84064 - BLUE COMPLETE E-HOLSTER

BOR84065 - RED BLUE COMPLETE E-HOLSTER

---

## CABINETS

---



### VISIO DEFIBRILLATOR METAL CABINET - DEF040

---

Metal Cabinet for Defibrillator with see through door. Made of Galvanized Carbon Steel PO2. Enameled by Epoxy Polyester white powder RAL9010 Heated at 180°. Transparent door with opening at 180° with a galvanized strong pivot. Strong Nickel plated zinc alloy closing device. Can be sealed or locked. Back holes for wall hanging. Dimensions: 425x425x160 mm. Weight: 6 kg. Supplied empty.



### CABINET - CAV605

---

Carbon steel cabinet, coated in epoxy polyester. Rounded edges. Three steel shelves. Door with 3 accessories holder. Locker with key. Dimensions: 402x202x483(h) mm. Weight: 6,60 kg. Supplied empty.

## ALUMINIUM CASES

---



### ALUMINIUM CASE TYPE A - VAL42000

---

Type A emergency medical case made of anodized aluminium. High resistance and long durability. Inside made in washable and disinfectable thermo-format with modular system. Arranged to accommodate oxygen cylinders (1 or 2 litres).

Dimensions: about 52x36x17(h) cm.

Weight: about 5,7 kg. Supplied empty.

---

VAL42000 - ALUMINIUM CASE TYPE A  
VAL42010 - 2 AMPOULE HOLDER RAIL  
DEPOT FOR CASE TYPE A

---



### ALUMINIUM CASE TYPE B - VAL43000

---

Type B emergency medical case, made of anodized aluminium. High resistance and long durability. Inside made in washable and disinfectable thermo-format with modular system. Arranged to accommodate oxygen cylinders (1 or 2 litres). Removable internal basin made in thermo-format with rail ampoules holder (max 60 pcs).

Dimensions: about 52x36x24(h) cm.

Weight: about 7,6 kg. Supplied empty.

---

VAL43000 - ALUMINIUM CASE TYPE B

---



### RAIL AND EASY RAIL VIALS HOLDER

---

We have available RAIL and EASY-RAIL vials holders made of painted metal with internal rubber, suitable for vials of any size. Available sizes indicated on request.

## CASES IN ABS

---



### SAN CASE IN ABS - VAL40130

---

San type case empty. Made in orange shockproof ABS material. Wall bracket anchorage system. 2 equal sides made of 2 adjustable inner compartments. 2 transparent dividers. Dimensions: about 32,5x12,5x21,5(h) cm. Weight: about 1,1 kg.



### MULTI CASE IN ABS - VAL40110

---

Multi type case empty. Made in orange shockproof ABS material. Wall bracket anchorage system. 2 equal sides made of 2 adjustable inner compartments. 2 transparent dividers. Dimensions: about 42x15x30(h) cm. Weight: about 3,3 kg.



### MEDIC CASE IN ABS

---

Medic 2 case empty with orange or transparent cover. Made in shockproof orange ABS material. Equipped with wall bracketing. Dimensions: about 39,5x13,5x26(h) cm. Weight: about 1,6 kg.

---

CAV399 - MEDIC 2 CASE TRANSPARENT COVER  
CAV401 - MEDIC 2 CASE ORANGE COVER

---



### FLEXI CASE IN ABS - CAV406


---

Flexi case empty. Made in shockproof orange polypropylene material. Equipped with wall bracketing. Dimensions: about 25,5x9x19,5 (h) cm. Weight: about 0,45 kg.



## WATER PROOF CASES

---



---

VAL43140 - SMALL BOSCAROL HIGH IMPACT CASE  
VAL43141 - MEDIUM BOSCAROL HIGH IMPACT CASE  
VAL43142 - LARGE BOSCAROL HIGH IMPACT CASE

---

### BOSCAROL HIGH IMPACT CASE

---

Boscarol High Impact small orange case with high resistant thick walls. Automatic pressure release valve. Easy to open double throw latches. Load tested soft grip handle. Holes for double padlocking. Holes for shoulder strap. Stackable interlocking system. Thick body material. Extra fixing points. Fully length hinges. Nylon pins, corrosion free. Pick N Pluck foams in the bottom and convoluted foam available on demand. Screen printing space on the lid. Rectangular resin label space on the lid.

Available in 3 sizes:

Small: about 33,6x30x14,8 cm.

Weight: 1,75 kg.

Medium: 46,4x36,6x17,6 cm.

Weight: 2,78 kg.

Large: 55,5x42,8x21,1 cm.

Weight: 3,86 kg.



### EXPLORER CASE

---

Emergency waterproof and watertight case up to 10 meters deep. With pressurization valve. Sturdy closing system with predisposition for a lead lock. Available in different sizes.

EXPLORER MINI

Dim.: 30,5x19,5x27h cm - Weight: 1,9 kg

EXPLORER P

Dim.: 41x20,5x34 cm - Weight: 3,1 kg

EXPLORER M

Dim.: 52x23x43,5 cm - Weight: 4,6 kg

EXPLORER

Dim.: 44,5x34,5x19 cm - Weight: 4 kg

EXPLORER G

Dim.: 65x24x51 cm - Weight: 7,1 kg

---

VAL43151 - EXPLORER MINI  
VAL43152 - EXPLORER P  
VAL43162 - EXPLORER M  
VAL43164 - EXPLORER  
VAL43172 - EXPLORER G  
VAL43174 - EXPLORER P WITH 4 POUCHES  
VAL43176 - EXPLORER M WITH 5 POUCHES  
VAL43178 - EXPLORER G WITH 6 POUCHES

---



## PRODUCTS FOR HEALTH EMERGENCY

---

Our immobilization and transport medical devices meet the standards in force when expected and were born from the desire to provide products dedicated to the specific sector of health emergency.

## MADE IN ITALY

---

The MADE IN ITALY manufacturing allows us to offer high quality medical devices, by verifying and tracing every phase of the production process.

## SUCTION VALVE

---



All our vacuum medical devices are equipped with a suction valve specifically designed and patented by Oscar Boscarol srl.

## VACUUM MATTRESSES

---



### SDM - IMM120012

---

Vacuum mattress SDM with 3 independent chambers. Perfect and absolute immobilization of traumatized patient. High resistance cloth. External and compact aspiration valve with ¼ on/off system. 4 transport handles on each side. 3 safety belts. 100% X-ray transparent. CAT scan compatible.

Dimensions: about 210x87 cm.

Weight: about 8,4 kg.



### PROFI - IMM120150

---

PROFI vacuum mattress with 3 independent chambers. Perfect and absolute immobilization of traumatized patient. High resistance material. Protective cloth under the mattress. On/off external, compact ¼ suction valve. 4 handles on each side obtained from the protective cloth itself. 2 head-handles and 2 foot-handles. 3 safety belts. 100% X-ray transparent. CAT scan compatible.

Dimensions: about 210x90 cm.

Weight: 10,5 kg.



### PROFI EVOLUTION - IMM120129

---

PROFI EVOLUTION vacuum mattress with 3 independent chambers. Provided with 2 loops designed for the mattress fixation onto the stretcher by using the stretcher fastening belts. Perfect and absolute immobilization of traumatized patient. High resistance material. Protective cloth under the mattress. On/off external, compact ¼ suction valve. 4 handles on each side obtained from the protective cloth itself. 2 head-handles and 2 foot-handles. 3 safety belts. 100% X-ray transparent. CAT scan compatible.

Dimensions: about 210x90 cm

Weight: 11 Kg



## VACUUM MATTRESSES

---



### DOLOMITES - IMM120154

---

Dolomites vacuum mattress with 3 independent chambers. Perfect and absolute immobilization of traumatized patient. High resistance material. Protection removable made in washable Cordura. On/off external, compact  $\frac{1}{4}$  suction valve. 4 handles on each side obtained from the protective cloth itself. 4 coloured cross safety foldaway belts for shoulder/chest and legs. 2 head straps. 3 standard safety belts. 100% X-ray transparent. CAT scan compatible. Dimensions: about 210x90 cm. Weight: about 10 kg.



### GENESIS - IMM120160

---

Genesis vacuum mattress with 4 independent chambers. Vacuum Head-Immobilisator integrated in the mattress itself with 2 chin-supports. Perfect and absolute immobilization of traumatized patient. The structure allows security belt passing through patient shoulders. High-resistance material. Cordura protection cover. On/off external, compact  $\frac{1}{4}$  suction valve. 3 safety belts. 2 additional coloured cross safety foldaway belts for shoulder/chest and legs in the case itself. 100% X-ray transparent. CAT scan compatible. Dimensions: about 210x97 cm. Weight: 10 kg.



### VACUUM MATTRESS TRANSPORT BAG - IMM120030

---

Transport bag for vacuum mattress. Tearproof impermeable material. Transport handles. Supplied with all the complete sets of mattresses, it can be used with all the vacuum mattresses we manufacture. Dimensions: about 89x24x72(h) cm.



## VACUUM MATTRESS AND SPLINT

---



### BABY - IMM120060

---

Vacuum Mattress BABY, for safe rescue and transport of newborn - pediatric patients. Complete body immobilization. High resistant material. CAT compatible and X-Ray transparent. Ventil aspiration valve. 4 transport handles. Dimensions: about 100x64 cm.



### VACUUM SPLINT

---

High resistant and washable material. Quick and painless immobilization of all kind of fractures (even exposed ones). Immobilizes without constrictions thus allowing blood circulation. External aspiration valve. 100% X-ray transparent (any metal parts included). Suitable for use on baby patients. CE medical device.

Dimension hand vacuum splint: 50,5x33 cm.  
Dimension arm vacuum splint: 75x47 cm.  
Dimension leg vacuum splint: 92x64 cm.  
Dimension leg 108 vacuum splint: 108x64 cm. Singularly or in set packaged, complete with transport bag or transport bag and suction pump.

- 
- IMM120172 - VACUUM SPLINT LEG 108CM
  - IMM120174 - VACUUM SPLINT HAND
  - IMM120176 - VACUUM SPLINT ARM
  - IMM120178 - VACUUM SPLINT LEG
  - IMM120200 - SET 2 VACUUM SPLINT WITH TRANSPORT BAG
  - IMM120180 - SET 2 VACUUM SPLINT COMPLETE (TRANSPORT BAG AND SUCTION PUMP)
  - IMM120190 - SET 3 VACUUM SPLINT WITH TRANSPORT BAG
  - IMM120184 - SET 3 VACUUM SPLINT COMPLETE (TRANSPORT BAG AND SUCTION PUMP)
  - IMM120218 - TRANSPORT BAG FOR VACUUM SPLINT
-

## ACCESSORIES FOR VACUUM DEVICES

---



### MAXI ALUMINIUM SUCTION PUMP - IMM120240

---

Suction pump for vacuum mattress and splints  
Made in aluminium with foot support for easy aspiration.

Dimensions: about 13x5,5x54(h) cm.

Weight: about 1,4 kg.



### MINI ALUMINIUM SUCTION PUMP - IMM120260

---

Suction pump for vacuum mattress and splints. Made in aluminium.

Dimensions: about 4x32 cm.

Weight: about 500 g.



### MINI HIGH EFFICIENCY SUCTION PUMP - IMM120262

---

Mini high efficiency aspiration pump.  
Powerful, practical, lightweight and small.  
Suitable for vacuum splints and for vacuum mattresses alike. Made in plastic with aluminium cylinder. Complete with pvc tube (about 1 mt). Dimensions: about 6,5x33 cm.  
Weight: about 400 g.



### REPAIR KIT - AIM5000

---

It contains five pieces of pvc fabric and a bottle of special glue: Make sure the damaged spot is not covered with oil or grease, clean it thoroughly and make it rough. Apply the special glue and let it dry (check with the back of the fingers). Apply the fabric with firm pressure to that area. Check the seal.

## RIGID VACUUM SPLINT AND SSPLINT

---



### RIGID VACUUM SPLINT

---

Rigid splint for upper/inferior limbs immobilization. Made in PVC and Nylon. Washable and X-ray transparent. Velcro closing.

Dimensions upper limbs:

about 51x27,5x3,5 cm.

Dimensions inferior limbs: 74x40x3,5 cm.

Available in set of 2 vacuum splints complete with transport bag.

---

IMM123100 - UPPER LIMBS RIGID SPLINT  
IMM123110 - INFERIOR LIMBS RIGID SPLINT  
IMMM123120 - SET 2 RIGID SPLINT WITH TRANSPORT BAG  
IMM123125 - TRANSPORT BAG FOR RIGID SPLINT

---



### SSPLINT

---

Splint system set with aluminium core for easy limbs immobilization in case of fracture or sprain. Made in neoprene material. Adaptable to fit each body frame. It does not cause tissue compressions. It maintains body temperature. It reduces congestion around injured areas. Velcro modelling closing. Lightweight and easily washable even in washing-machine. Least encumbrance.

Dim. leg ssplint: 58x21x0,6 cm

(not suitable for femur fractures).

Dim. arm ssplint: 60,5x17x0,6 cm.

Dim. forearm ssplint: 39x16x0,6 cm.

Dim. wrist ssplint: 32x15x0,6 cm.

Dim. ankle ssplint: 57x18x0,6 cm.

Available in set of 5 vacuum splints complete with transport bag.

---

IMM124030 - SSPLINT OB SET 5 PCS WITH TRANSPORT BAG  
IMM124031 - SSPLINT OB LEG  
IMM124032 - SSPLINT OB ARM  
IMM124033 - SSPLINT OB FOREARM  
IMM124034 - SSPLINT OB WRIST  
IMM124035 - SSPLINT OB ANKLE  
IMM123153 - RED TRANSPORT BAG FOR SSPLINT  
IMM123154 - BLUE TRANSPORT BAG FOR SSPLINT

---

## OB NEW ONE-PIECE CERVICAL COLLARS

One-piece cervical collar, made with injection mould system. Perfect immobilization of cervical section. Frontal opening for trachea check. Made of polypropylene X-ray transparent with non absorbent padding. CE medical device. Available in different sizes and in set of 2,3, 4 or 6 collars complete with transport bag.




**0 INFANT PINK - IMM124100**

Size: 0 - Infant pink



**1 PEDIAT. LIGHT BLUE - IMM124102**

Size: 1 - Pediatric light blue



**2 NECKLESS VIOLET - IMM124104**

Size: 2 - Neckless violet



**3 SHORT BLUE - IMM124106**

Size: 3 - Short blue



**4 REGULAR ORANGE - IMM124108**

Size: 4 - Regular orange



**5 TALL GREEN - IMM124110**

Size: 5 - Tall green



**SET OF 3 OR 4 COLLARS**

IMM124112: set of 3 OB NEW one-piece collars (size from 1 to 3) with transport bag.  
IMM124114: set of 4 OB NEW one-piece (size from 1 to 4) with transport bag.



**SET OF 6 COLLARS - IMM124116**

Set of 6 OB NEW one-piece collars (sizes from 0 to 5), with transport bag.  
IMM120656: transport bag.



**SET OF 2 COLLARS - IMM124118**

Set di 6 OB NEW one-piece collars (sizes 0 and 1), with transport bag.



## OBNECK TWO-PIECES CERVICAL COLLARS

Two-pieces cervical collar, made with injection mould system. Perfect immobilization of cervical section. Frontal opening for trachea check. Made of polypropylene X-ray transparent with non absorbent padding. CE medical device. Available in different sizes and in set of 2,3 or 6 collars complete with transport bag.



**0 INFANT PINK - IMM120782**

Size: 0 - Infant pink



**1 PEDIAT. LIGHT BLUE - IMM120784**

Size: 1 - Pediatric light blue



**2 NECKLESS VIOLET - IMM120786**

Size: 2 - Neckless violet



**3 SHORT BLUE - IMM120788**

Size: 3 – Short blue



**4 REGULAR ORANGE - IMM120790**

Size: 4 - Regular orange



**5 TALL GREEN - IMM120792**

Size: 5 - Tall green



**SET 3 COLLARS - IMM120794**

Set of 6 OBNECK two-pieces collars (sizes from 1 to 3), complete with transport bag.



**SET 6 COLLARS - IMM120796**

Set of 6 OBNECK two-pieces collars (sizes from 0 to 5), complete with transport bag.



**SET 2 COLLARS - IMM120798**

Set of 6 OBNECK two-pieces collars (sizes 0 and 1), complete with transport bag.

IMM120653 - TWO-PIECES COLLARS TRANSPORT BAG  
IMM120658 - FANTASY TWO-PIECES COLLARS TRANSPORT BAG

## ADJUSTABLE CERVICAL COLLARS

---

### CHAMPION ADJUSTABLE CERVICAL COLLAR

---

Double valve collar adjustable in different heights. Central hole for the control of the trachea. Flat two-piece design allows efficient placement in limited spaces. The end of the Velcro is in orange color for a better identification of the grip. Polyethylene structure, black LD padding, nylon buttons and Velcro strap latch. Can be worn during CAT scan and MRI. 100% Latex-Free. The device is disposable. However, it is possible to wash in water at 30°C with neutral soap.

---

9030 - ADULT CHAMPION COLLAR  
9031 - PEDIATRIC CHAMPION COLLAR

---

### FAST ADJUSTABLE CERVICAL COLLAR

---

Cervical collar adjustable to different heights. Central hole for the control of the trachea and the execution of a tracheotomy procedure.

Flat one-piece design allows efficient placement in limited spaces. Polyethylene structure, black LD padding, nylon buttons and Velcro strap latch. Latex-Free. The device is disposable. However, it is possible to wash in water at 30°C with neutral soap.

---

9015 - ADULT FAST COLLAR  
9015PED - PEDIATRIC ADULT FAST COLLAR

---

### PATRIOT ADJUSTABLE CERVICAL COLLAR

---

4 positions adjustable cervical collar Patriot: low, short, regular, tall. Simple and versatile usage. Excellent immobilization. Wide frontal opening allowing an adequate access to windpipe. Easily washable. X-Ray transparent.

---


PHP-PA - PATRIOT ADULT  
PHP-PP - PATRIOT PEDIATRIC  
PHP-PA-G - PATRIOT ADULT - OLIVE GREEN  
IMM120938 - SET 2 PATRIOT ADULT + PEDIAT.

---



## ADJUSTABLE CERVICAL COLLARS

---



### STIFNECK ADJUSTABLE CERVICAL COLLAR

---

4 positions adjustable cervical collar Stifneck: low, short, regular, tall. Easy and versatile. Perfect immobilization. Frontal opening for windpipe check. Easy to wash. X-ray transparent.

---

IMM120648 - STIFNECK ADULT COLLAR  
IMM120650 - STIFNECK PEDIATRIC COLLAR

---



### ECO CERVICAL COLLARS

---

Economic cervical collars. Realised in high density Elastolen/Veolene, non-deformable (closed cells polyethylene) resistant to external agents, totally non-toxic. Velcro closing. X-ray transparent.

---

IMM120770 - ECO COLLAR SMALL N. 0  
IMM120780 - ECO COLLAR TIGHT N. 1  
IMM120774 - ECO COLLAR MEDIUM N. 2  
IMM120778 - ECO COLLAR LARGE N. 3  
IMM120770 - SET 3 ECO COLLARS (0-2-3)

---



### SOFT CERVICAL COLLARS

---

Soft collars for cervical immobilization. Structure in high density foam rubber with tubular cotton jersey covering. Velcro closing.

---

IMM120830 - SOFT CERVICAL COLLARS - S  
IMM120832 - SOFT CERVICAL COLLARS - M  
IMM120834 - SOFT CERVICAL COLLARS - L  
IMM120836 - SOFT CERVICAL COLLARS - XL

---



### RESCUE CAN - 6400300

---

Light and imperishable. Printed in high density polyethylene without joints. Endowed with floating 186 cm rope and shoulder strap.

## IMMOBILIZATION DEVICES

---



### GREEN EED SPINAL IMMOBILIZER - IMM121964

---

Device for patients' immobilization and extraction. Immobilizes optimally the patient's torso, head and spinal column. Made of durable synthetic fabric and equipped with rigid bamboo internal slats. The immobilizer is provided with coloured fixing belts equipped with quick hook / unhook buckles. Completely x-Ray transparent. It can be washed with water-based detergents at room temperature. It is equipped with a transport case that can also contain cervical collars. CE medical device. Weight: about 3,3 kg.

---

IMM121964 - GREEN EED

IMM121664 - SET 2 CHIN-STRAPS FOR EED

---



### PATIENT TRANSFER SHEET - IMM120321

---

New transfer sheet, even more resistant. Made of red tearproof impermeable PVC material. 8 side handles. Easily washable and disinfectable. In accordance to EN1865 norm. Perimeter reinforced with a 4cm, double sewn band on both sides. In accordance to EN597-1 on devices flammability.

Dimensions: about 187x70 cm.

Weight: about 1200 g.

Loading: about 150 kg.



## SPINAL BOARDS



### OB SPINAL BOARD

High resistant plastic spinal board, suitable on difficult emergency activities. 4 buttonholes for quick and total fixing of head immobilizer and 2 hollows placed on leaning side for a traditional base blocking thus allowing strap belts not to get teared during use. Provided with 20 transport handles. X-ray translucent.

Maximum load: about 150 kg.

Dimensions: about 184x41x5,5 cm .

Weight: about 5,5 kg.

---

IMM121074 - OB SPINAL BOARD  
IMM121076 - OB SPINAL BOARD WITH SPIDER BELTS AND OB HEAD IMMOBILIZER  
IMM121075 - OB SPINAL BOARD WITH SPIDER BELTS AND WIND HEAD IMMOBILIZER

---



### PEDIATRIC SPINAL BOARD

Paediatric Spinal Board for children's immobilization and transport up to 30 kg in weight and 120 cm in height. Complete belt system for shoulder, chest, abdomen and legs, adjustable length allowing children's respiratory and circulatory system control. Inferior belt system allows independent legs and ankles treatment. Made in treated multi layer with blue and yellow Tyrol coating allowing radiology check. X-ray transparent. Expanded polyurethanes layer allowing patient comfort considerable improvement. Complete of 4 lifting and transport handles. Complete with transport bag. Dimensions: about 126x34x2(h) cm. Weight: about 5,1 kg.

Head immobilizer to be purchased separately.

---

IMM121046 - PEDIATRIC TYROPOL SPINAL BOARD  
IMM121048 - PEDIATRIC SPINAL BOARD WITH HEAD IMMOBILIZER  
IMM121668 - BABY HEAD IMMOBILIZER BLOCKS  
IMM121044/C - TRANSPORT BAG FOR SPINAL BOARD

---

## UNIVERSAL HEAD IMMOBILIZER

---



### OB HEAD IMMOBILIZER - IMM121628

---

Universal head immobilizer OB. Made in high density treated yellow foam material, water repellent. Featuring 1 basic pillow, two head support pillows with opening for ears access, 2 chin supports with central padding and 2 straps with Velcro already integrated to the basis for the positioning of the head Immobilizer on the stretcher. X-Ray translucent.

Dimensions: about 42,5x26,5x16(h) cm.

Weight: about 900 g.

---

IMM121628 - OB HEAD IMMOBILIZER

AIM4190 - BASE FOR HOB EAD IMMOBILIZER

AIM4191 - SET OF 2 CHIN STRAPS FOR OB HEAD IMMOBILIZER

AIM4192 - SET OF 2 HEAD-BLOCK FOR OB HEAD IMMOBILIZER

---



### WIND HEAD IMMOBILIZER - IMM121640

---

Universal head immobilizer WIND.

Realized in stuffed MDF and covered with washable, disinfectable, water repellent and thermo-welded material. 2 head-fixing-pillows with opening for access to the ears are integrated at the base and are part of the base itself if not used.

Allows fixed head immobilizer storage with stretcher. Complete of 2 chin straps and 2 Velcro straps, 8 belts with loop and 2 velcro bands for the fixing to the spinal board. X-ray translucent.

Dimensions: about 44,5x30x2,5 cm.

Weight: about 1,2 kg.

---

IMM121640 - WIND HEAD IMMOBILIZER

IMM121634 - SET OF 2 CHIN-STRAPS FOR WIND HEAD IMMOBILIZER


IMM121652 - SET OF 4 HEAD IMMOBILIZER BELTS

IMM121654 - SET OF 2 VELCRO STRAPS FOR WIND HEAD IMMOBILIZERS

---

## IMMOBILIZATION BELTS

---



### SPIDER IMMOBILIZATION BELTS

---

IMM121033 - Immobilization-restraint belts system. It can be used for all kind of stretcher and patient due to the mobility of the coloured horizontal belts alongside vertical one. Velcro fixing system ensuring an easy and fast use. Supplied with carrying bag.



### SET 2 METAL BUCKLE FIXING BELTS

---

IMM120282 - Set of 2 closed safety belts. Quick hook-unhook metal buckle. Dimensions: about 5x180 cm.



### SET 3 PLASTIC BUCKLE FIXING BELTS

---

IMM120280 - Set of 3 stretcher fixing belts. Quick hook-unhook plastic buckle. Complete with case with Velcro closing.



### SET OF 3 SAFETY BELT PAIRS

---

IMM122080 - Complete set of 3 pairs security belts (shoulder, chest, legs) for stretcher with loop. Rapid unhooked metal buckles. 2 stuffing passing through the shoulder belts for a better comfort. Adaptable to all type of stretchers.



### 2 SET 2 PLASTIC BUCKLE SAFETY BELTS

---

IMM120284 - Set of 2 2-pieces safety belts with loop. Quick hook-unhook plastic buckle for the use on automatic loading stretcher. Dimensions: about 5x176 cm.



### SET 2 2-PIECES METAL BUCKLE SAFETY BELTS


---

IMM120286 - Set of 2 2-pieces safety belts with loop. Quick hook-unhook metal buckle for the use on automatic loading stretcher. Dimensions: about 5x152 cm.

## CTS - CHILDREN TRANSPORT UNIT - IMM120050

### PEDIATRIC EMERGENCY TRANSPORT

CTS, innovative newborn's and children's containment and transport system, absolutely safe on ambulances. Fantasy fabric material, slightly stuffed.



### FROM 1 TO 24 MONTHS

The device is intended to be used for newborn and children transport, approximately from 1 to 24 months old (shoulder / crotch length from 35 to 57 cm).



### FIXED TO THE STRETCHER

Complete of safety-belts allowing to block the newborn in 3 positions, avoiding children usual dangerous and unexpected movements. 3 external safety belts allowing child further fastening. 3 belts allowing the fixing of the CTS to the stretcher. Complete of transport bag. CE medical device.  
Dimensions: about 110x70 cm.





# INFANT TRANSPORT UNIT - IMM120071

## NEWBORN EMERGENCY TRANSPORT

Newborn emergency transportation system. Grants a warm, protected and safe environment. Made in PVC material fantasy coloured. 3 transparent sides allowing complete and continuous monitoring of newborn by medical personnel.



## INNER FIXING BELTS

New edition PLUS complete with inner fixing belts. Internal mattress ensuring comfortable transport. Innovative folda-ble-bag system cordura material allowing, when not used, to be taken everywhere, always ready for use.




## ANCHORAGE TO THE STRETCHER

2 external handles for transportation. Arranged for the anchorage to the ambulance stretcher. Dimensions open: 90x38x34,5(h) cm. Dimensions closed: 90x38x15(h) cm.

## RESCUE STRETCHER

---



### SCOOP STRETCHER - 630

---

Aluminium scoop stretcher. X-ray translucent. Effective hook/unhook device for stretcher parting. Adjustable length with automatic block. Bag with n. 3 quick unhook safety belts for patient immobilization. Capacity: about 170 kg. Weight: about 10 kg. Height: about 7 cm. Width: about 43 cm.

Minimum Length: about 167 cm.

---

630 - SCOOP STRETCHER  
614 - HORIZONTAL STIRRUPS  
616 - VERTICAL STIRRUPS

---



### RESCUE STRETCHER TOBOGA - 16100

---

Rescue stretcher, extremely lightweight, suitable for patient transport in uncomfortable environment conditions. Made in polymeric ABS material with aluminium frame. Complete of expanded polyutherane mattress and adjustable footrest. 4 quick unhook security belts. 12 transport handles. Arranged for the anchorage to helicopters.

Dimensions: about 216x61x17h cm.  
Weight: 15 kg. Capacity: max 300 kg.

---

16100 - RESCUE STRETCHER TOBOGA  
16100-005 - LIFTING HARNESS  
16100-002 - FLOTATING SYSTEM

---



### 2 PIECES RESCUE STRETCHER- 16102

---

2 pieces rescue stretcher easy to store and transport. Made in polymeric ABS material with aluminium frame. Complete of expanded polyutherane mattress and adjustable footrest. 4 quick unhook security belts. 12 transport handles. Arranged for the anchorage to helicopters. Dimensions: about 216x61x17h cm. Weight: about 16 kg. Capacity: max 300 kg.

## RESCUE STRETCHER



### EMERGENCY STRETCHER X1

Emergency stretcher made of lightweight aluminium. Foldable in length with safety blocks on use position. Made in resistant fireproof, tearproof and washable material with central reinforce. Complete with 2 automatic unhook safety belts. Removable handles in plastic non-sliding material. Available with feet and wheels.

Dimensions (opened): 230x55,6x13,7 cm.

Dimensions (closed): 192x15x14,5 cm.

Weight: about 9 kg.

IMM121532 - X1 WITH 4 FEETS

IMM121534 - X1 WITH 4 WHEELS

IMM121530 - X1 WITH 2 FEETS AND 2 WHEELS



### EMERGENCY STRETCHER X2

Emergency stretcher made of lightweight aluminium. Foldable in length and width with safety blocks on use position. Made in resistant fireproof, tearproof and washable material with central reinforce. Complete with 2 automatic unhook safety belts. Removable handles in plastic non-sliding material. Available with feet and wheels.

Dimensions (opened): 230x55,6x13,7 cm.

Dimensions (closed): 104x15x14,5 cm.

Weight: about 9 kg.

IMM121538 - X2 WITH 4 FEETS

IMM121540 - X2 WITH 4 WHEELS

IMM121536 - X2 WITH 2 FEETS AND 2 WHEELS




### ANATOMIC PILLOW FOR STRETCHER - IMM122050

Anatomic pillow for stretcher. Made in antibacterial, fireproof, washable material. Completely thermo-welded.

Dimensions: about 40x27x7 cm.

Weight: about 300 g.



### CONFORMING PILLOW FOR STRETCHER - IMM122070

Conforming orange pillow for stretchers. Made in antibacterial, fireproof, washable material. Completely thermo-welded.


Dimensions: about 38x25x10/14h cm.

Weight: about 400 g.

## SEDAN CHAIR EXTRA-670/BR

### SEDAN CHAIR

EXTRA is the new stair chair with tracked system to descend stairs. Entirely made of lightweight aluminium alloy frame. The backrest and back carrying handles inclination, along with the unique frontal arms configuration, enable a perfect weights distribution, relieving operators' efforts.





### OPTIONAL AND INCLUDED

It is completed with two practical armrests and a telescopic rear bar covered by an high comfort padded sheath. Turnable footrest with magnetic locks and innovative bio-mechanical solutions, complete the characteristics of this new rescue device. As optional, a version with brake on the back wheels is available.

### COMFORT

Long and high-grip tracks allow the EXTRA stair-chair to be safely manoeuvred and handled on 3 stairs when descending with a seated patient. While manoeuvring on smooth and plane surfaces, the 200 mm rear and 150 mm front wheels give exceptional comfort.



Height: 101 cm. Width: 55 cm. Length: 81 cm. Max length: 139 cm. Folded thickness: 29 cm. Weight: 15 kg. Loading capacity: 200 kg.



## MERCURY LITE SELFLOADING STRETCHER - 7080/4RG PROOF

### 4 ROTATING WHEELS

This version is available for "Mercury Lite" with 4 rotating wheels. It grants extremely maneuverability. This product is ideal for narrow spaces. Operating an easy command the two frontal wheels can be unlocked or reset to their original configuration.



### FIXING SYSTEM

Fixing of stretcher inside the rescue vehicle is granted by a certified system. Supplied complete with mattress and belts.



### DIMENSIONS

Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789. Material: painted aluminium HDPE (High-density polyethylene).

Intermediate height: 74 cm. Rear part height: 85 cm. Max loading height: 66 cm. Min loading height: 64 cm. Folded height: 28 cm. Width: 57 cm. Length: 196 cm. Length with folded trolley: 162 cm. Weight: 36,5 kg. Loading capacity: 250 kg.



### ALL-IN-ONE RESUSCICATION BAGS

---

The resuscitation bags are available in the new ALL-IN-ONE version with a valve directly integrated into the bottom of the bag. This new feature allows the use of reservoirs not equipped with non-return valve. The model with external valve so far marketed is available only for pediatric and neonatal bags.

## SILICONE RESUSCITATOR BAGS

---



### ADULT

---

Transparent silicone resuscitator bag, autoclavable up to 135°C under the pressure of 2 bar. Volume 1600 ml. Complete with 360° rotating connection patient's valve. Adjustable safety valve to limit the overpressure applied by the patient. New all-in-one valve directly integrated into the bottom of the resuscitator bag, that allows the use devoid of a non-return valve (less bulky and lighter). Complete version supplied with mask no. 5, reservoir and oxygen tube. Latex free. CE certified in class IIa.

---

PAL34000 - ADULT SILICONE  
PAL34004 - ADULT SILICONE  
COMPLETE

---



### PEDIATRIC

---

Transparent silicone pediatric resuscitator bag, autoclavable up to 135°C under the pressure of 2 bar. Volume 500 ml. Complete with 360° rotating connection patient's valve. Adjustable safety valve to limit the overpressure applied by the patient. New all-in-one valve directly integrated into the bottom of the resuscitator bag, that allows the use devoid of a non-return valve (less bulky and lighter). Complete version supplied with mask no. 3, reservoir and oxygen tube. Latex free. CE certified in class IIa.

---

PAL34010 - PEDIATRIC SILICONE  
PAL34014 - PEDIATRIC SILICONE  
COMPLETE

---



### NEWBORN

---

Transparent silicone newborn resuscitator bag, autoclavable up to 135°C under the pressure of 2 bar. Volume 240 ml. Complete with 360° rotating connection patient's valve. Adjustable safety valve to limit the overpressure applied by the patient. New all-in-one valve directly integrated into the bottom of the resuscitator bag, that allows the use devoid of a non-return valve (less bulky and lighter). Complete version supplied with mask no. 1, reservoir and oxygen tube. Latex free. CE certified in class IIa.

---

PAL34020 - NEWBORN SILICONE  
PAL34024 - NEWBORN SILICONE  
COMPLETE  
PAL32020 - NEWBORN SILICONE  
WITH EXTERNAL VALVE  
PAL32024 - NEWBORN SILICONE  
COMPLETE WITH EXTERNAL VALVE

---

## DISPOSABLE RUBBER RESUSCITATOR BAGS



### ADULT

Adult non-autoclavable rubber resuscitator bag. Volume 1600 ml. Complete with 360° rotating connection patient's valve. Adjustable safety valve to limit the overpressure applied by the patient. New all-in-one valve directly integrated into the bottom of the resuscitator bag, that allows the use devoid of a non-return valve (less bulky and lighter). Complete version supplied with mask size 5 with reservoir and shockproof oxygen tube. Latex free. To be disinfected in bactericide solution. CE certified in class IIa.

PAL34200 - ADULT RUBBER  
PAL34204 - ADULT RUBBER  
COMPLETE



### PEDIATRIC

Pediatric non-autoclavable rubber resuscitator bag. Volume 500 ml. Complete with 360° rotating connection patient's valve. Adjustable safety valve to limit the overpressure applied by the patient. New all-in-one valve directly integrated into the bottom of the resuscitator bag, that allows the use devoid of a non-return valve (less bulky and lighter). Complete version supplied with mask size 3 with reservoir and shockproof oxygen tube. Latex free. To be disinfected in bactericide solution. CE certified in class IIa.

PAL34210 - PEDIATRIC RUBBER  
PAL34214 - PEDIATRIC RUBBER  
COMPLETE  
PAL32210 - PEDIATRIC RUBBER  
WITH EXTERNAL VALVE  
PAL32214 - PEDIATRIC RUBBER  
COMPLETE WITH EXTERNAL VALVE



### NEWBORN

Newborn non-autoclavable rubber resuscitator bag. Volume 240 ml. Complete with 360° rotating connection patient's valve. Adjustable safety valve to limit the overpressure applied by the patient. New all-in-one valve directly integrated into the bottom of the resuscitator bag, that allows the use devoid of a non-return valve (less bulky and lighter). Complete version supplied with mask size 1 with reservoir and shockproof oxygen tube. Latex free. To be disinfected in bactericide solution. CE certified in class IIa.

PAL34220 - NEWBORN RUBBER  
PAL34224 - NEWBORN RUBBER  
COMPLETE  
PAL32220 - NEWBORN RUBBER  
WITH EXTERNAL VALVE  
PAL32224 - NEWBORN RUBBER  
COMPLETE WITH EXTERNAL VALVE



## DISPOSABLE PVC RESUSCITATOR BAGS



### ADULT

Disposable economic resuscitator bag. Made in PVC material. Volume 1600 ml. Complete with 360° rotating connection patient's valve. Adjustable safety valve to limit the overpressure applied by the patient. New all-in-one valve directly integrated into the bottom of the resuscitator bag, that allows the use devoid of a non-return valve (less bulky and lighter). Complete version supplied with mask size 5 with reservoir and shockproof oxygen tube. Latex free. CE certified in class IIa.

PAL34100 - ADULT PVC  
PAL34104 - ADULT PVC COMPLETE



### PEDIATRIC

Disposable economic resuscitator bag. Made in PVC material. Volume 500 ml. Complete with 360° rotating connection patient's valve. Adjustable safety valve to limit the overpressure applied by the patient. New all-in-one valve directly integrated into the bottom of the resuscitator bag, that allows the use devoid of a non-return valve (less bulky and lighter). Complete version supplied with mask size 3 with reservoir and shockproof oxygen tube. Latex free. CE certified in class IIa.

PAL34110 - PEDIATRIC PVC  
PAL34114 - PEDIATRIC PVC COMPLETE  
PAL32110 - PEDIATRIC PVC WITH EXTERNAL VALVE  
PAL32114 - PEDIATRIC PVC COMPLETE WITH EXTERNAL VALVE



### NEWBORN

Disposable economic resuscitator bag. Made in PVC material. Volume 240 ml. Complete with 360° rotating connection patient's valve. Adjustable safety valve to limit the overpressure applied by the patient. New all-in-one valve directly integrated into the bottom of the resuscitator bag, that allows the use devoid of a non-return valve (less bulky and lighter). Complete version supplied with mask size 1 with reservoir and shockproof oxygen tube. Latex free. CE certified in class IIa.

PAL34120 - NEWBORN PVC  
PAL34124 - NEWBORN PVC COMPLETE  
PAL32120 - NEWBORN PVC WITH EXTERNAL VALVE

## VENTILATION MASKS

---



### SILICON MASKS

---

Transparent ventilation mask with silicone autoclavable lip-board ensuring maximum adherence to the face. The transparent shell allows patient's monitoring during ventilation.

---

PAL32050 - SILICONE MASK SIZE 0  
PAL32054 - SILICONE MASK SIZE 1  
PAL32058 - SILICONE MASK SIZE 2  
PAL32062 - SILICONE MASK SIZE 3  
PAL32066 - SILICONE MASK SIZE 4  
PAL32070 - SILICONE MASK SIZE 5

---



### DISPOSABLE RUBBER MASKS

---

Anatomic ventilation mask with blue rubber rim. Transparent shell allowing patient's monitoring during ventilation.

---

PAL32250 - RUBBER MASK SIZE 3  
PAL32254 - RUBBER MASK SIZE 4  
PAL32258 - RUBBER MASK SIZE 5

---



### DISPOSABLE PVC MASKS

---

Disposable ventilation mask with turning standard connection (22 mm). Shell shaped transparent plastic material and border/cushion in soft inflated plastic.

---

PAL32150 - PVC MASK SIZE 0  
PAL32154 - PVC MASK SIZE 1  
PAL32158 - PVC MASK SIZE 2  
PAL32162 - PVC MASK SIZE 3  
PAL32166 - PVC MASK SIZE 4  
PAL32170 - PVC MASK SIZE 5

---

## FILTER AND RESERVOIR

---



### PATIENT VALVE

---

Autoclavable patient valve for resuscitator bags OB 360° rotating connector. Adjustable safety over pressure valve.

---

PAL32370 - AUTOCLAVABLE PATIENT VALVE  
PAL32375 - MEMBRANE FOR POS. 15 VALVE

---



### DISPOSABLE RESERVOIR WITH EXTERNAL VALVE

---

Disposable reservoir for resuscitator bags OB with external valve.

---

PAL32350 - 2500 ml  
PAL32356 - 900 ml

---



### PATIENT VALVE WITH PAL34370 PEEP VALVE CONNECTION

---

Autoclavable patient valve for OB resuscitation bags 360° spinning connection. Adjustable over-pressure security valve. Equipped with PAL34371 Peep valve connection.

---

PAL34370 - PATIENT VALVE WITH PAL34370 PEEP VALVE CONNECTION

---



### DISPOSABLE RESERVOIR WITH INTEGRATED VALVE

---

Disposable reservoir for resuscitator bags OB with integrated valve.

---

PAL34350 - 2500 ml  
PAL34356 - 900 ml

---



### PEEP VALVE - PAL34371

---

Peep Valve manufactured for patient valve PAL34370. Built in polycarbonate, autoclavable with steam for 30 cycles. Adjustable from 2,5 to 20 cmH<sub>2</sub>O. It cannot be used on valves of other types or brands.



### DISPOSABLE FILTER - 3505879

---

Disposable filter for resuscitator bags and automatic respirators. Sterile package.

## OXYGEN THERAPY MASKS

---



### OXYGEN THERAPY MASK - 11914020304

---

Adult oxygen therapy mask. Manufactured in transparent, atoxic and soft material complete with adjustable elastic. Supplied with 210 cm crushproof tube . Singulary packaged.



### PEDIATRIC OXYGEN THERAPY MASK - 11914020302

---

Pediatric oxygen therapy mask made of transparent, atoxic and soft material complete with adjustable elastic. Supplied with 210 cm crushproof tube. Singulary packaged.



### ADULT OXYGEN THERAPY MASK WITH RESERVOIR - 11914020313

---

Oxygen therapy mask with reservoir and non-return valve. Made of transparent, atoxic and soft material complete with adjustable elastic. High concentration (60% Oxygen). Metal nose clip allowing a safe and comfortable appliance. Supplied with 210 cm crushproof tube and universal connector. Size: adult nr. 6.



### PEDIATRIC OXYGEN THERAPY MASK WITH RESERVOIR - 11914020311

---

Oxygen therapy mask with reservoir and non-return valve. Made of transparent, atoxic and soft material complete with adjustable elastic. High concentration (60% Oxygen). Metal nose clip allowing a safe and comfortable appliance. Supplied with 210 cm crushproof tube and universal connector. Singulary packaged.



## OXYGEN THERAPY MASKS

---




### ADULT VARIABLE CONCENTRATION MASK

- 11914020328

---

Mask for oxygen therapy with variable concentration. Made of transparent, non-toxic and soft material, complete with adjustable elastic. Metal nose clip for safe and comfortable application. Equipped with a 210 cm anti-crushing hose and an adjustable concentration connection. Size: adults.



### OXYGEN THERAPY NASAL CANULA 180CM

- 11914020101

---

Disposable oxygen therapy nasal canula. Manufactured in atoxic and soft material. It allows freedom of movement and do not interfere with feeding and vocal activities. Supplied with a 180 cm tube in singularly sealed package.



### CRUSHPROOF OXYGEN TUBE FOR RESUSCITATOR BAGS

- 11914020204

---

Oxygen 2 m crushproof tube extension complete with male/female connectors.

## VENTILATION SETS

---



### REANIMATION BAG

---

Reanimation bag with shoulder strap for transport, featuring: resuscitator bag with turning connection security valve, 2 ventilation masks, reservoir, connecting tube for oxygen, 5 Guedel airways.

Dimensions: about 28x15x16 cm.

Weight: about 1 kg.

---

VAL40910 - REANIMATION BAG - SILICONE

VAL40914 - REANIMATION BAG - RUBBER

VAL40916 - REANIMATION BAG - PVC

---



### RESUSCITATION CASE SILICONE OB - VAL40902

---

ABS resuscitation case with wall anchorage system. 1 autoclavable silicone resuscitator bag. Reservoir and oxygen tube. 3 silicone transparent ventilation masks. 5 Guedel airways. Dimensions: 39,5x13,5x26(h) cm. Weight: about 2,3 kg.



### RESUSCITATION CASE DISPOSABLE RUBBER OB - VAL40904

---

ABS resuscitation case with wall anchorage system. 1 rubber disposable resuscitator bag. Reservoir and oxygen tube. 3 disposable ventilation masks. 5 Guedel airways.

Dimensions: about 39,5x13,5x26(h) cm.

Weight: about 2,2 kg.



### RESUSCITATION CASE DISPOSABLE OB - VAL40905

---

ABS resuscitation case with wall anchorage system. 1 pvc disposable resuscitator bag. Reservoir and oxygen tube. 3 disposable ventilation masks. 5 Guedel airways.

Dimensions: about 39,5x13,5x26(h) cm.

Weight: about 2,2 kg.

## PEDIATRIC-NEWBORN VENTILATION SETS

---



### NEWBORN REANIMATION CASE - VAL40715

---

Shockproof ABS case with wall bracketing system. 1 newborn resuscitator bag, 1 paediatric resuscitator bag, 3 resuscitator masks, 4 Guedel airways, 2 mucus aspirators, 3 suction catheters, 1 paediatric laryngoscope with 2 blades, 1 stylet, 1 Magill forceps, 5 endotracheal tubes, 1 paediatric sphygmomanometer with cuff, 1 stethoscope, 1 diagnostic pen, 5 sterile syringes, 2 IV administration sets, 6 sterile disposable butterfly needles, 1 Alumed cloth, 6 Alumed compresses, 3 elastic bandages, 1 baby isothermal blanket, 1 medicated plaster, 1 plaster roll, 1 tourniquet, 1 Iodopovidone bottle, 4 disinfectant towelettes, 2 umbilical forceps, 1 anatomic forceps, 1 Mosquito forceps, 1 instrument case, 1 disposable scalpel, 1 pair of scissors, 1 haemostatic forceps, 1 medication package, 1 trash bag.

Dimensions: about 42x15x31(h) cm.

Weight: about 5,9 kg.



### PEDIATRIC/NEWBORN RESUSCITATOR BAGS SETS

---

Set of 2 silicone resuscitator bags (newborn 240ml, pediatric 500ml) featuring: 3 ventilation masks, 1 reservoir, 1 oxygen tube, 3 Guedel airways.  
Dimension: about 25x25x8(h) cm.

---

VAL44240 - SILICONE SET  
VAL44280 - PVC SET  
VAL44300 - RUBBER SET

---



## RELIABILITY

---

Diagnostic devices selected from the best manufacturers available on the market. All the devices we market are scrupulously checked according to documented procedures to guarantee quality and reliability.

## DAE MADE IN EUROPE

---

Our DAE defibrillator completes this range of devices, as the result of a collaboration with a well-known and specialized European manufacturer.



## EXTERNAL SEMI-AUTOMATIC DEFIBRILLATOR OBLS300

### MADE IN EUROPE

Made in Europe, available as semi-automatic defibrillators (ref. DAE300) or automatic (ref. DAE300A).

PADs shelf life from production date: 30 months.

Internal memory capacity: 2 hours of ECG signals /500 events (memory card).

Battery capacity: 180 shock/ 5 years shelf life. Daily or weekly self-test.



### MULTIPUSE BIOWAVE®

The defibrillator features the biphasic pulse defibrillation impulse Multipulse Biowave® patented on a worldwide level, which allows more efficient defibrillation. The defibrillator measures impedance and adjusts energy delivery accordingly.



### DISPLAY AND DIMENSIONS

High-resolution LCD display allows the user to notice and evaluate immediately possible faults. Besides, it permits to easily visualize text instructions concerning reanimation and defibrillation interventions.

Dimensions: 70x272x240 mm.

Weight: 1,5 kg.

DAE300 - OB LIFE SUPPORT 300 SEMI-AUTOMATIC  
DAE300A - OB LIFE SUPPORT 300A AUTOMATIC

## DEFIBRILLATOR ACCESSORIES

---



### ADULT PADS - F7956W

---

PADs for adults.  
Suitable for defibrillator OB Life Support 300. PADs with blue connector should be used for patients, whose weight is more than 25 kg.



### PEDIATRIC PADS - F7956P

---

PADs for children.  
Suitable for defibrillator OB Life Support 300. PADs with yellow connector should be used for children, whose weight is less than 25 kg (younger than 8 years of age).



### DEFI SET - MED98950

---


Ready to use set for early defibrillation, consisting of:  
Clothes-cutting scissors.  
Trichotomy razor. 1 pair of disposable gloves. 2 disinfectant wipes. 2 non-woven absorbing compresses 10x10 cm. Life key ventilation mask.



### TRAINER AED OBLS 300 - DAE290

---

Please see section "Other, accessories and didactics".



### SAFETY SIGN OF AED - ACC170551

---

Please see section "Other, accessories and didactics".



### VISIO DEFIBRILLATOR METAL CABINET - DEF040

---

Please see section "Emergency bags and rucksacks".

## STETHOSCOPE

---



### SINGLE HEAD STETHOSCOPE

---

Compact and light flat-head stethoscope. Colored “Y” tube made of a single PVC piece. Very sensitive membrane with excellent tone transmission. CE marked according to the latest legal provisions concerning medical devices.

DIA01012B - SINGLE HEAD BLUE  
DIA01012P - SINGLE HEAD VIOLET  
DIA01012PK - SINGLE HEAD PINK  
DIA01012R - SINGLE HEAD RED  
DIA01012Y - SINGLE HEAD YELLOW



### SINGLE HEAD SATINED MASTER STETH. - DIA01020B

---

Professional single-head stethoscope. Guarantees a highly sensitive level of auscultation. Head made of galvanized aluminum alloy with satin finish that guarantees lightness and durability. Excellent sound transmission thanks to the soft and flexible PVC tube. Excellent finish and ease of grip. Equipped with a flexible stainless-steel bow. The package includes a pair of rigid plastic cups and a soft PVC pair. CE marked according to the latest legal provisions concerning medical devices.



### CHICKEN STETHOSCOPE COVER - DIA106908

---

Washable soft and plush stethoscope cover suitable for standard size stethoscope.



### DUAL HEAD STETHOSCOPE

---

Compact and light dual head stethoscope. Colored “Y” tube made of a single PVC piece. Very sensitive membrane with excellent tone transmission. CE marked according to the latest legal provisions concerning medical devices.

DIA106746 - DUAL HEAD YELLOW  
DIA01016 - DUAL HEAD BLACK

## SPHYGMOMANOMETERS

---



### ANEROID SPHYGMOMANOMETER

---

Red or blue aneroid sphygmomanometer. High visibility pressure gauge with graduated scale up to 300 mm/Hg. Nylon cuff with “D” ring and Velcro closing (allows a rapid and fast fastening). Separate inflation bulb with pressure regulation and release valve. Offered with protection and transport case and user manual. Measuring accuracy:  $\pm 3$  mmHg. CE marked according to the latest legal provisions concerning medical devices.

DIA02003RE - RED ANEROID SPHYGMOM.  
DIA02003BL - BLUE ANEROID SPHYGMO.

---



### PEDIATRIC PALMAR SPHYGMOMANOMETER - DIA100070

---

Aneroid palmar sphygmomanometer. Big manometer with graduated scale up to 300 mm/Hg. Complete with 2 cuffs, newborn and paediatric in fantasy fabric material with velcro closing. One-handed operations allowed, accurate regulation valve. Supplied with fantasy fabric case. Manometer Ø 70 mm. CE marked according to the latest legal provisions concerning medical devices.



### BLUE PALMAR SPHYGMOMANOMETER - DIA02022

---

Palmar sphygmomanometer with blue cuff. Large pressure gauge with graduated scale up to 300 mm/Hg (diameter 65mm). Nylon cuff with Velcro closing for rapid use. Inflation bulb with pressure regulation valve. Measuring accuracy:  $\pm 3$  mmHg. Supplied with protection and transport case. CE marked according to the latest legal provisions concerning medical devices.



## SPHYGMOMANOMETERS



### TABLE AUTOMATIC SPHYGMOMANOMETER - DM490

Innovative design with easy-to-read LCD display (size: 4.6 x 6.3cm) which simultaneously displays the rapid and precise measurement of systolic and diastolic pressure, the heart rate, date/hour and detection of irregular IHD heart rate (Arrhythmia). One-user-only memory, 120 slots in total. Oscillometric reading method with indication according to the international standard (WHO). Armlet with Velcro closure with automatic inflation. Power: 4 batteries (AAA) 1.5 V. Size: length 14cm; depth 9.8cm; height 4,8 cm. Weight: 227g (without batteries).



### BLACK PALMAR SPHYGMOMANOMETER - DIA02008

Palmar sphygmomanometer with highly resistant pressure gauge in the event of a hit. Large pressure gauge with graduated scale up to 300 mm/Hg (diameter 65mm). Nylon cuff with Velcro closing for rapid use. Inflation bulb with pressure regulation valve. Measuring accuracy:  $\pm 3$  mmHg. Supplied with protection and transport case. CE marked according to the latest legal provisions concerning medical devices.



### PALMAR SPHYGMOMANOMETER DELUXE 1 TUBE- DIA02005A

Professional palmar sphygmomanometer with 1 tube. High visibility pressure gauge (65 mm) with graduated scale up to 300 mm/Hg. Nylon cuff with Velcro closing for a rapid use. Inflation bulb with pressure regulation and rapid release valve. Measuring accuracy:  $\pm 3$  mmHg. Right and left use (by rotating the gauge). Practical protection and transport case. CE marked according to the latest legal provisions concerning medical devices.

## SPHYGMOMANOMETERS

---



### SPHYGMOMANOMETER WITH STETHOSCOPE - DIA02002

---

Black aneroid sphygmomanometer complete with integrated stethoscope. Greater flexibility and ease of use. High visibility pressure gauge with graduated scale up to 300 mm/Hg. Black nylon cuff with "D" ring and Velcro closing. Separate inflation bulb with pressure regulation and release valve. Offered with protection and transport case and user manual. Measuring accuracy:  $\pm 3$  mmHg. CE marked according to the latest legal provisions concerning medical devices.



### TABLE ANEROID SPHYGMOMANOMETER - DIA02009

---

Circular wall and/or table aneroid sphygmomanometer (diameter 170mm). Scale from 0 to 300 mm/Hg. Measuring accuracy:  $\pm 3$  mmHg. Velcro-closing cuff and inflation bulb, spiral connecting tube. Wall anchor support with basket bracelet holder. CE marked according to the latest legal provisions concerning medical devices.

## OMRON SPHYGMOMANOMETERS



### WRIST SPHYGMOMANOMETER OMRON RS4 - OMRS4

Electronic wrist blood pressure monitor Omron RS4 HEM-6181-E is fully automatic.

The device stores up to 60 measurements in memory and calculates an average value based on the last three measurements taken within 10 minutes.

Dimensions of the wrist: 135-215 mm.

Big display to ensure better readability. The values of systole, diastole and pulse are displayed simultaneously. AP sensor detects correct placement of the cuff. "IntelliSense" technology.

Rigid storage case.



### OMRON M6 SPHYGMOMANOMETER - OMM6COMIT

Digital Automatic Blood Pressure Monitor with computer connection. Compact, fully automatic, user-friendly and quick. Measures blood pressure and pulse rate. Detects irregular heartbeats. Memory for 100 sets with date and time for each user (1 and 2). Digital LCD display. "IntelliSense" technology. Storage case. CE marked according to the latest legal provisions concerning medical devices.



### OMRON M2 SPHYGMOMANOMETER - OMM2EC

Digital Automatic Blood Pressure Monitor. Compact, fully automatic, user-friendly and quick. Measures blood pressure and pulse rate. "Detects irregular heartbeats. Memory for 30 sets. Digital LCD display. "IntelliSense" technology. Storage case. CE marked according to the latest legal provisions concerning medical devices.



### OMRON M3 SPHYGMOMANOMETER - OMM3BC

Digital Automatic Blood Pressure Monitor. Compact, fully automatic, user-friendly and quick. Measures blood pressure and pulse rate. Irregular heartbeats detection. Movement error detection. Average value function. Memory for 60 sets with date and time for each user (1 and 2). Digital LCD display. Storage case. CE marked according to the latest legal provisions concerning medical devices.

## PULSE OXIMETERS

---



### ADULT FINGER PULSE OXIMETER - FS20E

---

Device for monitoring blood oxygen saturation. Compact and lightweight, equipped with colour LED display. Detects oxygen saturation (SpO2) and pulse rate. Allows 4 display modes. Automatic switch-off if the device is not used. Battery level indicator. Supplied with 2 AAA batteries, allows up to 600 readings at 25° C. Dimensions: about 58x33x35 mm. Maximum weight: 100g (batteries included). CE medical device. Isolation class: II. Class of applied part: BF. Italian user manual. Measurements must be followed by further specific tests performed by competent practitioners.

---

FS20E - ADULT FINGER PULSE OXIMETER  
DIA109521 - CASE FOR FINGER OXIMETER

---



### HANDHELD PULSE OXIMETER - HS10A

---

Lightweight for carrying and Easy-To-Use. Silicon rubber shell protection and stable metal bracket for table usage. Adult, pediatric and neonate oximeter probe available. CE medical device. Isolation class: II. Protection against electric shock: type BF – applied part. Protection against the penetration of fluids and substances: IPX1. Visual and acoustic alarm available. Big size 2.4" LCD display with SPO2, pulse ratio, pulse bar-graph and plethysmogram. SPO2 resolution 1 %. Low Battery voltage indicator. Refreshing rate: 1 sec. The unit can work with alkaline batteries (2000 mAh) or rechargeable Lion battery pack (1600 mAh). User manual available in English, Italian and German languages.

---

HS10A - HANDHELD PULSE OXIMETER

---



## DIAGNOSTIC PENLIGHTS



### OB DIAGNOSTIC PENLIGHT

Rechargeable diagnostic bright penlight. Built in plastic material. Pocket clip. Complete with 2 batteries AAA.

DIA101010 - OB DIAGNOSTIC PENLIGHT



### DIAGNOSTIC METAL PENLIGHT

Rechargeable diagnostic bright penlight. Built in metal material. Pocket clip. Complete with two batteries AAA.

DIA101011 - DIAGNOSTIC METAL PENLIGHT



### DIAGNOSTIC PENLIGHT WITH LED

Rechargeable diagnostic bright penlight with LED light. Built in metal material. Pocket clip. Complete with 2 batteries AAA.

DIA101012 - DIAGNOSTIC PENLIGHT WITH LED



### PEN LIGHT CLIP

Tongue depressor clip for diagnostic pen, allows one-hand operation.

DIA101020 - PEN LIGHT CLIP



### PINK DIAGNOSTIC PENLIGHT

Rechargeable diagnostic bright penlight. Built in plastic material. Pocket clip. Complete with 2 batteries AAA.

DIA101016 - OB PINK



### SET 6 DIAGNOSTIC PENLIGHTS

Set of 6 disposable diagnostic pens with filament lamps. Made of plastic material and water resistant. To turn on the light simply press the clip on the pen itself. Disposable device, battery not replaceable.

DIA05004 - SET 6 PCS DISPOSABLE DIAGNOSTIC PENLIGHTS

## THERMOMETERS

---

### THERMOSCAN 7 THERMOMETER - 26IRT6520

---



Instantaneous auricular thermometer BRAUN THERMOSCAN® 7, with Age Precision® IRT6520. Just set the child's age, measure the temperature and the color code displayed will help you interpret the temperature value and act accordingly. The ExacTemp® positioning system confirms the right placement and the measurement precision by lighting an indicator light and producing an acoustic signal. Infrared heat detector through tympanum in one second. Complete with 20 pcs lens covers. CE medical device. Internal and external temperature difference is 0,5 °C.

---

26IRT6520 - THERMOSCAN 7 THERMOMETER  
DIA104725 - EAR TIPS BOX 40 PCS

---



### FLEXIBLE DIGITAL THERMOMETER JUMBO 30 SEC - RM310

---

Flexible digital thermometer.  
Range: 32°C÷42°C, precision: 0,1°C  
Self-shutting system. Waterproof. The internal temperature is higher than the one outside by +/- 0,5 degrees. CE medical device.  
Dimension: 139x27x12 mm.  
Weight: 10 g.



### MERCURY FREE THERMOMETER- TER174

---

Clinical mercury-free thermometer for body temperature detection. The liquid used is a "Gallio/Indio" alloy combination, totally non-toxic and safe for the environment. CE medical device. Detection range: 35-42° C, precision: 0,10° C  
Dimensions: about 125x12x8 mm.

## DIAGNOSTICS DEVICES AND ACCESSORIES

---



### FIBER-OPTIC OTOSCOPE DIMED - DMV520

---

High quality otoscope for pocket and doctor's bag. Unobstructed view, unit equipped with fibre-optic light. Robust and light weight plastic with metal inserts and 3x rotating magnifying glass. Handle requires 2 batteries (AA) and it is equipped with grooves for a better grip. Speculum input is rapid and stable. Provided with bag and with 4 additional reusable specula.



### MINI OTOSCOPE WITH CONVENTIONAL LIGHT - DMV540B

---

High quality blue otoscope for pocket and doctor's bag. Unobstructed view, unit equipped with halogen bulb. Robust and light weight plastic with metal inserts and 3x rotating magnifying glass. Handle requires 2 batteries (AA) and it is equipped with grooves for a better grip. Speculum input is rapid and stable. Provided with bag and with 4 additional reusable specula.

## DIAGNOSTICS DEVICES AND ACCESSORIES


---



### PRESSURE INFUSER WITH OB MANOMETER - DIA100390

---

Pressure infuser for medicaments contained in plastic material sachets. Complete with manometer for infusion monitoring and mono-ball with deflating valve for accurate dosage. Suitable for 500 ml infusion. CE medical device.



### BLACK OPHTHALMOSCOPE - DIA100415

---


Fiber optic 2,5V bulb. Ophthalmoscope: head with screw blocking; corrective lenses from 0 to +20 and 0 to -20 dioptres with 5 different diaphragm apertures. Plastic/metal handle with batteries; working with 2 batteries type AA (not included). Complete with shockproof plastic case. CE medical device.



### WOODEN TONGUE DEPRESSOR - 90554

---

Wooden tongue depressor. Separately packaged. Dimensions: 150x17 mm.



### CAMERA D'ARIA

---

Adult 1 tube bladder: 22x12 cm.  
Adult 2 tubes bladder: 22x12 cm.  
Newborn 2 tubes bladder: 6x12 cm.

---

DIA102100 - ADULT 1 TUBE BLADDER  
DIA102104 - ADULT 2 TUBES BLADDER  
DIA102108 - INFANT 2 TUBES BLADDER

---



## DIAGNOSTICS DEVICES AND ACCESSORIES

---



### ECG ELECTRODES

---

E.C.G. electrodes. Base made of foam material. Conductive paste. Optimal detection of the tracing.

Adult: diameter 55 mm.

Pediatric: diameter 43 mm.

50 pcs package.

---

F9060 - ADULT ECG ELECTRODES 50 pcs

F9060P - PEDIATRIC ECG ELECTRODES 50 pcs

---



### ONE-SHOULDER WITH BACK VALVE - DIA 102212

---

One-shoulder with back valve for standard sphygmomenometer.



### ONE-SHOULDER ARI TAP - DIA 102208

---

Air deflating tap for standard sphygmomanometer



## REANIMATION

---

We have a range of products for the most critical situations: medical devices for reanimation and specific kits. Our goal of guaranteeing reliability and safety for patients and users is achieved through constant quality controls.

## CONVENTIONAL LIGHT LARYNGOSCOPE

---



### MC-INTOSH LARYNGOSCOPE

---

Mc Intosh blades for conventional light laryngoscope, made in steel. Different sizes available. Sets with autoclavable laryngoscope available with 3 blades (sizes from 1 to 3) or 4 blades (sizes from 0 to 3). Sets are supplied with a practical roll case.

---

RIA53000 - 3 BLADES SET MC-INTOSH CONVENTIONAL LIGHT  
RIA53004 - 4 BLADES SET MC-INTOSH CONVENTIONAL LIGHT  
RIA53020 - MC-INTOSH CONVENTIONAL LIGHT BLADE SIZE 0  
RIA53022 - MC-INTOSH CONVENTIONAL LIGHT BLADE SIZE 1  
RIA53024 - MC-INTOSH CONVENTIONAL LIGHT BLADE SIZE 2  
RIA53026 - MC-INTOSH CONVENTIONAL LIGHT BLADE SIZE 3  
RIA53028 - MC-INTOSH CONVENTIONAL LIGHT BLADE SIZE 4

---



### MILLER LARYNGOSCOPE

---

Conventional light metal autoclavable laryngoscope with straight Miller blades in steel for conventional light laryngoscope. Blades are available in different sizes. Available in set with 3 blades (sizes from 0 to 2) and a practical roll case.

---

RIA53100 - 3 BLADES SET MILLER CONVENTIONAL LIGHT  
RIA53120 - MILLER CONVENTIONAL LIGHT BLADE SIZE 0  
RIA53122 - MILLER CONVENTIONAL LIGHT BLADE SIZE 1  
RIA53124 - MILLER CONVENTIONAL LIGHT BLADE SIZE 2  
RIA53126 - MILLER CONVENTIONAL LIGHT BLADE SIZE 3  
RIA53128 - MILLER CONVENTIONAL LIGHT BLADE SIZE 4

---

## FIBER OPTIC LARYNGOSCOPE

---



### MC-INTOSH LARYNGOSCOPE

---

Mc Intosh fiber optic blades made in steel for fiber optic laryngoscope. Different sizes available. 2,5V halogen lamp set inside the blade for an easy sterilization. Sets with laryngoscope available with 3 blades (sizes from 1 to 3) or 4 blades (sizes from 0 to 3). Sets are supplied with a practical roll case.

---

RIA53200 - 3 BLADES SET MC-INTOSH FIBER OPTIC  
RIA53204 - 4 BLADES SET MC-INTOSH FIBER OPTIC  
RIA53220 - BLADE MC-INTOSH FIBER OPTIC SIZE 0  
RIA53222 - BLADE MC-INTOSH FIBER OPTIC SIZE 1  
RIA53224 - BLADE MC-INTOSH FIBER OPTIC SIZE 2  
RIA53226 - BLADE MC-INTOSH FIBER OPTIC SIZE 3  
RIA53228 - BLADE MC-INTOSH FIBER OPTIC SIZE 4  
RIA53230 - BLADE MC-INTOSH FIBER OPTIC SIZE 5

---



### MILLER LARYNGOSCOPE

---

Miller fiber optic blades made in steel for fiber optic laryngoscope. Different sizes available. 2,5V halogen lamp set inside the blade for an easy sterilization. Sets with laryngoscope available with 2 blades (sizes from 0 to 1) or 3 blades (sizes from 0 to 2). Sets are supplied with a practical roll case.

---

RIA53300 - 2 BLADES SET MILLER FIBER OPTIC  
RIA53304 - 3 BLADES SET MILLER FIBER OPTIC  
RIA53320 - BLADE MILLER FIBER OPTIC SIZE 0  
RIA53322 - BLADE MILLER FIBER OPTIC SIZE 1  
RIA53324 - BLADE MILLER FIBER OPTIC SIZE 2  
RIA53326 - BLADE MILLER FIBER OPTIC SIZE 3

---



## FIBER OPTIC LARYNGOSCOPE WITH DISPOSABLE BLADES



### MC-INTOSH DISPOSABLE BLADES

Disposable fiber optic Mc Intosh laryngoscope blade, made in plastic material. Different sizes available.

Size 1: 92x15(h)

Size 2: 112x15(h)

Size 3: 132x25(h)

Size 4: 165x30(h)

RIA53444 - DISPOSABLE MC INTOSH BLADE SIZE 1

RIA53446 - DISPOSABLE MC INTOSH BLADE SIZE 2

RIA53448 - DISPOSABLE MC INTOSH BLADE SIZE 3

RIA53450 - DISPOSABLE MC INTOSH BLADE SIZE 4



### MILLER DISPOSABLE BLADES

Disposable fiber optic Miller laryngoscope blade, made in plastic material. Different sizes available.

Size 0: 75x8(h)

Size 1: 100x8(h)

Size 2: 132x9(h)

Size 3: 172x12(h)

RIA53432 - DISPOSABLE MILLER BLADE SIZE 0

RIA53434 - DISPOSABLE MILLER BLADE SIZE 1

RIA53436 - DISPOSABLE MILLER BLADE SIZE 2

RIA53438 - DISPOSABLE MILLER BLADE SIZE 3



### 3 FIBER OPTIC DISPOSABLE BLADES SET - RIA53410

Fiber optic laryngoscope complete with 3 Mc Intosh blades (2-3-4). Made in plastic material. Autoclavable handle. Supplied with a practical roll case.



### 6 FIBER OPTIC DISPOSABLE BLADES SET - RIA53400

Fiber optic laryngoscope complete with 6 blades (3 Mc Intosh, 3 Miller) and handle. Supplied with a practical case. Dimensions: about 28x4,5x25 cm.

## LARYNGOSCOPE ACCESSORIES



### CONVENT. LIGHT LARYNG. HANDLE

Conventional light laryngoscope handle made in stainless steel.  
Dimensions: Ø 30x148 mm.

RIA53030 - CL L. HANDLE  
RIA53182 - CL L. BLADE BULB



### FIBER OPTIC LARYNG. HANDLE

Fiber optic laryngoscope handle made in stainless steel.  
Dimensions: Ø 30x152 mm.

RIA53232 - FO L. HANDLE  
RIA53380 - FO L. BLADE BULB



### SMALL CONVENT. LIGHT LAR. HANDLE

Small conventional light laryngoscope handle made in stainless steel.  
Dimensions: Ø 20x155 mm.

RIA53130 - SMALL CL L. HANDLE  
RIA53180 - CL L. SMALL BLADE BULB



### SMALL FIBER OPTIC LARYNG. HANDLE

Small fiber optic laryngoscope handle made in stainless steel.  
Dimensions: Ø 20x152 mm.

RIA53330 - SMALL FO L. HANDLE  
RIA53382 - FO L. SMALL BLADE BULB



### DISPOSABLE FIBER OPTIC L. HANDLE

Fiber optic disposable laryngoscope handle, autoclavable.  
Dimensions: Ø 30x154 mm.

RIA53454 - DISPOSABLE FO L. HANDLE  
RIA53458 - FO DISPOSABLE L. BLADE BULB



### ROLL-CASE FOR LARYNGOSCHOPE

Roll-case for laryngoschope - red. 5 inner compartments. Made in washable cloth with inner stuffing.  
Dimensions : about 19x34 cm.

BOR84041 - ROLL CASE FOR LARYNGOSCHOPE

## AIRWAYS AND CATHETER MOUNT



### GUEDEL AIRWAY NO. 000

Pink - 40 mm  
Ref.: 127141001004



### GUEDEL AIRWAY NO. 00

Blue - 50 mm  
Ref.: 127141001006



### GUEDEL AIRWAY NO. 0

Black - 60 mm  
Ref.: 127141001008



### GUEDEL AIRWAY NO. 1

White - 70 mm  
Ref.: 127141001010



### GUEDEL AIRWAY NO. 2

Green - 80 mm  
Ref.: 127141001012



### GUEDEL AIRWAY NO. 3

Yellow - 90 mm  
Ref.: 127141001014



### GUEDEL AIRWAY NO. 4

Red - 100 mm  
Ref.: 127141001016



### GUEDEL AIRWAY NO. 5

Light blue - 110 mm  
Ref.: 127141001018

RIA50290 - SET OF 5 GUEDEL AIRWAYS (0/1/2/3/4)



### WENDL AIRWAY

Wendl nasopharyngeal airway, with ending pipe ring. Supplied singularly in sterile package. Different sizes available.

185200CH26 - WENDL AIRWAY CH 26-8,7 mm  
185200CH28 - WENDL AIRWAY CH 28-9,3 mm  
185200CH30 - WENDL AIRWAY CH 30-10 mm



### CATHETER MOUNT - 1515TG

Catheter Mount, lightweight, transparent, disposable and at low compliance. Supplied with perfect capacity connectors. Extremely flexible. Allows up to 360° angulation without any modification to inside volume, crushproof. Connectors in compliance to ISO norms. Dimensions: 20 cm.

## REANIMATION DEVICES

---



### STERILE ENDOTRACHEAL TUBES

---

Cuffed and uncuffed endotracheal tubes in plastic material. Unidirectional valve, standard connector in accordance to ISO norms. Reference notches, x-ray opaque, non-traumatic tip. Supplied in sterile package. Different sizes available.

---

11814102003	- CH12 3,0MM
11814102004	- CH14 3,5MM
11814102005	- CH16 4,0MM
11814102006	- CH18 4,5MM
11814102007	- CH20 5,0MM
11814102008	- CH22 5,5MM
11814102009	- CH24 6,0MM
11814102010	- CH26 6,5MM
11814102011	- CH28 7,0MM
11814102012	- CH30 7,5MM
11814102013	- CH32 8,0MM
11814102014	- CH34 8,5MM

---



### THOMAS HOLDER

---

Safety endotracheal tubes blocking system. Screw system, integrating a no-lockjaw surface which protects teeth and gums and avoids endotracheal tube crushing. It avoids use of plaster.

---

05-1650	- ADULT THOMAS HOLDER
05-1651	- PEDIATRIC THOMAS HOLDER

---



### HELICOIDAL MOUTH OPENER - RIA50302

---

Helicoidal mouth opener, manufactured in plastic material (atoxic polymer) completely autoclavable.

## REANIMATION DEVICES



### LIFE KEY MASK KEY-RING

Key-ring with mouth to mouth ventilation mask. Elastics for ear fixing. One way valve for a safe ventilation. Singulary packaged . Non-sterile.

RIA50247 - RED LIFE KEY MASK KEY-RING  
RIA50248 - BLUE LIFE KEY MASK KEY-RING  
FS-102 - LIFE KEY MASK



### QUICKTRACH

Urgent coniotomy sterile set quicktrach type. Complete with 1 plastic canula with intubation internal-blade stylet with diameter of 4 mm for adults and 2 mm for child. Complete with ISO joint for the connection to resuscitation bag. Quicktrach features a sharp conical-tip intubation stylet that avoids windpipe incision using classic scalpel and a "lock" that avoids to penetrate too much in depth. 2 lateral wings allowing airways safe fixing.

VBM30-04-0021 - CHILD QUICKTRACH  
VBM30-04-0041 - ADULT QUICKTRACH



### REMASK VENTILATION MASK - RIA50104

Remask transparent ventilation Mask, with pneumatic border. Mask adjustable to all kind of resuscitator bag. Disposable one-way safety valve. Blocking elastic. Connecting juncture to oxygen source. Standard connector 22mm. Supplied in a practical yellow case. CE certified in class IIa.



### INTUBATION STYLETS

Flexible and coated intubation stylets. Different sizes available.

RIA52213 - INTUBATION STYLET FR6 - 2,0 MM  
RIA52214 - INTUBATION STYLET FR10 - 3,3 MM  
RIA52215 - INTUBATION STYLET FR14 - 4,7 MM



## INTUBATION SETS - CONVENTIONAL LIGHT



### INTUBATION BAG COMPLETE - RIA50820

Intubation bag complete of: 1 conventional light laryngoscope with 3 blades (1-2-3), 1 Magill forceps, 1 Pean forceps, 1 scissors, 5 endotracheal tubes, 3 Guedel airways, 1 syringe, 1 mouth opener, 1 plaster roll, 2 Wendl airways, 1 stylet.

Dimensions: about 32x22x7,5(h) cm.

Weight: 1,2 kg.



### PEDIATRIC INTUBATION BAG COMPLETE - RIA50816

Intubation bag complete of: 1 pediatric conventional light laryngoscope with 3 blades (0-1-2), 1 Magill forceps, 1 Pean forceps, 1 scissors, 4 endotracheal tubes, 3 Guedel airways, 1 syringe, 1 mouth opener, 1 plaster roll, 2 Wendl airways, 1 stylet.

Dimensions: about 32x22x7,5(h) cm.

Weight: 1 kg.



### INTUBATION SET IN POUCH - VAL44410

Complete intubation set in pouch with transparent frontal, featuring: 1 laryngoscope pediatric handle, 1 C.L. Mc Intosh blade size 0, 1 C.L. Miller blade size 1, 1 C.L. Miller blade size 0, 1 Magill forceps, size 16, 1 Pean forceps, 4 endotracheal tubes, 4 Guedel airways sizes 000/00/0/1, 1 intubation stylet 1 plaster roll.

Dimensions: about 27x12x9(h) cm.

Weight: 600 g.



### PEDIATRIC INTUBATION SET IN POUCH - VAL44390

Complete intubation set in pouch with transparent frontal, featuring: 1 laryngoscope with 3 conventional light blades, 1 Magill forceps, 1 Pean forceps, 1 scissors, 5 endotracheal tubes, 3 Guedel airways, 1 syringe, 1 mouth opener, 1 plaster roll, 2 Wendl airways, 1 Intubation stylet.

Dimensions: about 27x12x9h cm.

Weight: 1kg.

## INTUBATION SETS - FIBER OPTIC

---

### INTUBATION BAG COMPLETE - RIA50822

---

Intubation bag complete of: 1 fiber optic laryngoscope with 3 blades (1-2-3), 1 Magill forceps, 1 Pean forceps, 1 scissors, 5 endotracheal tubes, 3 Guedel airways, 1 syringe, 1 mouth opener, 1 plaster roll, 2 Wendl airways, 1 stylet.

Dimensions: about 32x22x7,5(h) cm.

Weight: 1,3 kg.

### PEDIATRIC INTUBATION BAG COMPLETE - RIA50825

---

Intubation bag complete of: 1 fiber optic laryngoscope with 3 blades (0-1-2), 1 Magill forceps, 1 Pean forceps, 1 scissors, 4 endotracheal tubes, 3 Guedel airways, 1 syringe, 1 mouth opener, 1 plaster roll, 2 Wendl airways, 1 stylet.

Dimensions: about 32x22x7,5(h) cm.

Weight: 1 kg.

### DISPOSABLE INTUBATION BAG COMPLETE- RIA50826

---

Intubation bag complete of: 1 fiber optic laryngoscope with 6 blades, 1 Magill forceps, 1 Pean forceps, 1 scissors, 5 endotracheal tubes, 3 Guedel airways, 1 syringe, 1 mouth opener, 1 plaster roll, 2 Wendl airways, 1 stylet.

Dimensions: about 32x22x7,5(h) cm.

Weight: 1kg.





## SURGICAL MEDICATION AND INSTRUMENTS

---

We distribute surgical medication and instruments that are subjected to quality controls before being placed on the market.

## SPECIFIC MEDICATION KITS

---

We create specific kits to deal with emergencies of various entities and nature, for domestic, non-professional and professional use.

## FIRST AID CASES

---

We provide first aid boxes that meet the requirements of workplace safety.

## GAUZE COMPRESSES



### STERILE GAUZE COMPRESSES 2 PCS

Sterile medication gauze compresses in an easy-to open package, two 8-ply 10x10 compresses, pure cotton tit.12/8. Indicated for emergency and first aid activities. 2 pcs package. Ref.: MED96001.



### NON WOVEN STERILE GAUZE COMPRESSES 2 PCS

Sterile non woven compresses. Absorbing, non-adherent and resistant, two 8-ply. 10x10 compresses. 2 pcs package. Ref.: MED96101.



### STERILE GAUZE COMPRESS

Sterile medication gauze compresses in an easy-to open package, 1 18x40 cm compress, pure cotton tit.12/8. Indicated for emergency and first aid activities. Ref.: MED96008.



### NON WOVEN SELF-ADHESIVE

Non woven self-adhesive roll for medication fixing on all body surfaces. Dimensions: 10 m x15 cm. Ref.: 35CER334.



### 1 KG NON-STERILE GAUZE COMPRESS

Non-sterile medication gauze, pure cotton. 1 kg package. Different sizes available.

61610 - 1 kg GAUZE 10x10 cm  
61645 - 1 kg GAUZE 15x15 cm  
61670 - 1 kg GAUZE 20x20 cm  
61710 - 1 kg GAUZE 30x30 cm




### NON-STERILE GAUZE COMPRESSES – 100PCS

12-layer, non-sterile gauze compresses made of pure cotton. Package: 100 pcs. Different sizes available.

61050 - 100 pcs 5x5 cm  
61051 - 100 pcs 7,5x7,5 cm  
61052 - 100 pcs 10x10 cm



## ELASTIC TUBULAR NET, BANDAGES AND COTTON



### ELASTIC TUBULAR NET

Elastic tubular net for medications on every body's part. Easy and quick application.

75412 - SIZE 1 - 50 m  
75414 - SIZE 2 - 50 m  
75416 - SIZE 3 - 50 m  
75418 - SIZE 4 - 50 m  
75420 - SIZE 5 - 50 m  
75422 - SIZE 6 - 25 m



### ELASTIC TUBULAR NET

Elastic tubular net for medications. Easy and quick application. Length: about 40cm  
The size is indicated for hand, wrist, foot, ankle.

75466 - ELASTIC TUBULAR NET MEDIUM



### ELASTIC MEDICATION BANDAGE

Elastic medication bandage. Singularly packaged. Dimensions: about 4 m x 10 cm

MED96280 - BENDA ELASTICA



### EDGED BANDAGE

OB flexible medication bandage. Singularly packaged in sealed wrapping. Different sizes available.

FS-3570 - BANDAGE 5 cm  
FS-3572 - BANDAGE 10 cm  
FS-3554 - BANDAGE 8 cm



### COHESIVE BANDAGE

Self-fixing elastic bandage, lightweight and soft. Self adherent thanks to its atural rubber latex micro particles. Does not adhere to the skin and it is highly transpiring. Different sizes available, lenght 4 m.

370 - BENDA COESIVA 4 cm  
371 - BENDA COESIVA 6 cm  
372 - BENDA COESIVA 8 cm  
373 - BENDA COESIVA 10 cm



### COTTON

Absorbent cotton wool 20 or 50 g package. Not a medical device.

COT104 - COTTON 50 g  
COT279 - COTTON 20 g



## SOFT NEXT COHESIVE BANDAGE AND DETECTABLE DRESSINGS

### SOFT NEXT COHESIVE BANDAGE

Ideal for covering wounds especially at critical points, such as hands fingers. It is used directly on the wound. It is mouldable, autostick, elastic, water resistant and latex free. It absorbs the blood and does not stick to the skin, the hair and the wound itself. Skin or blue colour, different sizes available.

CER050 - SOFT NEXT 3x450 cm SKIN  
CER051 - SOFT NEXT 3x450 cm BLUE  
CER059 - SOFT NEXT 6x100 cm SKIN  
CER064 - SOFT NEXT 6x100 cm BLUE  
CER063 - SOFT NEXT 6x450 cm SKIN  
CER061 - SOFT NEXT 6x450 cm BLUE



### DETECTABLE BLUE DRESSINGS - CER049

Pack of 100 detectable sterile blue dressings, 5 assorted sizes. The metal band fitted in the internal part of the dressing is easy to find by the metal detectors. Contains: 20 pcs 9x72 mm, 40 pcs 30x72 mm, 32 pcs 55x18 mm, 4 pcs butterfly, 4 pcs 72x50 mm.



### EXTRA LONG 50PCS BLUE DETECTABLE DRESSING - CER143

The dressing length helps you to wrap your finger more than once, protecting the eventual cut or graze and avoiding accidental detachments, as it may happen for the smaller dressings. Dimensions: 18x2 cm. Pack of 50pcs.



## TAPE, SKIN SAVER AND STRIP PLASTERS

### NEURO-MUSCULAR KINESIOLOGY TAPE

Neuro-muscular tape, Thera-Band® Kinesiology Tape eliminates the risks of incorrect application thanks to the exclusive XactStretch™ Technology: our guide indicators allow you to always apply the tape with the correct tension. Relieves pain. Provides support for muscles and joints. Assists the treatment process of the edema. 98% Cotton, 2% Sapdex. Latex Free. Dimensions: 5 cm x 5 m.

305464 - KINESIOLOGY TAPE BLACK/GREY  
305465 - KINESIOLOGY TAPE WHITE/BALCK  
305467 - KINESIOLOGY TAPE BLUE  
305468 - KINESIOLOGY TAPE PINK/WHITE  
305469 - KINESIOLOGY TAPE RED/BLACK  
305470 - KINESIOLOGY TAPE GREEN/YELLOW



### TAPE PLASTER

Sport tape in white viscose with zinc oxide adhesive, suitable for rigid bandaging. Grainy fray-proof borders for an easy snatch. Available in two dimensions for a length of 14 m.

35CER692 - TAPE 3,8 cm  
35CER694 - TAPE 5 cm



### ELASTICATED STRIP PLASTER

Medicated strip plaster in elasticated cotton material, with transpiring gauze pad. Available in two dimensions for a width of 6 cm.

70051 - STRIP PLASTER 1 m  
70054 - STRIP PLASTER 5 m



### SKIN SAVER - MED97049

Thin and light layer made of polyurethane foam. Apply before bandaging. Dimensions: cm 7x27,5 m.



### POLYETHYLENE STRIP PLASTER - 70104

Skin coloured polyethylene medicated strip plaster. Waterproof, elastic, transpiring and hypoallergenic. High absorbing viscous gauze pad, does not stick to the wound. Dimensions: 5 m x6 cm. Singulary packaged.



## PLASTERS AND PLASTER ROLL



### NON WOVEN PLASTER ROLL

Non woven plaster roll.  
Available in two different  
sizes with a lenght of 9,14  
m.

00181 - ROLL 2,5 cm  
00183 - ROLL 5 cm



### NON WOVEN PLASTER

Adhesive, sterile medication  
with central absorbing pad  
non adhering to wound  
Hypoallergenic transpiring  
adhesive. Available in  
different sizes with a lenght  
of 10 cm.

35CER342 - PLASTER 8 cm  
35CER346 - PLASTER 15 cm  
35CER348 - PLASTER 20 cm



### SILK PLASTER ROLL

Adhesive silk plaster roll  
Available in two different  
sizes with a lenght of 9,14  
m.

35CER196 - SILK ROLL 2,5 cm  
00193 - SILK ROLL 5 cm  
02023173100000 - SOFFIX  
SILK ROLL 5 cm



### 20 PLASTERS - 35CER212

20 medicated plasters with  
gauze pad on adhesive  
support. Dimensions:  
about 19x76 mm  
Singlarily packaged.



### CLOTH PLASTER ROLL

Adhesive cloth plaster roll.  
Available in two different  
sizes with a lenght of 5 m.

02023132000000 -  
CLOTH PLASTER 2,5 cm  
02023133000000 -  
CLOTH PLASTER 5 cm



### 40 DIFFERENT SIZED PLASTERS - 35CER216

40 different sized  
medicated plasters with  
gauze pad on adhesive  
support. Various sizes.  
Singlarily packaged.

## PLASTERS


---



### OMNISTRIP SUTURE PLASTER - 540684

---

Sterile suture plaster skin colour, coated with an hypoallergenic adhesive layer. Suitable for small wound non-traumatic closing. 10 sterile strips. Dimensions: about 6x101 mm.



### ADHESIVE STERILE EYE BANDAGE - 35CER314

---

Sterile medication ocular compress with a high absorbent power on adhesive support. Hypoallergenic and transpiring. For secure and comfortable application on the eye.



### PLASTER FOR LINE-IN CANNULA NEEDLE - 35CER400

---

Non-woven plaster used to fasten and secure cannula needles. One package, 50 pcs. Size: 9x6 cm.

## IMMOBILIZATION AND MEDICATION SHEET

---



### TRIANGULAR IMMOBILIZATION SHEET

---

Triangular immobilization sheet, singularly packaged.

Available in cloth: 90x90x129 cm.

Available in non-woven: 90x90x127 cm.

---

FS-140 - CLOTH MATERIAL

FS-143 - NON-WOVEN MATERIAL

---



### STERILE MEDICATION CLOTH

---

Sterile medication cloth, ready to use for physiological treatment of wounds and burns. Hypoallergenic, transpiring with high. Absorbing capability. Allows transpiration of wound or burn and absorbs without sticking. Provided in a single sterile wrap.

---

62050 - STERIL CLOTH 40x60 cm

FS-148 - STERIL CLOTH 60x80 cm

62052 - STERIL CLOTH 80x120 cm

---



### HEMOSTATIC BANDAGE PAD - 1002001

---

Compressive bandaging for wounds and haemorrhages. It assures wound hemostasis thanks to the compression performed by elastic bandage on the hemostatic pad. The pad assures the hemostasis and moreover does not allow blood leakage. Supplied in sterile package. Dimensions: about 300x10 cm.



## MEDICATION PACKAGES

---



### MEDICAL PACKAGE D.I.N.

---

Sterile dressing package according to norm DIN 13151. Ready to use bandage and gauze. Hypoallergenic and with great transpiration power. Highly absorbent. Sterile packaging, immediately ready to use. Bandage length: 4 m. Different sizes available.

- 
- 62000 - DIN MEDICATION PACKAGE S - 6 cm
  - 62001 - DIN MEDICATION PACKAGE M - 8 cm
  - 62002 - DIN MEDICATION PACKAGE L - 10 cm
- 



### ALUMINIUM MEDICATION PACKAGES, ALUMINIUM COMPRESS AND ALUMINIUM CLOTHS

---

Sterile dressing for wounds and burns: aluminized medication package. Hypoallergenic, transpiring, anti-bacteria, highly absorbent and incohesive. Composed of a layer of aluminum and a layer of gauze with high absorbency that greatly reduces the infection. The aluminum film, in direct contact with the wound or with the burn, allows the transudate or blood to pass through with the aluminum layer not adhering to the wound. Promotes epithelialization and scar tissue formation. Thanks to the great absorbent power, but not adherent, it facilitates the change of the dressing.

- 
- 63050 - ALUMINIUM MEDICATION PACKAGE S
  - 63051 - ALUMINIUM MEDICATION PACKAGE M
  - 63052 - ALUMINIUM MEDICATION PACKAGE G
  - 63003 - ALUMINIUM COMPRESS 10x10 cm
  - 63100 - ALUMINIUM CLOTH 40x60cm
  - 63101 - ALUMINIUM CLOTH 60x80 cm
-

## DISPOSABLE GLOVES IN LATEX

---

### LATEX TALCUM POWDER

---



Disposable non-sterile gloves made of pure latex from natural rubber, lubricated with corn starch. Ambidexter with snatch-free wristband. High bio-compatibility. Food compatible (oil and grease excluded). 100 pcs package.

42576 - SIZE S 100 pcs  
42577 - SIZE M 100 pcs  
42578 - SIZE L 100 pcs

---

### LATEX POWDER FREE

---



Disposable latex gloves (natural rubber), powder-free with a double internal and external chlorination process. This process guarantees the absence of powder and reduces considerably protein elements. Not sterile. White, versatile ambidextrous with anti-rip wrist protection. Suitable for contact with food (oil and grease excluded). 100 pcs package.

39920 - SIZE S 100 pcs  
39921 - SIZE M 100 pcs  
39922 - SIZE L 100 pcs  
39923 - SIZE XL 100 pcs

---

### STERILE LATEX

---



Chirurgical sterile gloves made of pure latex from natural rubber, with powder with anatomic shape. Enlarged backside, lubricated with corn starch and provided with ring beam. Singularly packaged with double wrapping made of medical paper 1 pair.

CLASSIC65 - SIZE 6,5  
CLASSIC70 - SIZE 7,0  
CLASSIC75 - SIZE 7,5  
CLASSIC85 - SIZE 8,5  
40400DM - GLOVES BOXES  
WALL BRACKET

---

### STERILE LATEX SIZE L

---



1 pair of disposable sterile latex gloves. Singularly packaged. Size: L.

01BML0S0G - STERILE  
LATEX GLOVES SIZE L

---

## DISPOSABLE GLOVES IN NITRILE AND VYNIL



### NITRILE

Disposable non-sterile nitrile gloves internally oiled with bio-absorbing maize amid. "Latex Free" preserves from possible allergic reactions due to the natural rubber latex proteins. Ambidexter with snatch-free wristband. Micro-coarse surface granting the utmost sensibility. Highly resistant to chemicals and detergents, concentrated soaps and hydrocarbon. 100 pcs package.


39961S - SIZE S 100 pcs  
39961M - SIZE M 100 pcs  
39961L - SIZE L 100 pcs  
39961XL - SIZE XL 100 pcs



### NITRILE POWDER FREE 30 CM

Disposable non-sterile nitril gloves powder free, thicker than usual, for patients' examination. Hypoallergenic: latex free, non-allergenic. Blue, ambidextrous with snatch-free wristband, 30cm long over the arm. Thanks to a thicker made, they are suitable for jobs where a better protection and resilience are required. Food compatible. 100 pcs package.


39021 - SIZE S 100 pcs  
39022 - SIZE M 100 pcs  
39023 - SIZE L 100 pcs  
39024 - SIZE XL 100 pcs



### NITRILE POWDER FREE BLACK

Disposable non-sterile nitrile gloves powder free. Hypoallergenic: latex free, non-allergenic. Black, ambidextrous with snatch-free wristband and reinforced hems. The particular formula with which they are made, gives them an exceptional sensibility and softness and at the same time great flexibility. Micro-coarse surface granting the utmost sensibility. Food compatible. 100 pcs package.

42573S - SIZE S 100 pcs  
42573M - SIZE M 100 pcs  
42573L - SIZE L 100 pcs  
42573XL - SIZE XL 100 pcs



### VINILE POWDER FREE

Disposable non-sterile vinyl gloves powder free, non-adherent, coated with a waterproof polyurethane film. Hypoallergenic: latex free, non-allergenic natural rubber. Transparent, ambidextrous with snatch-free wristband, internal/external anti-glare matt finishing. Food-compatible. 100 pcs package.

39056 - SIZE S 100 pcs  
39057 - SIZE M 100 pcs  
39058 - SIZE L 100 pcs  
39059 - SIZE XL 100 pcs

## TOURNIQUET

---

### LATEX FREE RED - MED98210

---



Quick release latex free tourniquet, one hand operated. Easy to use, ensuring an excellent hemostasis without skin and tissues laceration. Quick opening and progressive/adjustable haemostasis reduction. Red colour.

### LATEX FREE BLUE - MED98212

---



Quick release latex free tourniquet, one hand operated. Easy to use, ensuring an excellent hemostasis without skin and tissues laceration. Quick opening and progressive/adjustable haemostasis reduction. Blue colour.

### LATEX FREE BABY - MED98208

---



Quick release latex free tourniquet, one hand operated. Easy to use, ensuring an excellent hemostasis without skin and tissues laceration. Quick opening and progressive/adjustable haemostasis reduction. Fantasy colour.

### OB HEMOSTATIC BAND - MED98230

---



Tourniquet OB type Martin Wide surface mm 75 x 1000 mm for an excellent haemostasis. Suitable for huge wounds and amputations. Supplied with regulation non-traumatic buttons. Latex free.

### FLAT TOURNIQUET - LAC190

---



Flat tourniquet. Singularly packaged. Latex free.

### CAT COMBAT APPLICATION TOURNIQUET - CAT-007

---



Original Cat Combat Application Tourniquet haemostatic tourniquet, in nylon material with micro metrics regulation.

## NEEDLES BOX



### MULTI-SAFE MINI 300

Container suitable for needles and sharps. With safety locking system. Made of polypropylene. Volume: 0,3 L. Dimensions: 7,8x4,5x16 cm.

MED98609 - MULTI-SAFE MINI 300



### MULTI-SAFE SANI 200

Container suitable for needles and sharps. With safety locking system. Made of polypropylene. Volume: 0,2 L. Dimensions: 7,3x4,5x13 cm.

MED98613 - MULTI-SAFE SANI 200



### MULTI-SAFE VARIO 2300

Container suitable for needles and sharps. With safety locking system. Made of polypropylene. Volume: 2,3 l. Dimensions: 16,1x10,6x24,4 cm.

MED98621 - MULTI-SAFE VARIO



### RIGIBOX EUROMATIC

Waste and needle box, specific for ambulances. With safety closure. Made of eco-friendly polypropylene with a red cap. Volume: 7 l. Size: about 17x17x 39 cm.

MED98631 - RIGIBOX EUROMATIC  
MED98623 - WALL HOLDER FOR EUROMATIC



### KONTAMED

Trash and needles box suitable for ambulance. Wall bracket for wall fixing. Complete of safety closure. Manufactured in yellow high resistant plastic material. Dimensions: about 170x75x225h mm.

2009067 - KONTAMED



### KONTAMED MINI

Small trash and needles box suitable for emergency bags and cases. Complete of safety closure. Manufactured in yellow high resistant plastic material. Dimensions: 45x115h mm.

2009068 - KONTAMED MINI



## SYRINGE AND NEEDLES

### CANNULA NEEDLE



2 FEB ways intravenous catheter needle with pliable wings injection valve. Best convenient use, utmost fixing to the skin. Sterile, disposable, non-toxic, x-ray opacity, nonpyrogenic.

50AC2S14 - ORAN. 14-45 mm  
50AC2S16 - GREY 16-45 mm  
50AC2S18 - GREEN 18-45 mm  
50AC2S20 - PINK 20-32 mm  
50AC2S22 - L. BLUE 20-25 mm  
50AC2S24 - YELL. 24-23 mm  
35CER400 - PLASTER FOR  
LINE-IN NEEDLE 50 PCS

### BUTTERFLY NEEDLES



Sterile epicranial butterfly needle, for intravenous infusion. Indicated for patients with capillary fragility, for children, and in the case a thin atraumatic needle is needed. Easily flexible.

90185 - GREEN 21G  
90187 - LIGHT BLUE 23 G

### CUP CANNULA NEEDLE - 90227



Disposable sterile protective cap for needle canula.

### 3-WAY TAP - 90232



3-way tap with spinning luer lock.

### STERILE SYRINGE



Sterile disposable syringes, singularly packaged.  
1 Package, 100 pcs.

90152 - SYRINGE 2,5 cc  
90153 - SYRINGE 5 cc  
90154 - SYRINGE 10 cc  
90155 - SYRINGE 20 cc

### FLUIDS ADMINISTRATION SET - 01D193E2N



Infusion set for administration of solutions. (Not for blood or blood products). Sterile infusion set for delivery of liquids and fluids in intravenous therapy and medical application. With an antibacterial drip chambers, lines and tubing with roller clamp. Disposable with protective end caps and filter.

## PROTECTION MASKS AND GLASSES

---



### PROTECTION MASK FFP2 WITH VALVE - MED98714

---

Disposable mask with FFP2 protection level. Rear fixing supports. 10 pcs package.



### PROTECTION MASK FFP3 - MED98712

---

Disposable mask with FFP3 protection level. Valve with protection filter. Rear fixing supports. 10 pcs package.



### PROTECTION MASKS - 40550

---

Disposable protection mask. Back fixing. 100 pcs package.



### MASK WITH PROTECTION VISOR - 12MUS4A4VIAM

---


Mask with protection visor. Midway filter. 4 layers fabric. Anti-mist shield. To be used upon standard eyeglasses. Support elastics.



### PROTECTION GLASSES - MED98704

---

Eye protection glasses. Transparent lens and frames. Scratchproof lens. Lateral protection on the bars. May be used upon standard eye glasses.



### PROTECTION SET EYE AND AIRWAYS - MED98708

---

Protection set composed by mask and transparent visor. Complete protection of eyes and airways. Water repellent, antiacid.

## EYE WASHING KITS

---



### EMERGENCY EYE WASHING KIT - CPS258

---

Recommended to all companies and on transport motor vehicles for emergency eye washing and its subsequent medication. Green case with wall support featuring: 1 empty eye washing bottle 500 ml; 1 eye washing sterile saline solution; wash 500 ml; 2 sterile eye pads; 2 sterile gauze packages 18x40cm; 1 pair of lister scissors 14,5cm; 1 plaster roll 1mx2cm; 1 towelette package; 1 hand mirror; 1 instruction manual.

Dimensions: about 25,5x9x19h cm.

Weight: about 1,2 kg.



### EYE WASHING KIT

---

Eye washing kit featuring 2 disposable and sterile saline solution 500 ml bottles. Rapid eye pad. Thermoformed support with transparent frontal and powder free closing, equipped with small mirror, screws and nogs for wall mounting.

Dimensions: about 32x8x41 cm.

Weight: about 1,8 kg.

---

CPS202 – EYE WASHING KIT – EYE-SHAPED  
TAP003 - EYE WASH CAP WITH ANTI-DUST COVER  
SOL002 - PHYSIOLOGICAL SOLUTION 250ML - NOT INJECTION  
SOL004 - PHYSIOLOGICAL SOLUTION 500ML - NOT INJECTION

---

## DISINFECTANT DEVICES - UMONIUM

---

### UMONIUM MEDICAL WET WIPES

---

Ultra-fast acting full-spectrum disinfectant, UMONIUM38 Medical Tissues are wet wipes which clean and disinfect, killing germs after 5 seconds of coming into contact with them. They can be used to clean and disinfect fragile materials such as optical fibre, polycarbonate, latex etc... Cleaning and disinfection of non-invasive medical devices. The impregnating solution is validated for bactericidal, fungicidal, mycobactericidal (including B.K.), and viricidal action. Aldehyde-free. Non-toxic. Do not irritate the skin. Dermatologically tested (occlusive dressing for 48 hours). Non-corrosive (pH 7), therefore compatible with optical fibre, polycarbonate (incubators), acrylic, glass, PVC, HDPE, PET, neoprene, latex, silicone, paint, steel, 410 stainless steel. Wipes contained in the box do not dry out. Stability = 36 months. CE medical device.

### UMONIUM MEDICAL SPRAY BOTTLE

---

Ultra-fast acting full-spectrum disinfectant, UMONIUM38 Medical Spray kills all germs present on the treated surface in just one spray (including BK, MRSA, Rotavirus and Pseudomonas). Its performance is particularly guaranteed in cases in which many disinfectants are ineffectual. Thanks to its neutral pH, it poses no danger to the user or to materials. It is the essential disinfectant for protecting the health of medical staff and their patients. Disinfection of the surfaces of non-invasive medical devices. Disinfection of invasive medical devices. Powerful disinfectant; bactericidal, fungicidal, mycobactericidal (including B.K.), and viricidal. Non-corrosive; can be used on optical fibre, polycarbonate (incubators), acrylic, glass, Pyrex, PU, PVC, HDPE, PET, neoprene, latex, silicone, paint, 410 stainless steel. Fast-acting: it kills germs from the first minute that it comes into contact with them. No aldehyde or peroxide. No toxic vapours and not harmful if it comes into contact with the skin. Non-irritating. pH7 (neutral). Product stability = 36 months. CE medical device.



---

PF10748 - UMONIUM MEDICAL WET WIPES 10 PADS pcs  
PF10760 - DISPENSER WET WIPES UMONIUM MEDICAL  
100 pcs

---



---

PF12203 - UMONIUM MEDICAL SPRAY BOTTLE 250 ml  
PF12209 - UMONIUM MEDICAL SPRAY BOTTLE 1 l

---

## DISINFECTANT DEVICES



### SEPTALDEIDE FOR DISINFECTION AND STERILIZATION

Sterilizing solution containing glutaric surfactant aldehyde, with lemon essence. Efficient against spores and viruses including HIV, LAV, HBV and HCV. Pure application for sterilization of instruments and other objects, including thermolabile ones.

S431 - DISINFECTANT 1 LT  
S432 - DISINFECTANT 5 LT



### PHARMA G DISINFECTANT

Disinfectant and detergent based on an association between didecylmethylammonium chloride and orthophenylphenol. Can be nebulized with ULTRA LITE or in a one-liter bottle with trigger gun (available separately). Pleasant Lemon flavour. Surgical Medical Device.

Q641 - PHARMA G DISINFECTANT 1 LT  
Q642 - PHARMA G DISINFECTANT 5 LT



### 1 LT DISINFECTANT IMPREDIS - CE07

Disinfectant detergent for surfaces, surgical instruments and dental impressions. Active ingredients: combination of quaternary ammonium and chlorhexidine salts. Colourless, slightly orange scented. Alcohol-free and suitable for delicate surfaces non compatible with alcohol. Effective against viruses: HIV and HBV.



### PHARMA STERIL DISINFECTING WIPES - OX02

Disinfectant wipes in a canister, 100 pcs. Suitable for disinfection of surfaces and instruments. Contains quaternary ammonium salts. Effective against viruses: HIV, HBV, HCV and H1N1.



## DISINFECTANT DEVICES



### PHARMAIODIO DISINFECTANT

Povidone iodine 1000 ml bottle 10% P.M.C. iodine. Pharmaiodio is bactericide and virucidal and is suitable for pre-injection disinfection, pre-surgical preparation and intact skin antiseptic.

PF08110 - PHARMAIODIO 1 l  
PF08111 - PHARMAIODIO 500 ml  
PF08173 - PHARMAIODIO 125 ml



### HYDROGEN PEROXIDE

1 lt. hydrogen peroxide bottle. Hydrogen peroxide, 3% solution – 10% volume. Stabilized with proper treatment. Suitable for skin cleaning.

MED98508 - HYDROGEN PEROXIDE 250 ml  
MED98510 - HYDROGEN PEROXIDE 1 l



### NON-INJECTABLE PHYSIOLOGICAL SOLUTION

0.9% non-injectable sodium chloride physiological solution. Sterile and apyrogenic solution in a recyclable polyethylene container. For wounds, suitable in-between and post-surgery. For bladder catheters and oral cavity cleaning. Cleaning and rinsing of surgical instruments and equipment.

195B230541 - PHYSIOLOGICAL SOL. 250 ml  
195B230551 - PHYSIOLOGICAL SOLUTION 1 l



### PHARMASIL

Rapid and persistent disinfectant for non-wounded skin. Mixture of didecyldimethylammonium chloride and orthophenyl. It rapidly gets rid of unwanted pathogenic microorganisms. Efficient against gram+ and gram- acid resistant bacteria. Scented, no burn, no stains. Surgical medical device.

B300 - PHARMASIL 250 ml  
S801 - PHARMASIL 1 l

## DISINFECTANT DEVICES



### SEPTAMAN GEL

Surgical Medical Device. Gel colourless alcohol solution (Ethyl alcohol 70%), ready to use, suitable for rapid and effective disinfection of intact skin. Kills bacteria, fungus, micro-bacteria and viruses. Furthermore, inactivates HIV, HBV, HCV and Rotaviruses. Contains no allergenic substances.

PF32459 - SEPTAMAN GEL 100 ml

PF32411 - SEPTAMAN GEL 500 ml

AN0425058 - WALL HOLDER SEPT. 500 ml



### MEDICAL SPRAY DISINFECTANT - PF011

150 ml bactericide disinfectant. Bactericide and fungicide action, suitable for environmental disinfection treatments with automatic supply. This system allows the distribution of the product in the environment thus efficiently disinfecting floors, walls, surfaces and air. It allows the complete and rapid disinfection of waiting, emergency, consulting and surgical rooms besides ambulances and classroom. Products diffusion is about 4 minutes long.



### DISINFECTANT TOWELETTE - TOW223

Disinfectant towelette: each sachet contains 3 ml of Pharmaderm.  
Dimensions of sachet: 65x95 mm circa.  
Dimensions of towelette: about 140x190 mm.



### 4% AMMONIA WIPE - TOW013

Ammonia wipe in water solution 4%.  
10 pcs package.

## COLD/HOT THERAPY

---



### DISPOSABLE INSTANT ICE PACK - FCT2000

---

Disposable ice with controlled temperature (0°C for about 30 minutes). Its immediate use during traumas avoids haematoma and swelling formation. It stops little haemorrhages and considerably reduces the pain. It can be kept without limits at ambient temperature.

Made in PVC material. Dimensions 17x13,5 cm. CE medical device.



### DISPOSABLE NON-WOVEN ICE PACK - 02SF1853EASYICE

---

Disposable synthetic ice that keeps 0°C temperature for about 30 minutes.

Does not contain ammonium nitrate.

Made in non-woven material.

Dimensions: 13,5x17 cm.



### SPRAY ICE

---

Spray ice bottle ready to use. Particularly suitable for sport activities, small traumas and contusions. CE medical device.

---

SP200DISPO24 - SPRAY 200 ml

SP400DISPO24 - SPRAY 400 ml

---

## COLD/HOT THERAPY

---



### HOT/COLD GEL - DISPOGEL

---

Reusable gelatin hot/cold Dispo GEL  
It is a hot/cold reusable accumulator from outside source. Completely non-toxic, suitable for medical and domestic use: easily modeled, does not harden becoming ice and it is reusable. Cold therapy: keep in freezer for at least 2 hours. Take out of the freezer and put in a cloth and apply on the affected part. Hot therapy: place the bag in the microwave for at most 1 minute setting it to 650 Watt or dip the bag in boiling water for at most 5 minutes, remove and dry; in both cases insert the bag in a cloth and apply on the affected part.

Dimensions: about 14x24 cm.



### DISPOSABLE INSTANT HOT PACK - 002SC1853HOT

---

Disposable instantaneous hot pack with a controlled temperature of 60°degrees lasting about 25 minutes.

Dimensions: about 13,5x17 cm.

## COLD/HOT THERAPY

---



### WARMING INFUSION BAG - FCT2310

---

Warming infusion bag, suitable for 2 sacks  
Inner cloth for the insertion of an instantaneous or reusable warm compress.  
Made in blue cordura material. Inner covering made in insulating thermo-refractive material. Lateral transparent window suitable for examination.  
Anterior pocket suitable for flow control devices storage.  
Dimensions: about 23x13x28h cm.  
Weight: 450 g.



### FIREPROOF BLANKET

---

Fiberglass, fireproof blanket, EN 1869 compliant. Suitable for extinguishing small fires, especially for small oil fires.  
Supplied with a case equipped with a ring for wall mounting and velcro opening.

---


COP208 - FIREPROOF BLANKET 120x180 cm  
COP215 - FIREPROOF BLANKET 150x100 cm

---



## ISOTHERMIC BLANKETS

---



---

FCT4000 - ADULT ISOTHERMIC BLANKET  
FCT3010 - PEDIATRIC ISOTHERMIC BLANKET

---

### ISOTHERMIC BLANKET

---

Protection and rescue isothermic blanket  
Prevents patients' hypothermia whether traumatized, burned or post-operated  
Double function: warming and cooling  
Silver side, if positioned toward the body, isolates and reflects, preserving about 80% of body warmth. Thus protects patient from hypothermia, coldness, freezing and shock.

Golden side: placed toward patient body and silver one, reflects sun ray maintaining body freshness. Non sterile.

Adult dimensions: about 160x210 cm.

Pediatric dimensions: about 100x80 cm.



### BLIZZARD EMS BLANKET ORANGE- FCT3060

---

Designed specifically for the Emergency Services market the EMS is made from 2 layer Reflexcell™ and is widely used by all emergency services worldwide. The Blizzard EMS blanket works as an effective tool in the prevention and treatment of hypothermia. The Blizzard EMS blanket is a high performing blanket with self-adhesive closure that re-opens for easy access to casualties. The Blizzard EMS blanket is available vacuum packed for ease of storage: Small packaging. Easy to open and roll out.

Size (vacuum packed): 12x11x5 cm

Weight: 250 g

Size (unpacked): 2,3x1,96 meters (fully stretched) 156.8cm (working)

# ABSORBING TRANSFER SHEET WRASFER - 610007

## TRANSFER SHEET

### 150 KG AND 3,5 LITERS

Wrasfer sheet in addition to being highly resistant (150 kg maximum carrying weight) allows the absorbency of fluids and exudates up to 3,5 liters, which guarantees the optimum containment of possible infections. Disposable device of the risk class I, made of 3-layer construction: PE-foil, nonwoven and cellulose fluff.

It has been specifically studied for the transfer assistance of patients and it has been designed for use in paramedic services, hospitals, hospital emergency wards, nursing homes, funeral chambers and for remains transportation in special conditions.



### EASY TO USE

After its use the sheet can be easily eliminated and the carriers or stretchers are ready for the next use.

Sheet dimensions: 220x100 cm. Weight: 350g.

The application of the transfer sheet should be adequate to provide assistance to two/four patients.



## WASTE SACKS, BASINS AND SHEETS

---



### VOMIT ABSORBING BAG - V60020

---

Vomit bag “100% Protected Hands. Fitted with an absorbent pad with Alveol-Technology, it absorbs the vomit and allows its elimination: no dispersion. New MDPE liner: more resistant. Does not stick to the skin. Hypoallergenic. Dimensions: 25x25 cm. Packs containing: 20 pcs. Latex Free, single use only. CE Medical device.



### WHITE VOMIT AND WASTE SACK - SAC052

---

White waste sack with hermetic closing. It takes little space and can be stored inside emergency bags and rescue rucksacks. Dimensions: 30x50cm. 100 pieces package. Available singulary packaged with ref. XME0530.



### ROLLED SHEETS - 005 P LU

---

Sheet rolls in pure cellulose, suitable for examination tables Embossed surface, two layers. Average veil weight: 17 gr/mq. Pack of 6 pieces. Colour: white. Dimensions: 70m x 60 cm.



### DISPOSABLE PILLOWCASE - FCT2466

---

Disposable non woven plasticized pillowcase, non-allergenic, reinforced, absorbing and impermeable. Dimensions: 40x40 cm. 25 pieces package.



### PATIENTS BASIN - 390470

---

White patients basin Made in sterilizable polypropylene. Singulary packaged. Dimensions: 27x52,3 cm handle included. Weight: 300 g. Max carrying: >130 kg.



### URINAL

---

Plastic white urinal. Made in sterilizable polypropylene. Singulary packaged.

---

MED98784 - URINAL WITH COVER

---

## INSTRUMENTS - SCISSORS

---



### FIRST AID SCISSORS CUT CLOTHING

---

Stainless steel first aid scissors cut clothing. Autoclavable plastic handle. Suitable for first aid bags and cases. Dimensions: 18 cm.

---

STR110010 - BLACK SCISSOR  
STR110010/A - ORANGE SCISSOR  
STR110010/B - BLUE SCISSOR  
STR110010/G - YELLOW SCISSOR  
STR110010/R - RED SCISSOR  
STR110010/V - GREEN SCISSOR

---



### LISTER SCISSORS

---

First aid cut-clothing stainless steel scissors. Suitable for emergency bags and cases .

---

STR110290 - LISTER SCISSORS 14 cm  
STR110300 - LISTER SCISSORS 16 cm  
STR110310 - LISTER SCISSORS 18 cm  
STR110320 - LISTER SCISSORS 20 cm

---



### LITTAUER LIGATURE SCISSOR - STR110330

---

Ligature scissors. Made in stainless steel  
Dimensions: 14 cm.

## INSTRUMENTS - SCISSORS

---



### FIRST AID SCISSORS - FOR078

---

First aid scissors with plastic handle.  
Stainless steel blades.  
Dimensions: 10 cm.



### STRAIGHT SCISSORS WITH BLUNT

---

Straight medication scissors with blunt tips.  
Made in stainless steel.

---

STR110080 - SIZE 14 cm  
STR110090 - SIZE 16 cm

---



### STRAIGHT SCISSORS WITH SHARP TIPS

---

Straight medication scissors with sharp tips.  
Made in stainless steel.

---

STR110040 - SIZE 14 cm  
STR110050 - SIZE 16 cm

---



### STRAIGHT SCISSORS WITH ALTERNATE TIPS

---

Straight medication scissors with alternate tips, blunt/sharp. Made in stainless steel.

---

STR110060 - SIZE 14 cm  
STR110070 - SIZE 16 cm

---



### SCISSORS WITH SHARP TIPS

---

Medication scissors with bent alternate tips, blunt/sharp. Made in stainless steel.

---

STR110030 - SIZE 14 cm  
STR110032 - SIZE 16 cm

---



### BENT SCISSORS WITH ALTERNATE TIPS

---

Medication scissors with bent alternate tips, blunt/sharp. Made in stainless steel.

---

STR110020 - SIZE 14 cm  
STR110035 - SIZE 16 cm

---



## INSTRUMENTS - SCISSORS

---



### MULTIPURPOSE EMERGENCY SCISSORS - FOR299

---

Multipurpose Emergency scissors  
Ideal in cases of emergency to break crystals, cut clothes, cut safety belts, cut thick objects (knife function). Made of metal that is not subject to oxidation, without plastic parts. You can open it completely, dividing it into two pieces, to wash it more easily. Completely autoclavable. Supplied with cover and accessories for wall mounting.



### ROBIN SCISSORS - 124900

---

Multi-function emergency scissors, Robin original. Made in Solingen steel. It easily cuts fasten belts. Its handle allows car glasses crush. Supplied with wall support. Weight: 400 g.

---

124900 - ROBIN SCISSORS

124901 - WALL HOLDER FOR ROBIN SCISSOR

---



### STERILE DISPOSABLE SCALPEL

---

Disposable sterile scalpel, singularly packaged. Handle in abs material. Stainless steel blade. Blade protection in polythene, easy to remove.

---

90466 - DISPOSABLE SCALPEL SIZE 10

90467 - DISPOSABLE SCALPEL SIZE 11

90468 - DISPOSABLE SCALPEL SIZE 12

---

## INSTRUMENTS

---



### REMOVAL PENCIL - MAT022

---

Small metallic splinters eye tool remover. One side you have a NEODIMIO magnetic point, which enable you to withdraw all metallic splinters avoiding any eventual dangerous eye approach. On the other side you have an non-abrasive nylon loop which is used for non metallic splinters. Main body and support made of anodised aluminium and the neoprene seals guaranty a better pencil hold. Supplied in a practical box. CE medical device.



### PERCUSSION HAMMER - STR110270

---

Percussion hammer. Zinc plated handle. Head in pvc. Complete with needle. Dimensions: about 18x6x1,5 cm.



### DISPOSABLE TRICHOTOMY RAZOR - RSA225

---

Trichotomy razor suitable for hairs removal in the area of the skin supposed to surgical operation so to avoid possible infections. Stainless steel blade. Disposable non-sterile single blade.

## INSTRUMENTS - FORCEPS

---



### MAGILL FORCEPS

---

Intubation Magill forceps. Made in autoclavable and disinfectable stainless steel.

---

STR110230 - MAGILL FORCEPS 16 cm  
STR110240 - MAGILL FORCEPS 20 cm  
STR110250 - MAGILL FORCEPS 25 cm


---



### CLAMP PLASTIC FORCEPS - STR110450

---

Clamp plastic forceps.  
Dimensions: 12 cm.



### SURGICAL AND ANATOMIC TWEEZERS

---

Anatomic tweezer made in stainless steel.

---

STR110100 - ANATOMIC TWEEZER THIN 14 cm  
STR110110 - ANATOMIC TWEEZER 14 cm  
STR110120 - ANATOMIC TWEEZER 16 cm  
STR110130 - SURGICAL TWEEZER 14 cm  
STR110140 - SURGICAL TWEEZER 16 cm

---



### PEAN, KOCHER AND WIRE CONDUCTOR FORCEPS

---

Tweezer made in stainless steel.

---

STR110150 - PEAN FORCEPS 14 cm  
STR110160 - PEAN FORCEPS 16 cm  
STR110170 - KOCHER FORCEPS 14 cm  
STR110180 - KOCHER FORCEPS 16 cm  
STR110190 - MATHIEU WIRE CONDUCTOR  
FORCEPS 14 cm  
STR110200 - MATHIEU WIRE CONDUCTOR  
FORCEPS 16 cm  
STR110220 - WIRE CONDUCTOR FORCEPS  
14 cm  
STR110280 - MOSQUITO FORCEPS

---

## INSTRUMENTS - FORCEPS

---

### STERILE UMBILICAL FORCEPS - 420140

---

Disposable sterile umbilical forceps. Made in non-toxic PVC material. Singularly packaged. Dimensions: about 5,5 cm



### NYLON FORCEPS

---

Nylon forcep singularly packed  
Sterile version: 10 cm.  
Non-sterile version: 11 cm.



---

PIN110 - STERILE FORCEP  
FUS0590 - NON-STERILE  
FORCEP

---

### RING CUTTER FORCEPS - STR110500

---

Stainless steel ring cutter forceps. Dimensions: 18 cm.



### PLASTIC TICK-REMOVING TWEEZERS - STR110462

---

Suitable for animals. The pen-shaped structure allows to capture the body of the tick in order to extract it with a rotatory movement. The correct extraction prevents infections. Dimensions: 9 cm. Colours availability: random.



### STAINLESS STEEL TONGUE FORCEPS - STR110210

---

Stainless steel tongue forceps. Dimensions: 17 cm.



### PLASTIC TONGUE FORCEPS - PIN201

---

Polycarbonate tongue forceps. Dimensions: 16 cm.



## INTRAOSSEUS - B.I.G. BONE INJECTION GUN

### BONE INJECTION GUN

B.I.G.-Bone Injection Gun is a new automatic intraosseus infusion device providing an easy, safe and rapid access. B.I.G. is the effective alternative answer to unsuccessful intravenous access during emergencies, which dramatically increases the rescuer's ability to save lives.

### SAFETY AND READY TO USE

It avoids delay during intravenous injection activities for the infusion of medicaments in emergency operations thus ensuring the perfect solution. Infusion established in less than 1 minute. Users friendly and relatively painless penetration. Suitable for blood transfusions. Rapid infusion of fluids and drugs. Extreme safety – no direct contact with patient's blood.



### PEDIATRIC B.I.G.

Pediatric B.I.G.-Bone Injection Gun is available.



MED98850 - ADULT B.I.G.  
MED98854 - PEDIATRIC B.I.G.



## INFECTIVE SELF-PROTECTION KITS

---



### EMERGENCY KIT INFECTIVES VIRUS EBOLA - VAL44079

---

The set includes all the necessary devices to follow the procedures for the management of the cases of suspected Ebola Virus disease, according to the operative indicative guidelines written into the bulletin of Italian Ministry of Health. Tearproof pouch with transparent front containing: protection glasses, 2 filtering face respirators, 2 pairs on unsterile nitrile disposable gloves, powder free and latex free, 2 singularly packaged tissues for disinfection, spunbond non-woven fabric protection-suit. Dimensions: about 30x12x13 cm.



### INFECTIVE SELF-PROTECTION INDIVIDUAL SET - VAL44080

---

Indicated as standard equipment on rescue vehicles, featuring: 1 protection mask with filter FFP3, 1 pair of protection glasses, 1 protection overalls with cap, 1 pair of galoshes, 3 disinfectant wipes, 1 pair of latex gloves, 2 trash bags, 1 coloured pouch with transparent frontal. Dimensions: about 27x12x9(h) cm.



### INFECTIVE SELF-PROTECTION INDIVIDUAL DOUBLE SET - VAL44084

---

Indicated as standard equipment on rescue vehicles, featuring: 2 protection masks with filter FFP3, 2 pairs of protection glasses, 2 protection overalls with cap, 2 pairs of galoshes, 3 disinfectant wipes, 2 pairs of latex gloves, 5 trash bags, 1 coloured pouch with transparent frontal. Dimensions: about 25x25x8(h) cm.

## EMERGENCY KITS



### PLEUROTHORAX SET - 250.200

Closed System comprising: secretion bag 2,0 ltr. with non-return valve, stopcock on the base and three-way stopcock (fixed). Transfusion Tube Luer-Lock (unscrewable). Puncture Cannula Diameter 2 mm; puncture length 75 mm. 50-ml-syringe Luer-Lock. Sterile double wrapped; Peel-Pack.



### PERIPHERAL VEINS CANNULATION SET - VAL44060

Peripheral veins cannulation set in pouches with transparent frontal, containing: 1 haemostatic belt, 5 canula needles (16/17/18/20/22), 3 butterfly needles (21/23/25), 1 silky plaster roll 2,5 cm, 4 disinfectant wipes, 2 IV administration sets. Dimensions : about 27x12x9(h) cm



### VALVE KIT TYPE HEIMLICH - 220000066910

Non-return valve for thoracic drainage (Heimlich type). Valve body in flexible and transparent polyurethane with natural rubber sleeve in "duck's bill" shape acting as a non-return valve. This device is connected directly onto the patient's pleural or thoracic drain (adult, paediatric), or indirectly, using connector code 800.01 and an extension tube for neonates and infants. This product allows thoracic drainage to continue without risk to the patient even when transported, as the non-return valve means that the drain does not have to be clamped.



### DISPOSABLE SUTURE SET - 28827

The kit contains the following material: 1 non woven sheet + PE 50x60cm. 1 non woven sheet + PE with hole and cut, adhesive. 2 adhesive non woven stripes with support 3x20cm. 5 compresses of surgical gauze, pure cotton 20 F.U. type, 10x10cm with 8 layers. 4 tampons of surgical gauze, 2 latex gloves, 1 haemostatic forceps, 1 plastic anatomic forceps, 1 forceps with thin point in steel and plastic, 1 scissors, 1 needle-carrier forceps, 1 pack of antiseptic single dose disinfectant, 25ml chlorhexidine based gluconate + cetrimide. 1 waste bag.

## EMERGENCY KITS



### BABY FIRST AID KIT - VAL44077

Pediatric first aid kit complete with: 2 small bags for minor burns treatment, 4 2-pieces non-woven sterile compress, 1 plaster roll, 1 plaster strip, 2 glove pairs, 3 disinfectant tissues, 1 nylon non sterile tweezers, 1 cohesive elastic bandage, 1 first aid scissors, 2 ammonia towels, 1 suture plasters package, 2 disposable standard instant ice, 40 plasters. Dimensions: about 23x14x11 cm.



### MOUNTAIN RESCUE FIRST AID KIT- VAL44082

Mountain rescue first aid kit complete with: 1 emergency isothermal blanket, 2 small bags for minor burns treatment, 5 2-pieces non woven compresses, 1 sterile sheet, 1 bag of aluminized dressing, 1 striped medicated plaster, 2 pairs of non-sterile gloves, 4 disinfecting wipes, 1 small first aid scissors, 1 sterile nylon tweezers, 1 life saver wipe, 2 napkins, ammonia 4%, 1 pack of suture plasters, 1 pack of standard disposable ice, 40 assorted plasters, 1 instant hot compress, 250ml of non-injectable saline solution. Dimensions: about 23x14x11cm.



### TRAVEL KIT - VAL44081

Travel first aid kit complete with: 2 small bags for minor burns treatment, 2 2-pieces non-woven sterile compress, 1 plaster strip, 2 glove pairs in non-sterile package, 3 disinfectant tissues, 1 nylon non sterile tweezers, 1 cohesive elastic bandage, 1 first aid scissors, 2 ammonia towels, 1 suture plasters package, 1 disposable standard instant ice, 40 plasters. Dimensions: about 15x12x8,5 cm.



### PET KIT - VAL44076

First aid kit ideal during small emergencies involving pets, designed especially for dogs. Small handbag in washable PVC made of coloured waterproof and tearproof fabric with side Velcro strip for anchoring to backpacks or rescue bags. The handbag contains the following material: Plaster roll. Cohesive elastic dressing bandage. DIN medication package. 2x sterile gauze (2 pcs). 2x pairs of gloves. 2x disinfectant wipes. Small first aid scissors. Isothermal emergency blanket. Sterile single-use syringe. Nylon tweezers. Tick-removing tweezers. Dimensions: about 10x7x13(h) cm.

## EMERGENCY KITS

### HEMOCARE FIRST AID KIT - VAL44083

Homecare first aid kit complete with: 1 emergency isothermal blanket, 2 small bags for minor burns treatment, 5 2- pieces non-woven sterile compress, 1 medicated plaster strip, 3 glove pairs in non-sterile package, 4 disinfectant tissues, 1 nylon sterile tweezers, 40 plasters, 1 first aid scissors, 3 ammonia towels, 1 suture plasters package, 2 disposable standard instant ice, 250 ml non-injectable physiologic solution, 20 gm. hydrophilic cotton, 250 ml. Peroxide, 2 cohesive elastic bandages, 40 cm tubular medication net, 80 ml bottle for minor burns treatment, 1 digital thermometer.  
Dimensions: about 25,5x9x19,5(h) cm.

### DIN13164 PACKAGE - CPS685

Red nylon handbag. Content complies with what stated by norm DIN 13164 (mandatory on automotive vehicles in certain European Nations). 8 adhesive bandages 10x6cm. Plaster roll 5 m x 2,5 cm. 2 gauze bandages 6 cm x 4 m. 3 gauze bandages 8 cm x 4 m. 1 pair of bandage-cutting scissors. Package of 4 vinyl gloves. Non-woven fabric sheet 60x80 cm. 2 non-woven fabric sterile sheets 60x40 cm. 3 sterile swab bandages 80x100 mm. Sterile swab bandage 100x120 mm. 6 non-woven fabric sterile swabs 10x10 cm. 2 non-woven fabric triangular sheets 96x96x136 cm. Isothermal blanket.  
Dimensions: 25x15x7,5 cm. Weight: 0,45 kg.

### EURO DIN 13164 - CPS222

Polystyrene case in compliance with DIN 13164B norm. 8 adhesive bandages 10x6,1 cm. 1 adhesive plaster roll m 5x2,5 cm. 2 gauze bandages cm 6x4 m. 3 gauze bandages cm 8x4 m. 1 pair of bandage cutter scissors. 4 vinyl gloves. 1 sterile non woven sheet 60x80 cm. 2 sterile non woven sheets 60x40cm. 3 sterile tampon bandages 80x100 mm. 1 sterile tampon bandage 100x120 mm. 6 non woven sterile tampons 10x10cm. 2 non woven triangular sheets 96x96x136cm. 1 isothermal blanket. First aid instructions.  
Dimensions: about 26,3x17x8,3 cm.





## EMERGENCY KITS

---



### FIRST AID INDIVIDUAL KIT - VAL44078

---

Complete first aid and medication kit:  
Isothermic blanket, 3,5 ml Idrogel sachet,  
2 sterile non woven compresses, Sterile  
medication cloth 60x80 cm, 2 Alumed  
medication pack size M, Transparent plaster  
roll 2,5 cm, Cotton medicated plaster 10x6  
cm, 2 pairs of gloves in non-sterile package,  
2 disinfectant wipes, Small first aid scissors,  
Nylon sterile tweezer, Life saver, Ammonia  
4% wipe, First aid instructions.  
Dimensions: about 15x12x8,5 cm.  
Weight: about 200 g.



### CAR DIN 13164 SET INCREASED - VAL44362

---

Practical first aid bag for cars: 1 plaster  
roll, 8 strip plaster, 1 wound medication  
package L, 3 wound medication packages  
M, 2 immobilization triangular clothes, 1  
sterile medication cloth for burn, 2 sterile  
medication clothes for burns, 3 big sized  
medication bandage gauzes, 2 medium  
sized bandages, 6 sterile gauze compresses,  
1 angular medication scissors, 1 emergency  
isothermic blanket, 4 vinyl protection gloves,  
8 disinfectant wipes, 3,5 ml Idrogel burns  
sachet, 2 4% ammonia wipes for insect  
stings.  
Dimensions: about 27x12x9(h) cm.  
Weight : about 0,5 kg.



### FIRST AID KIT-MOTORBIKE DIN 13167 - 17008

---

First Aid Kit for motorbikes with content  
according to DIN 13167. 1 roll of adhesive  
plaster, 14-piece adhesive dressing  
assortment: 4 wound dressings, 2 fingertip  
dressings, 2 finger bandages, 2 plasters, 4  
plasters. 1 first aid dressing, sterile, 1 first  
aid dressing, sterile, 1 burn dressing, sterile,  
1 rescue blanket, 1 bandage scissors, 2  
pairs of disposable gloves, 2 wet cleansing  
wipes for unwounded skin, 1 first aid  
instructions, 1 list of contents.



## REPLANTATION KITS

---



### REPLANTATION SET - VAL40954

---

Replantation set in a tearproof impermeable equipped bag with 4 coloured pouches with handles transparent frontal, complete of:

1 sterile replacement/transport set for finger,  
1 sterile replacement/transport set for hand,  
1 sterile replacement/transport sack for arm,  
1 sterile replacement/transport sack for leg,  
Hemostatic belt/lace, 1 tourniquet type Martin,  
10 sterile gauze packs 10x10cm,  
1 plaster roll 2,5cm, 1 emergency pad bandage,  
1 Alumed cloth 40x60cm, 1 anatomic forceps,  
1 surgery forceps, 1 Pean forceps.

Dimensions: about 48x16x26(h) cm.

Weight: about 5,5 kg.

---

MED98820 - HAND REPLANTATION SET - 300x250 mm  
MED98824 - ARM REPLANTATION SET - 700x250 mm  
MED98828 - LEG REPLANTATION SET - 1100x400 mm  
MED98832 - FINGER REPLANTATION SET - 200x200 mm  
MED98810 - HAND REPLANTATION BAG - 300x200 mm  
MED98811 - FINGER REPLANTATION BAG - 200x200 mm  
MED98812 - ARM REPLANTATION BAG - 700x200 mm  
MED98813 - LEG REPLANTATION BAG 1100x400 mm  
VAL40954 - REPLANTATION SET

---

### REPLANTATION LIMBS SET CASE - VAL40950

---

Replantation limbs set case, made in shockproof orange ABS, featuring:

1 sterile replacement/transport set for finger,  
1 sterile replacement/transport set for hand,  
1 sterile replacement/transport sack for arm,  
1 sterile replacement/transport sack for leg,  
1 haemostatic lace/belt, 10 sterile gauze packs 10x10cm,  
1 plaster roll 2,5cm, 1 emergency pad bandage,  
1 Alumed cloth 40x60cm, 1 anatomic forceps,  
1 surgery forceps, 1 Pean forceps.

Dimensions: about 42x15x31(h) cm.

Weight: about 7,5 kg.





## WATER-JEL BURN DRESSING

---

Water-Jel Burn Dressings provide fast, one-step relief for emergency first aid when someone has been burned. The dressings stop burn progression and cool the skin, relieve pain and protect against airborne contamination. Water-Jel Burn Dressings are made with a scientifically formulated cooling gel and utilize a medical-grade nonwoven pad that is non-adherent, so they will not stick to burn injuries. The dressings are soaked in a water-based, water-soluble gel and sterilized to ensure a contamination-free atmosphere for emergency treatment of burns. With Water-Jel Burn Dressings, it takes no special training to give fast and effective emergency burn care (4-6 hours) for any type of burn injury. All Water-Jel Burn Dressings have a 5-year shelf life from date of manufacture.



### WATER-JEL BURN DRESSING

---

B0206 - BURN DRESSING 5x15 cm  
B0404 - BURN DRESSING 10x10 cm  
B0416 - BURN DRESSING 10x40 cm  
B0818 - BURN DRESSING 20x45 cm  
B0820 - HAND BURN DRESSING 20x55 cm  
B1216 - FACE BURN DRESSING 30x40 cm

---



### WATER-JEL IN GEL

---

GBJ140 - BURN JEL 4.0 gr  
GBJ80 - BURN JEL 80 ml

---



### WATER-JEL FIRE BLANKET

---

Fire blanket Water-Jel provide fast, one-step relief for emergency first aid when someone has been burned. The dressings stop burn progression and cool the skin, relieve pain and protect against airborne contamination. It takes no special training to give fast and effective emergency burn care. It can be used to extinguish flames and smoldering clothes.

G3630P - FIRE BLANKET 91x76 cm  
G7260C - TUBE FIRE BLANKET 183x152 cm

---

## WATER-JEL BURN DRESSING KITS



### MINI BURN KIT - VAL40515

Mini burn kit for burns treatment contained in a tearproof water-resistant orange pouch, featuring: 2 WATER-JEL burn medication gauze 10x10cm, 1 WATER-JEL burn medication gauze 20x20cm, 2 small towelettes for minor burns, 3,5ml each, 1 pair of sterile gloves, size M, 1 elastic medication bandage, 1 emergency scissors with plastic handles.

Dimensions: about 25x7x20(h) cm.

Weight: about 0,7kg.



### BURN KIT - VAL40511

Burn kit, burn treatment set in equipped tear-proof impermeable material, featuring: 1 burn cover 91x76 cm, 1 face cloth WATER-JEL 30x40 cm, 2 burn dressing WATER-JEL 10x10 cm, 1 WATER-JEL 80 ml soft burns bottle, 1 burn dressing WATER-JEL 20x45 cm, 1 hand burn dressing WATER-JEL 20x55 cm, 12 – 4 gm WATER-JEL sachets, 2 pairs of sterile gloves, 1 Isothermic blanket, 2 elastic medication bandages, 1 emergency scissor, plastic handle, 2 instant ice packs, 2 non-injectable sodium chloride 250 ml bottles. Dimensions: 48x16x26(h) cm. Weight 5,8 kg.



### WATER-JEL AMBULANCE KIT - B-BK-LG-WJE

Complete ambulance kit in pouch featuring: 1 Water-Jel wrap pouch 91x76 cm, 1 Water-Jel face burn dressing 30x40 cm, 1 Water-Jel burn dressing 20x45 cm, 2 Water-Jel burn dressing 10x40 cm, 1 Water-Jel hand burn dressing 20x55 cm, 4 bandages, 1 scissors.



### WATER-JEL CRITICAL KIT - B-BKCK-SM-WJE

Complete critical kit in pouch featuring: 1 x water-jel face burn dressing 30x40 cm, 2 x water-jel hand burn dressing 20x55 cm, 2 x water-jel burn dressing 10x10 cm, 2 bandages, 1 scissors.




## ITALY MEDICATION CASE D.M. 388/2003 - D.LGS. 81/2008 ANNEX 1



### ITALY MEDICATION CASE ANNEX 1

Medication case complete with:  
5 pairs of disposable sterile gloves; 1 anti-squirt visor; 1 litre of cutaneous Povidone-iodine with 10% iodine; 1.5 litre of saline solution; 10 packages of sterile gauzes 10 x 10 cm – single packaging; 2 packages of sterile gauzes 18 x 40 cm – single packaging; 2 disposable sterile sheets; 2 disposable sterile medical forceps; 1 medium-sized elastic bandage; 1 package of hydrophilic cotton; 2 packages of various sized plasters ready-to-use; 2 rolls of high plaster 2,5cm; 1 pair scissors; 3 tourniquets; 2 packages of instant ice; 2 disposable sanitary waste bags; 1 thermometer; 1 sphygmomanometer; 1 stethoscope.  
Dimensions: about 39,5x13,5x26(h) cm.  
Weight: about 5,8 kg.

VAL45180 - MEDICATION CASE D.L. 81/2008 ANN. 1  
VAL45184 - MEDICATION CABINET D.L. 81/2008 ANN. 1  
VAL45150 - MEDICATION BAG 33 D.L. 81/2008 ANN. 1  
MED98878 - MEDICATION PACKAGE D.L.81/2008 ANN. 1 BASE  
MED98874 - MEDICATION PACKAGE D.L.81/2008 ANN. 1 COMPLETE





## D.LGS. 81/2008 ANNEX 2



## ITALY MEDICATION CASE ANNEX 2

Medication case complete with:  
2 pairs of disposable sterile gloves; 125 ml of cutaneous Povidone Iodine with 10% iodine; 250 ml of saline solution; 3 packages of sterile gauzes 10 x 10 cm – single packaging; 1 packages of sterile gauzes 18 x 40 cm – single packaging; 1 disposable sterile medical forceps; 1 package of hydrophilic cotton; 1 package of various sized plasters ready-to-use; 1 roll of high plaster 2,5 cm; 1 roll of elastic edged bandage – height 10 cm; 1 pair scissors; 1 tourniquet; 1 package of instant ice; 1 disposable sanitary waste bags;  
Dimensions: 25x9x19,5 C cm.  
Weight : about 1,3 kg.

VAL45170 - MEDICATION CASE D.L. 81/2008 ANN. 2  
VAL45174 - MADICATION POUCHE D.L. 81/2008 ANN. 2  
MED98870 - MEDICATION PACKAGE D.L.81/2008 ANN.2



## BOLZANO MEDICATION CASE D.P.P. 25/2005 ANNEX 1

---

### BOLZANO MEDICATION CASE ANNEX 1

---



Medication case complete with:  
5 pairs of disposable sterile gloves; 1 anti-squirt visor; 1 litre of cutaneous Povidone Iodine with 10% iodine; 1.5 litre of saline solution; 10 packages of sterile gauzes 10 x 10 cm – single packaging; 2 packages of sterile gauzes 18 x 40 cm – single packaging; 2 disposable sterile sheets; 2 disposable sterile medical forceps; 1 medium-sized elastic bandage; 1 package of hydrophilic cotton; 2 packages of various sized plasters ready-to-use; 2 rolls of high plaster 2,5cm; 1 pair scissors; 3 tourniquets; 2 packages of instant ice; 2 disposable sanitary waste bags; 1 thermometer; 1 sphygmomanometer; 3 medication packages - 8 cm; 3 medication packages - 10 cm; 3 elastic edged bandage - size 8 cm; 3 elastic edged bandage - size 10 cm; 1 isothermic blanket; 1 protection handkerchief for artificial ventilation; 250 ml bottle of disinfectant; 1 stethoscope.  
Dimension: about 39,5x13,5x26 cm.  
Weight: about 6,2 kg.

---

VAL45186 - BZ MEDICATION CASE D.L. 81/2008 ANN. 1  
VAL45188 - BZ MEDICATION CABINET D.L.81/2008 ANN. 1  
VAL45152 - BZ MEDICATION BAG 33 D.L. 81/2008 ANN. 1  
MED98868 - BZ MEDICATION PACKAGE D.L.81/2008 BASE ANN. 1  
MED98866 - BZ MEDICATION PACKAGE D.L.81/2008 COMPLETE ANN. 1

---



## BOLZANO MEDICATION CASE D.P.P. 25/2005 ANNEX 2

---

### BOLZANO MEDICATION CASE ANNEX 2

---



Medication case complete with:  
2 pairs of disposable sterile gloves; 125 ml of cutaneous Povidone Iodine with 10% iodine; 250 ml of saline solution; 3 packages of sterile gauzes 10 x 10 cm – single packaging; 1 packages of sterile gauzes 18 x 40 cm – single packaging; 1 disposable sterile medical forceps; 1 package of hydrophilic cotton; 1 package of various sized plasters ready-to-use; 1 roll of high plaster 2,5 cm; 1 roll of elastic edged bandage – height 10 cm; 1 pair scissors; 1 tourniquet; 1 package of instant ice; 1 disposable sanitary waste bags; 2 medication packages - 8 cm; 2 medication packages - 10 cm; 2 elastic edged bandage - size 8 cm; 2 elastic edged bandage - size 10 cm; 1 isothermic blanket; 1 protection handkerchief for artificial ventilation.  
Dimensions: about 25,5x9x19,5 cm.  
Weight: about 1,5 kg.

---

VAL45176 - BZ MEDICATION CASE D.L. 81/2008 ANN.2  
VAL45178 - BZ MEDICATION POUCHE D.L. 81/2008 ANN.2  
MED98864 - BZ MEDICATION PACKAGE D.L.81/2008 ANN.2

---








A wide range of reflective and adhesive emblems and coat of arms is available, only few are showed below. Contact our commercial department for any further information.



**REFLEX CRI VDS - STE140110**

Round reflex emblem CROCE ROSSA - VOLONTARIO DEL SOCCORSO. Without velcro. Dimensions: Ø 7 cm.



**REFLEX CRI PIONIERI - STE140115**

Round reflex emblem CROCE ROSSA - PIONIERI. Without velcro. Dimensions: Ø 7 cm.



**REFLEX P.C. P.A. - STE140140**

Triangular reflex emblem Protezione Civile – Pubbliche Assistenze. Without velcro. Dimensions: 11x11x11 cm.




**INTERNATIONAL SYMBOLS - STE141500**

Adhesive reflex international symbolisms. 12 different pieces packages. Dimensions: 4x4 cm.



**REFLEX PROTEZ. CIVILE - STE140150**

Round reflex emblem PROTEZIONE CIVILE NAZIONALE. Without velcro. Dimensions: Ø 7 cm.



**REFLEX CRI - STE140220**

Round reflex emblem Croce Rossa. Without velcro. Dimensions: Ø 22 cm.

# OB LIFE SUPPORT 300 TRAINER DEFIBRILLATOR

## DEFIBRILLATOR TRAINER

The defibrillator OB Life Support 300® TRAINER is a simple cost-effective simulator of the defibrillator OB Life Support 300 that is designed to be used exclusively for training. There is no high voltage in the device, and training is risk free.



## 8 SCENARIOS

The device provides several simulations or scenarios (from 0 to 8) that aim to help to familiarize with the OB Life Support 300 and shows that the users have the knowledge required to use it in case of emergency. The device also allows the trainer to work using an infrared remote control to select a scenario or to intervene in real time. OB Life Support 300® TRAINER allows the simulated delivery of electric shocks. This device can be used on all types of mannequins, as it does not deliver an electric shock. It may only be used with the training electrodes commercialized by Oscar Boscarol Srl.



## ACCESSORIES

The basic equipment of the OB Life Support 300® TRAINER consists of:  
3 sets of disposable training electrodes  
1 reusable cable with the connector between the electrodes and the device. The connector of the cable is fitted with a magnet allowing the device to recognize the presence of the electrodes in the device without any electrical contact.  
1 remote control: it is programmed at delivery and ready for use.



DAE290 - OB LIFE SUPPORT TRAINER  
DAE312 - AED TRAINER PADS  
0-21-0006 - REUSABLE CABLE FOR TRAINING ELECTRODES



## CARDIO FIRST ANGEL AND MANIKIN

---



### CARDIO FIRST ANGEL - 2-62-0000

---

The Cardio First Angel allows everyone to give an efficient cardiac massage. It requires no electricity or battery as it works purely mechanically. Its applicability is shown in uncomplicated and international comprehensible icons. The specific design of the Cardio First Angel guarantees the correct positioning on the patient's breast, so the cardiac massage works exactly where it is needed. A special foam creates a soft contact surface and absorbs sudor. A special spring signals the perfect pressure of 41 kg with a clear click-sound, also audible with the following release of the spring. If the first-aider is geared to these two click-sounds he will automatically reach the required frequency of round about 100 compressions per minute.



### CPR ADULT MANIKIN - 2-19-0004

---

An adult, realistic, lightweight and inexpensive manikin; perfect for an high quality CPR training and allows an immediate check of the RCP practiced maneuvers, (frequency and chest compression control). The control of lights allows the student to practice the compressions at a rate of 100 times per minute. Click sound for the control of the depth: in addition to CPR monitor, students can hear a "click" whenever the chest is pressed to the correct depth, allowing you to experience the needed strength to practice real chest compressions. Quick and easy to prepare: thanks to the chest design with a shell opening it is possible to replace easily and quickly the airways. Supplied with a comfortable carrying bag that can be used as a mat during training. Contents: a torso, a head, 10 face protection, two AA batteries (already installed) and a carrying bag.

## MANIKINS



### LITTLE ANNE - 120-01050

A realistic airway construction means that the airway naturally remains obstructed without the use of proper head tilt, chin lift or jaw thrust. This leads to more realistic training as chest rise will be seen when correct ventilations are made. Little Anne's anatomically correct torso and sternal notch allow the student to practice identification of all anatomical landmarks relevant to conducting adult CPR training. Realistic oral and nasal passages allow for mouth-to-mouth ventilation. Due to the highly realistic chest movement and compression resistance feedback devices, such as the CPRmeter, can be used to enhance training and ensure real time feedback and quality debriefing is provided. Audible feedback reinforces correct compression depth. An optional "clicker" feature signals when the correct compression depth is achieved. Economical disposable airways for quick and easy clean-up.



### RESUSCI ANNE FIRST AID - 170-00150

High-quality, durable construction enable training in a variety of environments. Modular design allows various limbs to be attached to produce a more realistic training experience. Easy-to-use and maintain faces and disposable airways give you a cost-efficient solution for First Aid training. Natural obstruction of the airway allows students to learn the important technique of opening the airway. Head tilt/chin lift and jaw thrust allows students to correctly practice airway maneuvers. Occluded airway with hyperextension stresses proper head position. Realistic resistance for chest compressions. Articulating arms, legs, and head allow realistic weight and handling of an adult during obstructed airway maneuvers. Carotid pulse simulation to realistically check for pulse. Disposable non-rebreathing airways are suitable for use by more than one student during class and are quick and easy to change after each training session.

## OTHERS



### FACE PROTECTIVE SHEETS - DID150198

Ventilation's handkerchief with filter protection. It allows the mouth to mouth ventilation in absolute safety. Transparent PVC sheet, single-use, non-toxic. Suitable for training on manikins.

Dimensions: about 19,5x25 cm.  
Package of 50 pieces in roll.




### PAIR OVERSHOES IN TNT - 18COP100

Overshoes with elasticated ankle  
Made in green TNT.  
Set of 2 pieces.



### DISPOSABLE SUIT PROTECTION - 18TUT001

Suit with TNT hood in white spunbond, not sterile and with elastic cuffs, ankles and waist  
Ideal to cover the operating room and department staff: barrier effect to liquids and bacteria. High drapability and comfort and total absence of latex. Compliant with UNI EN 556.




### TRIAGE TAGS - PC190000

Triage cards made of cardboard. The triage card is composed not only from the standard part reporting the state and the practiced therapy but also by 4 stubs, perforated tear-off, with 4 different colors, indicating the severity of the patient. Package of 25 pieces  
Dimensions: 10,5x20,5 cm.



### SAFETY SIGN OF FIRST AID - ACC170550

Safety Sign of First Aid (white cross on a green background).  
Made of aluminum.  
Dimensions: 25x25 cm.



### SAFETY SIGN OF AED - ACC170551

Safety Sign AED.  
Made of aluminum.  
Dimensions: 25x25 cm.

## TORCHES AND FORCE-OPEN SET

---



### SNAPLIGHT LIGHTSTICK - ACC170044

---

Bright torch 15 cm long. Ideal emergency light for every occasion. Different usages in the night and in total darkness. Available in brilliant red light. 12 hours duration.



### TORCH "TRAFFIC LED BATON" 12 LED

---

Electronic warning torch with blinking light in red or yellow LED. More than any conventional torch for high quality and for the autonomy of 400 hours. Waterproof and very resistant to impacts. Dimensions: about 4x54 cm.

---

ACC170117 - 12 LED TORCH RED

ACC170111 - 12 LED TORCH YELLOW

---



### FORCE-OPEN SET FOR ROAD CRUSH - ACC170151

---

Force-open set for road crushes in tearless bag with inner shockproof basin. Complete with 3 transparent pockets, handle and shoulder belt for transport. Featuring: 1 crowbar, 1 hatchet with plastic handle, 1 quality shear/ripper for battery cables cutting, 1 pair of accident-proof gloves, 1 Lightstick flashlight.

Dimensions: about 48x26x16 cm.

Weight: about 7,6 kg.



## OSCAR BOSCAROL SRL

Via Enzo Ferrari 29, 39100 Bolzano - ITALY  
VAT NUMBER: 01458460217

**PH.** +39 0471932893

**FAX** +39 0257760140

**@** info@boscarol.it

**@** info@pec.boscarol.it

**WWW.** boscarol.it



OscarBoscarolSrl

## WORKING SCHEDULE

Monday - Tuesday

09:00-12:30 and 14:00-18:00

Friday: 09:00-13:00



## SELLING CONDITIONS

---

### 1. OVERALL

General selling conditions listed hereupon are to be considered valid for every purchasing made to our company Oscar Boscarol S.r.l. and make a fundamental requirement for every purchase order. Possible variations to such clauses rather than particular conditions must be explicitly agreed in written form.

### 2. ORDERS

Purchase orders are to be sent to the following: fax 0039 0257760140; email [info@boscarol.it](mailto:info@boscarol.it) or on-line via our web site [www.boscarol.it](http://www.boscarol.it) (against registration). All orders received will be processed by our commercial department and must be approved by the Oscar Boscarol S.r.l. through a regular order confirmation. Orders will have to be compliant to the packaging details stated on the price list.

### 3. TRANSPORT

Transport occurs to the customer's account with the amount stated on the invoice or as assigned port (the customer arrange the shipment with his own forwarder). The customer has to check the packaging status at the moment of delivery: should any anomaly be noticed the goods should be accepted with note to the forwarder itself and to the Oscar Boscarol S.r.l.

### 4. DELIVERY

Delivery terms are stated on our order confirmations and might be subject to variations which are communicated to the customer. Order processing is strictly subject to its written confirmation from the customer.

### 5. PRICES

Price listed hereupon are to be considered in Euro, VAT and tax excluded. OSCAR BOSCAROL S.r.l. owns the right to perform changes on the prices in case raw materials and/or components price increase occurred.

### 6. RETURN POLICY AND COMPLAIN

Possible complains must be sent in written form to Oscar Boscarol S.r.l. within 8 days from goods delivery. The return of the goods must be previously authorized by Oscar Boscarol S.r.l.

The authorized returning goods must be addressed to the Oscar Boscarol S.r.l. warehouse at the conditions stated on the return authorization and needs to be shipped together with a transport document stating the reasons of the return together with purchasing references. Return goods will be replaced and/or credited after check-in inspection.

### 7. SIZE AND DRAWINGS

Indications related to dimensions and drawings on each publication issued by Oscar Boscarol srl have a mere informative purpose and can be modified without notice.

### 8. PRODUCTS STOCK AVAILABILITY:

All the products marked with (\*\*\*) are to be intended as AVAILABLE UNTIL STOCK AVAILABILITY.

### 9. PAYMENT

Payment must be performed in due course and for the amount agreed. Failing to do so might result in the temporarily suspension of supplies and start bank interests appliance.

### 10. Refer to Bolzano forum for every controversy that might occur.

Please be advised that all the information featured in this catalogue are merely indicative and subject to continuous changes. For any updated information please take a look at our website [www.boscarol.it](http://www.boscarol.it) or contact our commercial department.

Oscar Boscarol company reserves the right to update what stated in the technical and advertising files without any prior notice. Please contact our offices via telephone +39 0471 932893 or via email to double check such information. By communicating your data to Oscar Boscarol you are automatically subscribed to our newsletter from which you can cancel at any time by clicking the "Cancel subscription" link at the bottom of every newsletter.







OSCAR BOSCAROL SRL  
Via Enzo Ferrari 29, 39100 Bolzano ITALY - Tel. +39 0471932893 - Fax +39 0257760140  
info@boscarol.it - www.boscarol.it - www.emergencymedicalproduct.com