



**BOR87040**

### **OXIBAG BLAU\***

- Padded oxygen cylinder case
- Waterproof and tear-proof high-strength PVC fabric
- The shape and padding allow the oxygen cylinder to be safely contained, protected and transported
- Can hold 1 or 2 l oxygen cylinders
- Dimensions: approx. 54x23x24 cm
- Weight: approx. 1.1 kg

*DIR/QM – Rev. 01\_26/07/2023*

Für weitere Informationen u/o technische Daten zum Produkt wenden sie sich an Oscar Boscarol srl