



## ETAN45C02

### TRACHEALTUBUS MIT CUFF, 4,5 MM, ORAL/NASAL

- Disposable, cuffed, oro-nasal tracheal tube
- Made of thermosensitive polyvinyl chloride (PVC)
- Intended for the following uses:
  - positive pressure ventilation;
  - induction of anaesthesia if it is not possible to use a face mask;
  - removal/protection of secretions from the trachea and bronchi during suction.
- 15M connector
- Graduated with depth marking to check positioning in relation to the vocal cords
- Radiopaque line for clear identification under X-rays
- Equipped with:
  - Low pressure, high volume, atraumatic cuff
  - Blue indicator balloon with Luer fitting and one-way valve
  - Atraumatic, rounded edge Murphy eye
  - Atraumatic, pre-curved Magill tip in compliance with applicable ISO/EN standards.

DIR/QM – Rev. 01\_31/03/2025

Für weitere Informationen u/o technische Daten zum Produkt wenden sie sich an Oscar Boscarol srl