



**BOR80290**

## **MINI RESCUE BAG RED EMPTY**

- Compact and practical emergency bag
- Made in tearproof and waterproof material
- Equipped with:
  - handles and shoulder strap
  - 2 windows for the insertion of identity cards/logos
  - lower part with inner shockproof basin
  - upper part with elastics
- Size: about 39 x 18 x 25h cm
- Weight: about 1 kg
- Supplied empty

*DIR/QM – Rev. 02\_08/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.