



BOR80290

MINI RESCUE BAG RED EMPTY

- Compact and practical emergency bag
- Made in tearproof and waterproof material
- Equipped with:
 - handles and shoulder strap
 - 2 windows for the insertion of identity cards/logos
 - lower part with inner shockproof basin
 - upper part with elastics
- Size: about 39 x 18 x 25h cm
- Weight: about 1 kg
- Supplied empty

DIR/QM – Rev. 02_08/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.