



BOR84065

E-HOLSTER RED COMPLETE*

- Complete belt E-Holster
- Made in tearproof and water-resistant red material
- 5 compartments with pockets for gloves and for Pen Light
- Belt loop of about 5,5 cm
- Dimensions: about 14x3x19(h) cm
- Complete with:
- 1 first aid scissors
- 1 anatomic forceps
- · Cocher forceps
- 1 Pen Light
- 1 Tourniquet

DIR/QM - Rev. 02_08/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.