



**BOR87020N**

## **EMERGENCY RUCKSACK MASTERMED 2.0**

- Mastermed Rucksack with Edition 2019 internal structure made of red nylon
- The tear-proof and waterproof fabric allows safe use even in precarious weather conditions
- The two front pockets are made with transparent internal compartments that are easy to open
- The front pocket is equipped with elastic bands for the fixing medical devices of different sizes
- The shoulder straps are padded and equipped with reflective strips along their entire length
- Rucksack with 8 coloured inner bags with transparent front, 4 of which with handle
- The part in contact with the back of the user is equipped with cushions
- Dimensions: about 35x30x53cm
- Volume: about 43 litres
- Weight: about 4.3kg
- Maximum weight capacity: 40kg
- It can be washed by hand at a maximum temperature of 30° C
- Supplied empty but can be equipped with medical devices upon customer request
- Mastermed Rucksack complies with REACH regulations but is not a medical device

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.