



**BOR87020N**

## **EMERGENCY RUCKSACK MASTERMED 2.0**

- Mastermed Rucksack Edition 2019 with internal structure made of red nylon
- Tear-proof and waterproof fabric allows safe use even in precarious weather conditions
- Equipped with:
  - two front pockets made with transparent internal compartments that are easy to open
  - front pocket with elastic bands for the fixing medical devices of different sizes
  - padded shoulder straps with reflective strips along their entire length
  - 8 coloured inner bags with transparent front, 4 of which with handle
  - shock-absorbing cushions on the part in contact with the back
- Dimensions: about 35x30x53cm
- Volume: about 43 litres
- Weight: about 4kg
- Maximum weight capacity: 40kg
- It can be washed by hand at a maximum temperature of 30° C
- Supplied empty but can be equipped with medical devices upon customer request

*DIR/QM – Rev. 02\_08/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.