



**BOR89010**

## **PROFI BAG EMPTY**

- Emergency empty Profi Bag
- Made of impermeable and tearproof red fabric material
- Equipped with:
  - handle and shoulder straps
  - modular interior
  - 1 frontal pocket
  - 2 wide lateral external pockets
  - 1 inside pocket with transparent portion on top.
- Arranged to hold 1-2 lts oxygen cylinder
- Dimensions: about 60x25x27(h) cm
- Weight: about 1,4 kg

*DIR/QM – Rev. 02\_13/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.