



BOR89060

TRAUMA BAG 30 EMPTY

- Functional, capacious, resistant and compact emergency bag
- Made of tearproof and impermeable red fabric material
- Totally removable cover
- Easy to carry either by hand, on shoulder or as a backpack
- Equipped with:
 - 4 refractive straps
 - 2 external lateral pockets
 - 3 external lateral pockets with elastics and additional pockets
 - 5 inner modular spaces
 - 2 inner removable pouches with transparent frontal
- Dimensions: about 63x37x31(h) cm
- Weight: about 2,3 kg
- Supplied empty

DIR/QM – Rev. 02_13/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.