



**BOR89060**

## **TRAUMA BAG 30 EMPTY**

- Functional, capacious, resistant and compact emergency bag
- Made of tearproof and impermeable red fabric material
- 4 refractive straps
- Totally removable cover
- Easy to carry either by hand, or on shoulder
- 2 external lateral pockets
- 3 external lateral pockets with elastics and additional pockets
- 5 inner modular spaces
- 2 inner removable pouches with transparent frontal
- Dimensions: about 63x37x31(h) cm
- Weight: about 2,3 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.