



BOR89300

INTUBATION BAG EMPTY

- Bag for the intubation
- Interior logic structure with elastics to arrange the intubation set
- Made of blue colour tear-proof impermeable material
- Equipped with:
 - 1 anterior refractive strip
 - 2 posterior Velcro strips for anchorage on rucksacks or rescue bag
- Dimensions: about 32x22x7,5(h) cm
- Weight: about 300 g

DIR/QM – Rev. 02_13/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.