



## VAL42000

### EMERGENCY MEDICAL CASE - TYPE A\*

- Made of anodized aluminium
- High resistance and long durability
- Inside made in washable and disinfectable thermo-format with modular system
- Arranged to accommodate oxygen cylinders (1 or 2 litres)
- Dimensions: about 52x36x17(h) cm
- Weight: about 5,7 Kg
- Supplied empty

*DIR/QM – Rev. 02\_20/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.