



**VAL43000**

### **ALUMINIUM EMERGENCY CASE - TYPE B\***

- Type B emergency medical case empty in aluminium
- Made of anodized aluminium
- High resistance and long durability
- Inside made in washable and disinfectable thermo-format with modular system
- Arranged to accommodate oxygen cylinders (1 or 2 litres)
- Removable internal basin made in thermo-format with rail ampoules holder (max 60 pcs)
- Dimensions: about 52x36x24(h) cm
- Weight: about 7,6 Kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.