



VAL43141

BOSCAROL HIGH IMPACT CASE M

- Boscarol High Impact medium orange case with very thick walls for impact resistance
- Watertight polypropylene case
- Hermetic seal and automatic pressurisation valve
- Double-step safety latch for easy opening
- Soft handle made of a special compound to improve grip and comfort
- Double padlock predisposition
- Shoulder belt predisposition
- Stacking system
- Zip along the entire profile of the case
- Nylon pins, no corrosion problems
- Cubed and embossed sponges available on request
- Arrangement for resin-coated label
- Inside dimensions: 42.6x29x15.9 cm
- External dimensions: 46.4x36.6x17.6 cm
- Weight: 2.78 kg

DIR/QM – Rev. 02_20/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.