



VAL43142

BOSCAROL HIGH IMPACT CASE L

- Boscarol High Impact large orange case with very thick walls for impact resistance
- Watertight polypropylene case
- Hermetic seal and automatic pressurisation valve
- Double-step safety latch for easy opening
- Soft handle made of a special compound to improve grip and comfort
- Double padlock predisposition
- Shoulder belt predisposition
- Stacking system
- Zip along the entire profile of the case
- Nylon pins, no corrosion problems
- Cubed and embossed sponges available on request
- Arrangement for resin-coated label
- Inside dimensions: 50x35x19.4 cm
- External dimensions: 55.5x42.8x21.1 cm
- Weight: 3.86 kg

DIR/QM – Rev. 02_20/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.