



VAL43151

EXPLORER CASE MINI 30.5X19.5X27

- Emergency waterproof and watertight case up to 10 meters deep.
- With pressurization valve
- Corrosion resistant
- Sturdy closing system with predisposition for a lead lock
- External dimensions: about 30,5x19,5x27h cm
- Internal dimensions: about 27,6x17x20h cm
- Volume: 9,3 L
- Weight: about 1,9 kg
- Supplied empty

DIR/QM – Rev. 02_20/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.