



**VAL43152**

### **EXPLORER CASE 410X330X205**

- Watertight rescue case up to 10 m depth
- Pressurisation valve
- Corrosion resistant
- Robust snap fasteners with padlock or lead-tightness
- External dimensions: 41x20.5x34 cm
- Internal dimensions: 38x18x27 cm
- Volume: 18.4 litres
- Weight: 3.1 kg
- Supplied empty

*DIR/QM - Rev. 02\_20/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.