



**VAL43162**

### **EXPLORER CASE 520X435X230**

- Watertight rescue case up to 10 m depth
- Pressurisation valve
- Corrosion resistant
- Robust snap fasteners with provision for padlock or lead
- External dimensions: 52x23x43.5 cm
- Internal dimensions: 48x20.5x37 cm
- Volume: 35.5 litres
- Weight: 4.6 kg
- Supplied empty

*DIR/QM - Rev. 02\_20/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.