



**VAL43172**

### **EXPLORER CASE 650X510X24**

- Waterproof rescue case watertight to a depth of 10 m
- Pressure valve
- Corrosion resistant
- Huge click closure suitable for optional plumb lock
- External dimensions: about 65x24x51(h) cm
- Inner dimensions: about 58x22x44(h) cm
- Volume: 56,1 litres
- Weight: about 7,1 kg
- Supplied empty

*DIR/QM – Rev. 02\_20/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.