



**VAL43174**

### **CASE EXPLORER S WITH 4 EMPTY POUCHES**

- Case Explorer P with 4 empty pouches
- Watertight to a depth of 10 m
- Pressurization valve
- Resistant to corrosive agents
- Strong click closing predisposed for padlock or lead
- 4 coloured pouches with transparent frontal and handle
- External dimensions: about 41x20,5x34(h) cm
- Internal dimensions: about 38x18x27(h) cm
- Volume: 18,4 l
- Weight: about 3,4 kg

*DIR/QM – Rev. 02\_20/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.