



**VAL43176**

### **CASE EXPLORER M WITH 5 EMPTY POUCHES**

- Case Explorer M with 5 empty pouches
- Water repellent up to 10m depth
- Pressurization valve
- Resistant to corrosive agents
- Strong click closing predisposed for padlock or lead
- 5 coloured pouches with transparent frontal and handle
- External dimensions: about 52x23x43,5(h) cm
- Internal dimensions: about 48x20,5x37(h) cm
- Volume: 35,5 l
- Weight: about 5 kg

*DIR/QM – Rev. 02\_20/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscrol Company.