



VAL43178

CASE EXPLORER G WITH 6 EMPTY POUCHES

- Case Explorer G with 6 empty pouches
- Water repellent up to 10m depth
- Pressurization valve
- Resistant to corrosive agents
- Strong click closing predisposed for padlock or lead
- 6 coloured pouches with transparent frontal and handle
- External dimensions: about 65x24x51 (h) cm
- Internal dimensions: about 58x22x44(h) cm
- Volume: 56,1 l
- Weight: about 7,6 kg

DIR/QM – Rev. 02_20/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.