



BOR87010N

EMERGENCY RUCKSACK MAXIMED 2.0

- Maximed Rucksack 2019 Edition made of red nylon
- Tearproof and waterproof material, allowing the use of the rucksack during rainy conditions
- Equipped with:
 - reflective stripes on the padded shoulder straps for better visibility, as well as on the front and side pockets
 - shock-absorbing cushions on the back area
 - 2 front pockets
 - 1 empty outer perimeter pocket and 1 outer perimeter pocket with 4 removable coloured pouches. The 4 coloured inner pouches have a transparent front and are removable
 - straps in the centre for secure fastening of the oxygen bottle
- Dimensions 48x30x52cm
- Capacity: about 54 litres
- Weight: about 4.5kg
- To be hand washed at a maximum water temperature of 30°C
- Provided empty, but with the possibility to be equipped with medical devices on request

DIR/QM – Rev. 02_08/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.