



BOR87000N-S

EMPTY VIPER RESCUE BACKPACK 2.0

- Rescue Rucksack Viper made of red nylon
- The tear-proof and waterproof fabric allows safe use even in precarious weather conditions
- Equipped with:
 - double handle (upper or lateral) for easier transport
 - front, side and rear reflective strips (on padded shoulder straps)
 - padded back
 - large front opening pocket with 2 internal elastic bands for the fixing of material
 - 3 external perimeter pockets: two little ones at the ends of the second handle and a big one on the opposite side.
 - 8 removable coloured pouches with transparent front and Velcro strips
 - internal central partition with zipper, complete with elastic bands for the fixing medical devices
- Reinforced bottom and rubber feet
- Dimensions: about 46x26x50cm
- Weight: about 3.8kg
- Maximum weight capacity: 40kg
- It can be washed by hand at a maximum temperature of 30° C
- Supplied empty but can be equipped with medical devices upon customer request

DIR/QM – Rev. 02_08/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.