



BOR87300R

BAG 33 EMPTY

- Bag 33 empty
- Made of red tearproof and impermeable fabric material
- Equipped with:
 - handle and shoulder strap
 - 1 capacious outer pocket
 - 5 inner adjustable spaces
 - Inside provided with elastics and a documents pocket
 - 1 transparent frontal pocket for association logo
- Dimensions: about 36x18x35(h) cm
- Weight: about 1,6 kg

DIR/QM – Rev. 02_13/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.