



**BOR89000R**

## **PROFI BAG MINI EMPTY**

- Emergency empty Profi Bag
- Made of impermeable and tearproof red fabric material
- Equipped with:
  - handle and shoulder straps
  - modular interior
  - 2 wide lateral external pockets
  - transparent inside pocket on top.
- Dimensions: about 36x26x20(h) cm
- Weight: about 0,7 kg

*DIR/QM – Rev. 02\_13/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.