



**BOR89040PR**

## **TRAUMA BAG 10 PEDIATRIC EMPTY**

- Compact but extremely large rescue bag
- Made of tear-proof and waterproof nylon
- The fabric has bright colours and textures suitable for use in the paediatric sector
- Completely removable and washable according to the specified conditions
- Equipped with 4 reflective strips placed on the sides of the bag
- Bag complete with handles, belt and shoulder straps for comfortable and safe transport
- The 2 external side pockets have internal spaces and elastic bands for the fixing of medical devices
- The inside of the bag is divided into 4 separate spaces that can be removed
- In the upper part there is a removable bag with a transparent front
- Dimensions: about 44x30x29cm
- Weight: about 2.1kg
- The bag is sold empty but can be equipped with medical devices upon customer request
- The Trauma Bag 10 Paediatric complies with the REACH regulations but is not a medical device

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.