



**BOR87006**

## **EMPTY PVC VIPER RESCUE BACKPACK 2.0**

- Viper rescue rucksack made of red PVC
- The tear-resistant and waterproof fabric allows for safe use even in adverse weather conditions.
- double handle (top or side) for easier transport
- reflective strips on the front, sides, and back (on the padded shoulder straps)
- cushioning pads on the part in contact with the back
- large front pocket with full opening equipped with 2 internal elastic strips for securing materials
- 3 external perimeter pockets: 2 smaller ones at the ends of the second handle on one side
- and 1 large one on the other side
- 8 removable colored pouches with transparent front and Velcro-type strips
- internal central divider with zipper, complete with elastic bands for securing medical devices
- reinforced bottom and rubber feet
- Dimensions: approx. 48x30x51 cm
- Weight: approx. 3.65 kg
- Maximum transportable weight 25 kg (please note DL.81)
- Can be hand washed at a maximum temperature of 30° C
- Supplied empty but can be equipped with medical devices at the customer's request

*DIR/QM – Rev. 01\_10/02/2025*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.