



BOR87003

EMERGENCY RUCKSACK MINI VIPER PVC

- Red PVC tearproof and water-proof material
- Equipped with:
 - double upper and side handle for easier transport
 - large front opening pocket with 2 internal elastic bands for the fixing of material
 - 2 small side pockets and a large one on the outside
 - front, side and rear reflective strips (on padded shoulder straps)
 - shock-absorbing cushions on the back side
 - 4 removable coloured pouches with transparent front and Velcro strips
 - reinforced bottom and rubber feet
- Dimensions: about 44x26x44cm
- Weight: about 2.7kg
- Maximum weight capacity: 20kg
- It can be washed by hand at a maximum temperature of 30° C
- Supplied empty but can be equipped with medical devices upon customer request

DIR/QM – Rev. 01_23/09/2025

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.