



126140107191

YANKAUER RIGID SUCTION PROBE COMPLETE*

- Yankauer rigid suction probe, complete of 180 cm long connecting tube with female connection
- Patient aspiration using only one hand
- Rigid cannula instantly ready and easy to use
- Vacuum regulation by the thumb of the hand
- Suitable for every kind of medical secretion suction unit
- Singularly sterile packaged
- EC medical device

DIR/QM - Rev. 02_18/07/2023

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.