



VAL43000

ALUMINIUM EMERGENCY CASE - TYPE B*

- Made of anodized aluminium
- High resistance and long durability
- Inside made in washable and disinfectable thermo-format with modular system
- Arranged to accommodate oxygen cylinders (1 or 2 litres)
- Removable internal basin made in thermo-format with rail ampoules holder (max 60 pcs)
- Dimensions: about 52x36x24(h) cm
- Weight: about 7,6 Kg
- Supplied empty

DIR/QM – Rev. 02_20/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.