



**BOR81146N**

## **EMERGENCY RUCKSACK ATLS MODULE 2.0**

- ATLS modular rescue backpack
- Made of tear-proof and waterproof Nylon
- Equipped with:
  - reflective strip on the padded shoulder straps for improved visibility
  - transparent front window suitable for displaying of a logo
  - large front pocket for safe storage of even large medical devices
  - soft back pads to reduce back fatigue
  - shoulder straps with sternal closure for easier carrying
  - rigid, shockproof polypropylene tray with 4 rigid subdivisions for storing bags
  - 4 removable coloured inner pouches with transparent front
- Dimensions: approx. 32x25x53 cm
- Volume: approx. 35 l
- Weight: approx. 3.3 kg
- Maximum loadable weight 30 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.



- Sold empty but can be equipped according to specific customer requests

DIR/QM – Rev. 02\_08/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.