



BOR89015

TRAUMA BAG 25 PVC*

- Functional, compact and capacious rescue bag
- Made in tear-proof and waterproof red PVC fabric
- Equipped with:
 - 5 yellow reflective strips
 - handles and shoulder straps
 - 3 external side pockets
 - 1 exterior front pocket equipped with elastic bands and additional pockets, and with an external lodgement for the association logo
 - 1 exterior top pocket with additional logo slot
 - 5 internal modular spaces
 - 1 inner bag with transparent front
- Reinforced bottom with support feet
- Dimensions: approx. 48x28x25 cm
- Weight: approx. 2.1 kg
- Supplied empty

DIR/QM – Rev. 02_13/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.