



**BOR89015**

## **TRAUMA BAG 25 PVC**

- Functional, compact and capacious rescue bag
- Made in tear-proof and waterproof red PVC fabric
- Equipped with:
  - 5 yellow reflective strips
  - handles and shoulder straps
  - 3 external side pockets
  - 1 exterior front pocket equipped with elastic bands and additional pockets, and with an external lodgement for the association logo
  - 1 exterior top pocket with additional logo slot
  - 5 internal modular spaces
  - 1 inner bag with transparent front
- Reinforced bottom with support feet
- Dimensions: approx. 48x28x25 cm
- Weight: approx. 2.1 kg
- Supplied empty

*DIR/QM – Rev. 02\_13/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.