



VAL40905

NEWBORN RESCUSC. CASE*

- ABS resuscitation case with wall anchorage system containing:
 - 1 pvc disposable resuscitator bag
 - 1 reservoir and 1 oxygen tube
 - 3 disposable ventilation masks
 - 5 Guedel airways
- Dimensions: about 39,5x13,5x26(h) cm
- Weight: about 2,2 kg

DIR/QM – Rev. 02_30/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.