



IMM120214

ARM VACUUM SPLINT BASIC+TRANSPORT BAG

- Anatomic conformed vacuum splint suitable for leg and ankle
- Fast and painless immobilisation of all types of fractures, even compound fractures
- Immobilises but does not compress and allows normal blood circulation
- High-strength fabric, easily washable
- Equipped with:
 - plastic suction valve positioned on the outer side
 - Tear-proof and waterproof fabric protection bag with quick-release fastening system
 - Radio-transparent
 - Dimensions: approx. 68x38x53 cm
 - Weight: approx. 0.8 kg
- EC Medical device

DIR/QM – Rev. 02_21/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.