



1333124

SILICONE RESUSCITATOR ADULT SIZE

- Transparent silicone manual resuscitator, autoclavable up to 121°C
- Volume 1700 ml
- Complete with 360° rotating connection patient's valve
- Safety valve for overpressure regulation
- Additional inlet for dynamometer
- PEEP valve inlet
- Supplied with mask no. 5, reservoir and oxygen tube
- CE certified in class IIa

User manual available on the [eIFU](#) webpage

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.