



### **BOR87043**

### **OXIBAG CORDURA RED 3\***

- Padded oxygen cylinder case
- High-strength, tear-proof nylon fabric
- The shape and padding allow the oxygen cylinder to be contained, protected and transported with absolute safety
- Can hold 2 or 3 l oxygen cylinders
- Dimensions: approx. 54x23x24 cm
- Weight: approx. 0.91 kg

*DIR/QM – Rev. 01\_26/07/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.