



RIA50816

PAEDIATRIC INTUBATION BAG COMPLETE*

- Blue fabric intubation bag
- Complete with:
 - Paediatric conventional light laryngoscope with 3 blades, 0-1-2
 - Magill's forceps
 - Pean forceps
 - Scissors
 - Endotracheal tubes
 - 3 Güdel cannulas
 - Syringe
 - Helical mouth opener
 - Bandage roll
 - Spindle
- Dimensions: approx. 32x22x7.5h cm
- Set consisting of EC medical devices

DIR/QM – Rev. 02_11/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.