



**RIA50816**

### **PAEDIATRIC INTUBATION BAG COMPLETE\***

- Blue fabric intubation bag
- Complete with:
  - Paediatric conventional light laryngoscope with 3 blades, 0-1-2
  - Magill's forceps
  - Pean forceps
  - Scissors
  - Endotracheal tubes
  - 3 Güdel cannulas
  - Syringe
  - Helical mouth opener
  - Bandage roll
  - Spindle
- Dimensions: approx. 32x22x7.5h cm
- Set consisting of EC medical devices

*DIR/QM – Rev. 02\_11/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.