



RIA50820

INTUBATION BAG CONVENTIONAL LIGHT*

- Intubation bag complete of:
 - 1 conventional light laryngoscope with 3 blades (1-2-3)
 - 1 Magill forceps
 - 1 Pean forceps
 - 1 scissors
 - 5 endotracheal tubes
 - 3 Guedel airways
 - 2 Wendl airways
 - 1 syringe
 - 1 mouth opener
 - 1 plaster roll
 - 1 stylet
- Dimensions: about 32x22x7,5(h) cm
- Set consisting of EC-marked medical devices

DIR/QM – Rev. 02_11/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.