



RIA50825

PEDIATRIC INTUBATION BAG OPT. FIB. COMP*

- Blue fabric intubation bag containing:
 - 1 paediatric fibre-optic laryngoscope with 3 blades, sizes 0-1-2
 - 1 Magill's forceps
 - 1 Pean forceps
 - 1 scissor
 - 4 endotracheal tubes
 - 3 Güdel cannulas
 - 1 syringe
 - 1 helical mouth opener
 - 1 roll of plaster
 - 1 mandrel
- Dimensions: approx. 32x22x7.5 cm
- Set consisting of CE-marked medical devices

DIR/QM - Rev. 02_11/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.