



**RIA53000**

### **C.L. LARYNGOSCOPE SET 3 MC INTOSH BLADES**

- Metal conventional light Mc Intosh Laryngoscope
- Consisting of 3 Mc Intosh-type blades (1-2-3), standard handle, finely finished
- Autoclavable
- Supplied with a practical roll case
- The handle is designed to contain 2 C batteries (not included)

*DIR/QM – Rev. 02\_11/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.