



RIA53200

LARYNGOSCOPE O.F. 3 BLADES MC 1/2/3*

- Mc Intosh fiber optic laryngoscope
- Made in stainless steel
- 2,5V halogen lamp set inside the blade for an easy sterilization
- · Accurately finished off
- 3 Mc Intosh F.O. blades (1-2-3), standard handle
- Supplied with a practical roll case
- The handle is designed to contain 2 C batteries (not included)

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.