



**RIA53200**

**LARYNGOSCOPE O.F. 3 BLADES MC 1/2/3\***

- Mc Intosh fiber optic laryngoscope made in stainless steel
- Consisting of 3 F.O. Mc Intosh blades (1-2-3), standard handle, finely finished
- 2,5V halogen lamp set inside the blade for an easy sterilization
- Supplied with a practical roll case
- The handle is designed to contain 2 C batteries (not included)

*DIR/QM – Rev. 02\_11/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.