



**RIA53204**

### **LARYNGO. O.F. MC BLADES SIZE 1/2/3/4**

- Mc Intosh fibre optic laryngoscope made of stainless steel
- Consisting of 4 Mc Intosh F.O. blades (1-2-3-4), standard handle, finely finished
- 2,5V halogen lamp set inside the blade for an easy sterilization
- Supplied with a practical roll case
- The handle is designed to contain 2 C batteries (not included)

*DIR/QM – Rev. 02\_11/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.