



**RIA53304**

**O.F. LARYNGO. PED. MILLER BLADE 0/1/2\***

- Miller pediatric fiber optic laryngoscope, made of stain steel
- Consisting of 3 F.O. blades (0-1-2), small handle, finely finished
- 2,5V halogen lamp set inside the handle for an easy sterilization
- Supplied with a practical roll case
- The handle is designed to contain 2 AA batteries (not included)

*DIR/QM – Rev. 02\_11/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.