



**VBM30-04-0021**

## **CHILD QUICKTRACH**

- Urgent CONIOTOMY sterile set QUICKTRACH type
- Complete with:
  - 1 plastic canula with intubation internal-blade stylet with diameter of 2 mm
  - 1 ISO joint for the connection to resuscitation bag
  - 1 sharp conical-tip intubation stylet that avoids windpipe incision using classic scalpel and a “lock” that avoids penetrating too much in depth
  - 2 lateral wings allowing airways safe fixing

*DIR/QM – Rev. 02\_11/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.